

# ESC Swab

ESC (Easy Surface Checking) swab range includes a new generation of sampling devices for microbiological examination of surfaces, especially recommended for food and pharmaceutical industries.



## Product features:

- prefilled swab device
- rayon bud
- wide range of diluents
- leak-proof closure
- transport and storage at room temperature

Ref.	Description	Packaging
85623	ESC swab BPW + Neutralizing 1 mL	100/pack
85615	ESC swab BPW + Neutralizing 4 mL	100/pack
85603	ESC swab Buffered Peptone Water 10 mL	100/pack
85604	ESC swab Buffered Peptone Water 5 mL	100/pack
85605	ESC swab D/E Neutralizing Broth 10 mL	100/pack
85606	ESC swab D/E Neutralizing Broth 5 mL	100/pack
85614	ESC swab Demi Fraser Broth 10 mL	100/pack
85621	ESC swab Letheen Broth 1 mL	100/pack
85607	ESC swab Letheen Broth 10 mL	100/pack
85608	ESC swab Letheen Broth 5 mL	100/pack
85609	ESC swab Maximum Recovery Diluent 10 mL	100/pack
85610	ESC swab Maximum Recovery Diluent 5 mL	100/pack
85612	ESC swab Mycoplasma Transport Broth 5 mL	100/pack
85622	ESC swab Neutralizing Rinse Solution 1 mL	100/pack
85601	ESC swab Neutralizing Rinse Solution 10 mL	100/pack
85602	ESC swab Neutralizing Rinse Solution 5 mL	100/pack



## References:

- ISO 18593:2018. Microbiology of the food chain – Horizontal method for surface sampling.
- ISO 17604:2015. Microbiology of the food chain – Carcass sampling for microbiological analysis.
- ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination. Part 1: General rules for the preparation of the initial suspension and decimal dilutions.



## Liofilchem srl

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy  
Tel +39 0858930745  
[www.liofilchem.com](http://www.liofilchem.com)  
Headquarters, Manufacturing Site, International distribution

## Liofilchem, Inc.

465 Waverley Oaks Rd. Suite 317  
Waltham, MA 02452, USA  
Tel 781-902-0312  
US Distribution Center

ref. 6550219  
Rev. 2 - 04.05.2023  
© Liofilchem 2023  
printed in Italy