

Liofilchem®

Product Catalogue 2023 Clinical and Industrial Microbiology

The screenshot shows the Liofilchem website interface. At the top left is the Liofilchem logo. To its right is a search bar with the placeholder text "Search products...". Below the logo is a navigation menu with the following items: Home, Company (with a dropdown arrow), Solutions, Technical Library (with a dropdown arrow), News (with a dropdown arrow), Contact us, and Customer login.

The main content area features a large image of a laboratory setting. On the left, a person wearing blue gloves is using a pipette to transfer liquid into a multi-well plate. In the center, there is a large petri dish containing a green agar medium. To the right, another petri dish shows several circular colonies on a yellow agar medium, with labels "ZCAL 20 PR" and "ZCTL 15 PR" in red. Further right, a multi-well plate is shown with yellow liquid in the wells, and a label that reads "ComASP Collagen 0.25-16" and "Lot. 102118001 Exp. 2018/06/14".

Overlaid on the bottom left of the main image is the text "Solutions for Clinical microbiology" in bold black font, with a "Learn more" button below it.

Below the main image is a row of four service tiles:

- Download CoA**: Includes a circular icon with a document and a "Download" button.
- Solutions**: Includes an icon of three circles with a microscope and a "Read More" button.
- Literature**: Includes an icon of an open book and a "Read More" button.
- Sign up for our newsletter**: Includes a calendar icon and a "Read More" button.

The footer section is divided into four columns:

- Liofilchem S.r.l.**: Via Scozia, Zona Industriale 64026, Roseto degli Abruzzi (TE) Italy
- Company**: About us, Quality System
- Products**: Products, Literature
- Service**: Technical Sheet and MSDS, Quality Control Certificates

At the bottom left of the footer is a "Sign up for our newsletter" link. At the bottom right is a circular button with an upward-pointing arrow.



Microbiology Products
Est. 1983

In memory of

Dr. Santa Ruggieri
1962-2011

Solutions



Visit www.liofilchem.com to be up to date with our product updates, latest news, special offers and information on products and services.

The Liofilchem website offers a convenient and eco-friendly complete on-line Product Catalogue, fliers, brochures and video featured product pages.

Liofilchem Quality Control Certificates, Technical Sheets, MSDS and much more are available for download at the on-line Technical Library.

The Liofilchem on-line technical documentation service is convenient, quick and available 24/7.



LIOFILCHEM® s.r.l.

Via Scozia, Zona Industriale
64026 Roseto degli Abruzzi (Te)
Italy
Tel.: +39 0858930745
Fax: +39 0858930330
www.liofilchem.com



Since 1983 Liofilchem® s.r.l. has committed its Research and Development department to the production of reliable and high quality bacteriology products. Our dedication to Customer satisfaction and exceptional manufacturing versatility have granted Liofilchem® continuous growth, as a result of the increasing quality of our products and the worldwide distribution of the Liofilchem® brand name. Liofilchem® is present in over 130 countries with agents and distributors placing our products in clinical and industrial laboratories around the world.

Our catalogue is comprised of approximately 6000 items basically classified into two main categories, clinical and industrial microbiology. Liofilchem® manufacturing offers a wide range of items so that our customers can find in us a complete one stop supplier that is always ready to respond to the needs of the ever evolving laboratory market.



Liofilchem® produces microbial identification and susceptibility testing systems, MIC Test Strip, antibiotic discs in cartridges, dehydrated culture media, selective supplements, ready to use culture media in petri dish, tube, bottle and dip-slide, systems for the chemical and microbial screening in food, bio-indicators for sterilization process control systems. Liofilchem® also performs OEM and Private Label productions.

The complete manufacturing processes and quality control are carried out in house by Liofilchem® employees at our facilities in Roseto degli Abruzzi of total over 7000 m² area, which meet state-of-the-art standards for the diagnostic manufacturing.

The new Liofilchem®'s facility was built according to the latest international recommendations for manufacturers of In-Vitro diagnostics. Moreover, since many of our customers are Pharmaceutical Industries from all over the world, our new production processes closely resemble the highest standard of pharmaceutical companies. As a supplier of pharmaceutical industries, Liofilchem® is constantly subjected to their strictest audits and keen to satisfy their requests.

Liofilchem® is certified ISO 9001 for the Quality Management and ISO 13425 as an IVD manufacturer.

In 2016 Liofilchem began the FDA approval process for sales of MIC Test Strips in the USA and Liofilchem, Inc. has been established in Waltham, Massachusetts as sales and distribution center for the US market.



Liofilchem® s.r.l.

Via Scozia, Zona Industriale
64026 Roseto degli Abruzzi (Te)
Italy

tel: +39 0858930745

fax: +39 0858930330

www.liofilchem.com

VAT no.: IT00530130673

R.E.A. TERAMO no.79331 (CCIAA registration no.)

Bank:

BANCA DELL'ADRIATICO

Fil. di Morro d' Oro (TE)

ABI: 05748 CAB: 76940

Checking Account: 1000 00001380

IBAN: IT74 H057 4876 9401 0000 0001 380

SWIFT: IBSPIT3PI48

BIC: IBSPIT3P

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.



Continuous Improvement

Our quality system

Liofilchem is certified by TÜV Süd to ISO 9001 for the quality management, to ISO 13485 as manufacturer of IVD (in-vitro diagnostics medical device) and to the MDSAP scheme (Medical Device Single Audit Program) for the sale of IVD (in-vitro diagnostic medical devices) in Australia, Brazil, Canada, Japan and the United States.

Liofilchem has obtained the EU Quality Management System Certificate (IVDR) pursuant to Regulation (EU) 2017/746 on in-vitro Diagnostic Medical Devices, Annex IX Chapters I and III (Class C and B Devices excluding self/near-patient-testing and Companion Diagnostics).

EU Quality Management System Certificate (IVDR)



The Liofilchem's Quality Management System is certified by TÜV Süd to the Regulation (EU) 2017/746 on in-vitro Diagnostic Medical Devices for the manufacturing of Class C and B in-vitro Diagnostic Medical Devices.

ISO 13485:2016 MDSAP



Liofilchem's Quality Management System is certified by Tuv Süd America Inc. according to the MDSAP scheme for design, development, manufacture and distribution of IVD (In Vitro Diagnostic Medical Devices) in Australia, Brazil, Canada, Japan and United States.

UNI EN ISO 9001:2015



Since 1997 the implementation of the Quality System Management permits to satisfy Customer expectations, optimize resources, qualify and consolidate the company position in national and international market.

UNI EN ISO 13485:2016



Liofilchem manufactures In Vitro Medical Diagnostic Devices conform to the use destination provided for and current regulation requirements.

ISO 45001:2018




Liofilchem has implemented an "Occupational Health and Safety Management System" which fulfills the requirements of the standard ISO 45001, to control risks in this regard and to improve its performance.

UNI EN ISO 9001:2015





Liofilchem's production of hydrogen peroxide and other sanitizers is certified to ISO 9001 to make sure that quality meets the customer's expectations and standards. .

ComASP® Colistin 0.25-16 µg/mL



Compact Antimicrobial Susceptibility Panel containing colistin (polymyxin E) in 7 two-fold dilutions, that allows for testing up to four samples per panel.

Ref. 75001
Packaging: 4 panels



Chromatic CRE

Chromogenic medium for the detection of carbapenem-resistant Enterobacteriaceae directly from clinical specimens.

500 g - ref. 611619
CRE supplement - ref. 81088

20 plates - ref. 11619

Contam Swab

A new product range of swabs for the microbial screening directly from work surfaces.

Easy handling, results in 24 hours, incubation at 37°C.

Available for *Listeria*, *Salmonella*, *E.coli*-coliforms, total count, yeasts&moulds, MRSA.

AF Genital System

System for ID, count and AST of urogenital mycoplasma, presumptive ID of bacteria and *Candida* spp.

Ref. 74156
Packaging 20 tests

EnteroPluri-Test

System with twelve Agar media for Enterobacteriaceae identification.

Ref. 78618
10 tests/pack
Ref. 78619
25 tests/pack



ESC swab


ESC (Easy Surface Checking) swab is a device used for the microbiological monitoring of surfaces by swab sampling method. Available with several transport media, in accordance to (ISO 18593, ISO 17604).

100 tests/pack

Clean Test

Swab for detecting organic residues for the hygiene control of surfaces.

Ref. 86001
50 tests/pack



VTM

VTM is suitable for collection and preservation of clinical specimens containing viruses, including SARS-CoV-2 (COVID-19), *Chlamydia* spp., *Mycoplasma* spp. and *Ureaplasma* spp. The transport medium can be used with swabs and aspirates and is available in a plastic, screw cap tube.

ref. 26490
100 tubes/pack

Made in Italy

LIOFILCHEM MANUFACTURING

Our catalogue is comprised of approximately 6000 items basically classified into two main categories, clinical and industrial microbiology. Liofilchem manufacturing offers a wide range of items so that our customers can find in us a complete one stop supplier that is always ready to respond to the needs of the ever evolving laboratory market.

Liofilchem® produces systems for:

- microbial identification,
- antimicrobial susceptibility testing,
- sterilization process control,
- chemical and microbial screening in food.

Liofilchem® also manufactures ready to use culture media in petri dish, tube, bottle and dip-slide, dehydrated culture media, selective supplements.



Contents

Microbial identification



- Systems for ID and AST directly from clinical specimens
- Systems for AST
- Systems for ID and AST
- Systems for ID
- Accessories and reagent kits for microbial ID & AST systems
- Microbial staining and fixatives
- Rapid tests for sugar fermentation (dried medium in tube)
- Tests for microbial identification

Antimicrobial Susceptibility Testing



- Antibiotic discs in cartridges
- Antibiotic discs in cartridges for Veterinary use
- Antimycotic discs in cartridges
- Multodiscs with eight antibiotics
- Antimicrobial resistance detection, Agar Dilution, NP tests
- Accessories for susceptibility testing: disc dispenser, densitometer

Quality Control



- CultiControl freeze-dried microorganisms



- Quanti-CultiControl freeze-dried microorganisms

Dip-Slide



- Dip-Slide for the screening of urinary infection
- Dip-Slide for the screening of urogenital infection
- Dip-Slide for the screening of dermatological infection
- Flexible contact slide for bacteriological screening of surfaces

Food safety and environmental monitoring



- Systems for microbial identification directly from food and surfaces
- Coliforms in water testing
- Pesticide residue detection in water, food and drinks
- Antibiotic residue detection
- Swabs w/ and w/o transport medium
- Swabs for organic residues detection on surfaces, microbial detection,
- Surface swabs for microbial detection
- Easy Surface Checking, sample collection, with or without transport medium.

Contents

	Biological indicators for sterilization processes validation and monitoring	
	Self-contained for steam sterilization	
	Ampoules for steam sterilization	
	Strips for steam sterilization	
	Spores suspension for steam sterilization	
	Steel coupon for steam sterilization	
	D-value calculation	
	Strips for irradiation sterilization	
	Strips for dry heat or ethylene oxide sterilization	
	Spores suspension for dry heat or ethylene oxide sterilization	
	Self-contained for ethylene oxide sterilization	
	Steel Coupons for hydrogen peroxide sterilization	
	Poliflex coupons for hydrogen peroxide sterilization	
	Strips for hydrogen peroxide sterilization	
	Spores suspension for hydrogen peroxide sterilization	
	PVC coupons for hydrogen peroxide sterilization	
	Aluminum coupons for hydrogen peroxide sterilization	
	Glass coupons for hydrogen peroxide sterilization	
	Ready to use culture media	
	Ready culture media in 55 mm contact plate	
	Ready culture media in 60 mm petri dishes	
	Easy Dry media for water testing	
	Ready Bags	
	90 mm ready-to-use agar plates	
	Irradiated 90 mm ready-to-use agar plates	
	140 mm ready-to-use agar plates	
	Ready culture media in petri dishes with two sectors	
	Ready-to-use culture media in tubes and bottles	
Gaspak		
	Dehydrated culture media, ingredients and supplements	
	Dehydrated culture media	
	Ingredients for culture media	
	Freeze-dried selective and growth supplements	
	Liquid supplements	
	Animal blood derivates	
	Saline solutions and sputum treatment	

Contents



Liofilchem® IMMUNOLOGY

Incubators - Air sampler

Radial Immuno Diffusion plates

Latex tests

Rapid tests on Card



Turbidimetry

Coagulation

Bacterial suspension

Glucose beverage

Sterile disposable materials

Systems for presumptive ID and AST directly from clinical specimens		
Description	Packaging	Ref.
<p>A.F. Genital System System for detection, count and susceptibility testing of pathogenic urogenital microorganisms.</p> 	20 test 4 test	74156 79156
<p>A.F. Genital System – REAGENT Additional biochemical reagents for A.F. Genital System.</p>	100 test	80258
<p>Mycoplasma Transport Broth Liquid medium for transporting and storing clinical specimen for mycoplasma and other pathogenic microorganisms detection.</p>	20 vials	20158
<p>Copro System System for detection and identification of intestinal pathogenic microorganisms.</p>	40 test 8 test	71670 79670
<p>Copro System Plus System for detection and identification of intestinal pathogenic microorganisms.</p>	20 test 4 test	71675 79675
<p>Copro System AST System for the detection, identification and susceptibility testing of intestinal pathogenic bacteria.</p>	20 test 4 test	71671 79671
<p>Mycoplasma System Plus System for detection, count and susceptibility testing of urogenital mycoplasma.</p>	20 test 4 test	72592 79592
<p>Mycoplasma Transport Broth Liquid medium for transporting and storing clinical specimen for mycoplasma and other pathogenic microorganisms detection.</p>	20 vials	20158
<p>Urea-Arginine Screen Freeze-dried selective medium for the screening of urogenital mycoplasma.</p>	20 test	81079
<p>Mycoplasma System VET System for detection, count and susceptibility testing of urogenital mycoplasma of veterinary relevance.</p>	20 test	72593
<p>Pathogenic System 18 wells system for the presumptive identification of pathogenic microorganisms from pharyngeal, auricular, ophthalmic swabs, expectorate and other clinical specimens.</p>	20 test	71679
<p>Pathogenic System Double 12 wells system for the presumptive identification of pathogenic microorganisms from pharyngeal, auricular, ophthalmic swabs, expectorate and other clinical specimens.</p>	40 test 8 test	71678 79678
<p>Pathogenic System AST 32 wells system for the direct identification and susceptibility testing of pathogenic microorganisms from pharyngeal, auricular, ophthalmic swabs, expectorate and other clinical specimens.</p> 	20 test 4 test	71681 79681
<p>Urin System Chrom System for count, ID and AST directly from urine specimens.</p> 	20 test 4 test	74161 79161
<p>Urin System Plus System for enumeration, identification and susceptibility testing of urinary microorganisms.</p>	20 test 4 test	74160 79160

Systems for AST		
Description	Packaging	Ref.
SensiQuattro Gram-negative EU System for susceptibility testing of Gram-negative bacteria, 10 antimicrobial agents x 4 concentrations.	20 test 4 test	76031 79031
SensiQuattro Gram-positive EU System for susceptibility testing of Gram-positive bacteria, 10 antimicrobial agents x 4 concentrations.	20 test 4 test	76032 79032
SensiQuattro Candida EU System for susceptibility testing of yeasts, 8 antifungal agents x 4 concentrations.	20 test 4 test	76033 79033
SensiTest Gram-negative System for susceptibility testing of Gram-negative bacteria, 23 antimicrobial agents in single concentration.	20 test 4 test	76010 79010
SensiTest Gram-positive System for susceptibility testing of Gram-positive bacteria, 23 antimicrobial agents in single concentration.	20 test 4 test	76020 79020
ComASP® Antifungal Compact Antimicrobial Susceptibility Panel containing antifungal agents in various concentration ranges, that allows for testing one sample per panel.	4 panels (4 samples)	75101
ComASP® Cefiderocol (0.008-128 µg/mL) Compact Antimicrobial Susceptibility Panel containing cefiderocol in 15 two-fold dilutions, that allows for testing up to two samples per panel.	4 panels (8 samples)	75009
ComASP® Colistin (0.25-16 µg/mL) Compact Antimicrobial Susceptibility Panel containing colistin (polymyxin E) in 7 two-fold dilutions, that allows for testing up to four samples per panel.	4 panels (16 samples)	75001
ComASP® Piperacillin-tazobactam (0.008/4-128/4 µg/mL) Compact Antimicrobial Susceptibility Panel containing piperacillin-tazobactam in 15 two-fold dilutions, that allows for testing up to two samples per panel.	4 panels (8 samples)	75002
ComASP® Colistin / Piperacillin-tazobactam Compact Antimicrobial Susceptibility Panel containing the antibiotics in 15 two-fold dilutions, that allows for testing one sample per panel to colistin (0.008-128 µg/mL) and piperacillin-tazobactam (0.008/4-128/4 µg/mL).	4 panels (4 samples)	75003
ComASP® Ceftolozane-tazobactam / Ceftazidime-avibactam Compact Antimicrobial Susceptibility Panel containing the antibiotics in 15 two-fold dilutions, that allows for testing one sample per panel to ceftolozane-tazobactam (0.008/4-128/4 µg/mL) and ceftazidime-avibactam (0.008/4-128/4 µg/mL).	4 panels (4 samples)	75004
ComASP® Ceftolozane-tazobactam (0.008/4-128/4 µg/mL) Compact Antimicrobial Susceptibility Panel containing ceftolozane-tazobactam in 15 two-fold dilutions, that allows for testing up to two samples per panel.	4 panels (8 samples)	75006
ComASP® Oritavancin (0.001-16 µg/mL) Compact Antimicrobial Susceptibility Panel containing oritavancin in 15 two-fold dilutions, that allows for testing up to two samples per panel.	4 panels (8 samples)	75010
ComASP® Vancomycin / Teicoplanin Compact Antimicrobial Susceptibility Panel containing the antibiotics in 15 two-fold dilutions, that allows for testing one sample per panel to vancomycin (0.008-128 µg/mL) and teicoplanin (0.008-128 µg/mL).	4 panels (4 samples)	75005
ComASP® Vancomycin (0.008-128 µg/mL) Compact Antimicrobial Susceptibility Panel containing vancomycin in 15 two-fold dilutions, that allows for testing up to two samples per panel.	4 panels (8 samples)	75007
Customized versions System of 24, 32 or 40 wells with a customized set of antibiotic concentrations.	<i>minimum order 1000 panels</i>	<i>to be defined</i>



Systems for ID and AST		
Description	Packaging	Ref.
Copro System AST System for the detection, identification and susceptibility testing of intestinal pathogenic bacteria.	20 test 4 test	71671 79671
Integral System ENTEROBACTERIA System for <i>Enterobacteriaceae</i> identification and susceptibility testing.	20 test 4 test	71714 79714
Integral System Gardnerella System for <i>Gardnerella vaginalis</i> identification and susceptibility testing.	20 test 4 test	71724 79724
Integral System STAPHYLOCOCCI System for staphylococci identification and susceptibility testing.	20 test 4 test	71718 79718
Integral System STREPTOCOCCI System for streptococci identification and susceptibility testing.	20 test 4 test	71720 79720
Integral System YEASTS Plus System for yeasts identification and susceptibility testing.	20 test 4 test	71822 79822




Diagram of the Integral System YEASTS Plus test strip. It shows 24 wells arranged in three rows of eight. The substrates/reagents are: Row 1: 1. GLU, 2. MAL, 3. SAC, 4. LAC, 5. GAL, 6. MEL, 7. CEL, 8. IND; Row 2: 9. XYL, 10. RAF, 11. TRE, 12. DUL, 13. CHH, 14. NY, 15. AMB, 16. FIC; Row 3: 17. EON, 18. MCA, 19. CLO, 20. MIO, 21. TRT, 22. VOR, 23. FEU, 24. GEM. Below the wells is a label 'INTEGRAL SYSTEM YEASTS Plus' with fields for 'Sample' and 'Date'.

Systems for ID		
Description	Packaging	Ref.
EnteroPluri-Test System with 12 Agar media by simultaneous inoculum for Enterobacteriaceae identification.	10 test 25 test	78618 78619
 VP (KOH) Reagents VP (KOH) Reagents Droppers KOH and alpha-naphthol solutions for Voges Proskauer test. KOVAC'S - Reagent Reagent solution for indole test. EnteroPluri-Test Code Book Bacterial Code Book for Enterobacteriaceae identification using EnteroPluri-Test.	10 x 10 mL 20 x 2.5 mL 800/1000 test 20 x 2.5 mL dropper 1 book	80281 87007 80271 87001 71709
Oxi/Ferm Pluri Test System for the identification of Gram negative, oxidase positive bacteria. Oxi/Ferm Pluri Test Code Book Bacterial Code Book for Enterobacteriaceae identification using Oxi/Ferm Pluri Test.	10 test 25 test 1 book	78620 78621 71708
Anaerobe System 24 well system for the biochemical identification of anaerobic microorganisms. Nitrate Test KOVAC'S - Reagent H ₂ O ₂ Reagent	20 test	71620 88009 80271 80057
Entero System 18R System with 18 desiccated biochemical tests in wells for Enterobacteriaceae identification. Entero System 18R – Reagent Additional biochemical reagents for Entero System 18R. Entero System 18R Code Book Bacterial Code Book for enterobacteria identification using Entero System 18R.	20 test 4 test 100/200 test 1 book	71618 79618 80252 71710
Entero System 24R System with 24 desiccated biochemical tests in wells for Enterobacteriaceae identification.	20 test	71619
Listeria System 18R System with 18 desiccated biochemical tests in wells for <i>Listeria</i> spp. identification. Listeria System 18R - Reagent Additional biochemical reagents for Listeria System 18R.	20 test 4 test 100/200 test	71640 79640 80257
Staf System 18R System with 18 desiccated biochemical tests in wells for staphylococci identification. IDENTIFICATION System – Reagent Additional biochemical reagents for identification systems.	20 test 4 test 100/200 test	71630 79630 80260
Strepto System 12R System with 12 desiccated biochemical tests in wells for streptococci identification.	40 test 8 test	72560 79560

Accessories and reagent kits for microbial ID & AST systems		
Description	Packaging	Ref.
KOVAC'S - Reagent Reagent solution for indole test.	4 bottles x 25 mL 20 x 2.5 mL	80271 87001
McFarland STANDARD SET 0.5 - 1.0 - 2.0 - 3.0 - 4.0 Turbidity standard for preparing suspensions of microorganisms.	5 tubes	80405
McFarland 0.5 Barium Sulphate Standard	1 vial	80400
McFarland 1 Barium Sulphate Standard	1 vial	80401
McFarland 2 Barium Sulphate Standard	1 vial	80402
McFarland 3 Barium Sulphate Standard	1 vial	80403
McFarland 4 Barium Sulphate Standard	1 vial	80404
Physiological Solution Sterile physiological solution for diagnostic use.	20x7 mL	20095
VP (NaOH) Reagent	10x 10 mL	80280
VP (NaOH) Reagent Droppers	20 x 2.5 mL	87002
Vaseline Oil Supplementary reagent for identification systems.	4x50 mL 10 x 10 mL	80279 87006

Microbial staining and fixatives		
Description	Packaging	Ref.
ALBERT COLOR KIT Complete kit for <i>Corynebacterium diphtheriae</i> staining according to Albert.	2x250 mL	80294
Crystal Violet Solution (Gram staining) Crystal violet solution for Gram staining.	1000 mL	80299
DECOLOURIZING Solution (Gram staining) Decolourizing solution for Gram staining.	1000 mL	80295
GRAM COLOR KIT Complete kit for bacteria Gram staining.	4x250 mL 4x2x12 mL	80293 87101
KINYOUN COLOR KIT Complete kit for cold staining alcohol-acid resistant microorganisms.	2x250 mL	80380
LUGOL PVP Solution (Gram staining) Lugol solution for Gram staining.	1000 mL	80296
METHYLENE BLUE 0.3% Methylene blue solution for Ziehl-Neelsen staining.	250 mL	80277
MIF COLOR KIT Solution for eggs and fecal parasites detection.	50x12 mL	80275
SAFRANINE Solution (Gram staining) Safranine solution for Gram staining.	1000 mL	80290
ZIEHL – NEELSEN COLOR KIT Complete kit for hot staining alcohol-acid resistant microorganisms.	3x250 mL	80276

Rapid test for Sugar Fermentation (dried medium in tube)		
Description	Packaging	Ref.
Adonitol Test Rapid test for adonitol fermentation.	30/pack	88208
Arabinose Test Rapid test for arabinose fermentation.	30/pack	88209
Arabitol Test Rapid test for arabitol fermentation.	30/pack	88207
Dulcitol Test Rapid test for dulcitol fermentation.	30/pack	88210
Galactose Test Rapid test for galactose fermentation.	30/pack	88201
Glucose Test Rapid test for glucose fermentation.	30/pack	88202
Inositol Test Rapid test for inositol fermentation.	30/pack	88211
Inulin Test Rapid test for inulin fermentation.	30/pack	88212
Lactose Test Rapid test for lactose fermentation.	30/pack	88203
Levulose Test Rapid test for levulose fermentation.	30/pack	88213
Maltose Test Rapid test for maltose fermentation.	30/pack	88204
Mannitol Test Rapid test for mannitol fermentation.	30/pack	88214
Mannose Test Rapid test for mannose fermentation.	30/pack	88215
Raffinose Test Rapid test for raffinose fermentation.	30/pack	88205
Rhamnose Test Rapid test for rhamnose fermentation.	30/pack	88216
Salicin Test Rapid test for salicin fermentation.	30/pack	88217
Sorbitol Test Rapid test for sorbitol fermentation.	30/pack	88218
Sucrose Test Rapid test for sucrose fermentation.	30/pack	88206
Trehalose Test Rapid test for trehalose fermentation.	30/pack	88219
Xylose Test Rapid test for xylose fermentation.	30 tubes 30 discs	88220

Tests for microbial identification		
Description	Packaging	Ref.
Aesculin Bile Test Rapid test for aesculin hydrolysis.	30 tubes	88008
Arginine Decarboxylase Test Rapid test for arginine decarboxylation.	30 tubes	88016
Bacitracin Test Bacitracin 0.04 IU impregnated diagnostic discs for group A streptococci identification.	100 discs	9502
Beta Lactamase Sticks Sticks impregnated with Nitrocefin, used for the detection of β-Lactamase producing microorganisms.	30 Sticks	88033
CAMP Test – R Rapid test for listeriae identification.	30 discs	88027
CAMP Test – S Rapid test for listeriae identification.	30 discs	88021
Catalase Reagent Solution of hydrogen peroxide for microbial identification through the rapid detection of catalase enzyme.	20 x 1 mL dropper	87003
Catalase / Oxy Test Rapid test for catalase and peroxidase enzyme determination.	30 tubes	88023
H₂O₂ Reagent Reagent for catalase test.	5 mL	80057
Citrate Disc for the identification of Enterobacteriaceae and non-fermenting gram-negative bacteria.	30 test	88042
Coagulase Test Citrated freeze-dried rabbit plasma for coagulase test in tube.	5 vials (40 Tests)	88030
C-390 Discs containing antibiotic C-390 for the identification of <i>Pseudomonas aeruginosa</i> .	50 discs	88040
E. coli Test Rapid test for <i>Escherichia coli</i> identification (MUG, ONPG, Indole).	30 tubes 3 mL Kovac's	88006
Gram Test Stick Rapid test for gram-negative and gram-positive germs identification.	30 Sticks 3 mL DMAC	88031
H₂S Rapid Test Rapid test for H ₂ S production.	30 tubes	88013
Hippurate Test Rapid test for hippurate hydrolysis.	30 tubes	88007
Indole Test Stick Rapid test for indole production.	30 Sticks 3 mL DMAC	88032
Indole Test Rapid test for indole production.	30 tubes 3 mL Kovac's	88017
Lactophenol Cotton Blue Droppers Blue stain for microscopic examination of fungi, in dropper bottles.	10 bottles	87008
Listeria-mono Test Rapid test for <i>Listeria monocytogenes</i> identification.	20 tubes 20 discs 3 mL reagent	88010
Lysine Decarboxylase Test Rapid test for lysine decarboxylation.	30 tubes	88014
Metronidazole Test Discs for <i>Gardnerella vaginalis</i> identification.	100 discs	9508
Ninhydrin 7% Reagent for Hippurate Test.	10 mL bottle	80273
Nitrate Test Rapid test for nitrates to nitrites reduction.	30 tubes 3 mL α-Naphthylamine 3 mL Sulphanilic acid	88009
O129 discs Test for differentiating <i>Vibro</i> spp. from other Gram-negative bacteria.	50 discs 150 µg 50 discs 10 µg	88043 88044
ONPG Discs Rapid test for detecting the beta-galactosidase activity.	5 x 50 discs	88105
ONPG Test Rapid test for Lactose fermenting microorganisms differentiation.	30 tubes	88005
Optochine Test Diagnostic discs for pneumococci identification.	100 discs	9501
Ornithine Decarboxylase Test Rapid test for ornithine decarboxylation.	30 tubes	88015
Oxidase Test Discs Rapid test for cytochrome-oxidase enzyme detection.	30 discs	88004
Oxidase Test Sticks Rapid test for cytochrome-oxidase enzyme detection.	50 Sticks	88029
Oxidase Test Sticks Rapid test for detection of cytochrome oxidase enzymatic activity from pure colony.	50 Sticks	88029N

Tests for microbial identification		
Description	Packaging	Ref.
Oxidase Test Swabs Rapid test for cytochrome-oxidase enzyme detection.	30 swabs	88003
Peptidase A Sticks Rapid test for group A and group D streptococci identification.	30 Sticks 5 mL Pyr reagent	88034
Phenylalanine Reagent Droppers Solution of ferric chloride for microbial identification through the detection of phenylpyruvic acid.	20 x 2.5 mL	87004
S.F. Rapid Test Rapid test for group D streptococci determination.	30 tubes	88020
Sulphonamide Diagnostic discs for <i>Gardnerella vaginalis</i> identification.	100 discs	9511
Urea / Indole Test Rapid test for urea and indole detection.	30 tubes 3 mL Kovac's	88024
Urea Rapid Test Rapid test for urease determination.	30 tubes	88011
V Factor Test Diagnostic discs for <i>Haemophilus</i> spp. identification.	100 discs	9504
X Factor Test Diagnostic discs for <i>Haemophilus</i> spp. identification.	100 discs	9503
X + V Factor Test Diagnostic discs for <i>Haemophilus</i> spp. identification.	100 discs	9505

Antibiotic discs in cartridge					
Description		µg	CLSI ¹	EUCAST ^{3,4}	Ref.*
Amikacin	AK	30	✓	✓	9004
Amoxicillin	AML	2			9151
Amoxicillin	AML	10		✓	9133
Amoxicillin	AML	25			9179
Amoxicillin	AML	30			9005
Amoxicillin-clavulanic acid	AUG	3 (2/1)		✓	9191
Amoxicillin-clavulanic acid	AUG	7.5			9255
Amoxicillin-clavulanic acid	AUG	30 (20/10)	✓	✓	9048
Amoxicillin 10 + Clavulanic acid 0.1	AC	10.1 (10/0.1)			9278 ♦
Amoxicillin10 + Clavulanic acid 0.5	AC	10.5 (10/0.5)			9279 ♦
Amoxicillin 10 + Clavulanic acid 1	AC	11 (10/1)			9280 ♦
Ampicillin	AMP	2		✓	9115
Ampicillin	AMP	10	✓	✓	9006
Ampicillin-sulbactam	AMS	20 (10/10)	✓	✓	9031
Ampliclox (Ampicillin + Cloxacillin)	ACL	30 (25/5)			9122
Azithromycin	AZM	15	✓		9105
Azlocillin	AZL	75	✓		9007
Aztreonam	ATM	30	✓	✓	9008
Bacitracin	BA	10 units			9051
Carbenicillin	CAR	100	✓		9009
Cefaclor	CEC	30	✓	✓	9010
Cefadroxil	CDX	30		✓	9052
Cefamandole	MA	30	✓		9014
Cefazolin	KZ	30	✓		9015
Cefepime	FEP	5			9219
Cefepime	FEP	10			9220
Cefepime	FEP	30	✓	✓	9104
Cefepime + Clavulanic acid	FEL	40 (30/10)		✓ ⁹	9143
Cefiderocol	FDC	15			9265
Cefiderocol	FDC	30	✓		9266
Cefixime	CFM	5	✓	✓	9089
Cefoperazone	CFP	10			9210 ♦
Cefoperazone	CFP	30			9016
Cefoperazone	CFP	75	✓		9108
Cefotaxime	CTX	5		✓	9152
Cefotaxime	CTX	30	✓		9017
Cefotaxime + Clavulanic acid	CTL	15 (5/10)		✓ ⁸	9257 ♦
Cefotaxime + Clavulanic acid	CTL	40 (30/10)		✓ ⁹	9182
Cefotaxime + Clavulanic acid + Cloxacillin	CTLC			✓ ⁹	9203
Cefotaxime + Cloxacillin	CTC	230 (30/200)		✓ ⁹	9224
Cefotetan	CTT	30	✓		9081
Cefoxitin	FOX	30	✓	✓	9018
Cefoxitin + Cloxacillin	FOC	230 (30/200)			9144
Cefpirome	CR	30			9185
Cefpodoxime	PX	10	✓	✓	9064
Cefpodoxime + Clavulanic acid	PXL	11 (10/1)			9190
Cefprozil	CPR	30	✓		9112
Cefsulodin	CSD	30			9053
Ceftaroline	CPT	5		✓	9195
Ceftaroline	CPT	30	✓		9198
Ceftazidime	CAZ	10		✓	9153
Ceftazidime	CAZ	30	✓		9019
Ceftazidime-avibactam	CZA	14 (10/4)		✓	9206
Ceftazidime-avibactam	CZA	50 (30/20)	✓		9205
Ceftazidime + Clavulanic acid	CAL	20 (20/10)		✓ ⁸	9258 ♦
Ceftazidime + Clavulanic acid	CAL	40 (30/10)		✓ ⁹	9145
Ceftazidime + Clavulanic acid + Cloxacillin	CALC			✓ ⁹	9204
Ceftazidime + Cloxacillin	CAC			✓ ⁹	9225
Ceftibuten	CTB	30	✓	✓	9101

Antibiotic discs in cartridge						
Description		µg		CLSI ¹	EUCAST ^{3,4}	Ref.*
Ceftizoxime	CZX	30		✓		9054
Ceftobiprole	BPR	5			✓	9242
Ceftolozane-tazobactam	C/T	40	(30/10)	✓	✓	9246
Ceftriaxone	CRO	5				9235 ♦
Ceftriaxone	CRO	30		✓	✓	9020
Cefuroxime	CXM	1				9232
Cefuroxime	CXM	5				9236 ♦
Cefuroxime	CXM	30		✓	✓	9021
Cephalexin	CL	30			✓	9011
Cephalothin	KF	30		✓		9013
Cephadrine	CE	30				9055
Chloramphenicol	C	10				9128
Chloramphenicol	C	30		✓	✓	9022
Cinoxacin	CIN	100		✓		9057
Ciprofloxacin	CIP	5		✓	✓	9056
Clarithromycin	CLR	15		✓		9098
Clavulanic acid	CLA	1				9229 ♦
Clavulanic acid	CLA	2				9228 ♦
Clavulanic acid	CLA	5				9230 ♦
Clavulanic acid	CLA	10				9231 ♦
Clindamycin	CD	2		✓	✓	9047
Clindamycin	CD	10				9146
Cloxacillin	CX	5				9058
Colistin sulfate	CS	10		✓		9023
Colistin sulfate	CS	25				9184
Colistin sulfate	CS	30	units			9141
Daptomycin (includes Ca ²⁺)	DAP	30				9090
Dicloxacillin	DCX	1				9093
Dipicolinic acid	DP					9194
Doripenem	DOR	10		✓		9154
Doxycycline	DXT	30		✓		9059
EDTA	ED					9087
Eravacycline	ERV	20		✓		9238 ♦
Eravacycline	ERV	50				9237 ♦
Ertapenem	ETP	10		✓	✓	9061
Ertapenem + Cloxacillin	ET+CL					9199
Ertapenem + Phenylboronic acid	ET+BO					9202
Erythromycin	E	2				9180
Erythromycin	E	15		✓	✓	9024
Fosfomicin	FOS	50				9025
Fosfomicin (includes G6P 50)	FOS	100	(100/50)			9121
Fosfomicin (includes G6P 50)	FOS	200	(200/50)	✓	✓	9109
Fosfomicin 200 + G6P 100	FGP	100	(200/100)			9213 ♦
Fosfomicin 200 + G6P 200	FGP	200	(200/200)			9214 ♦
Fosfomicin 200 + G6P 400	FGP	400	(200/400)			9215 ♦
Furazolidon	FR	50				9099
Fusidic acid	FC	10		✓	✓	9049
Fusidic acid	FC	30				9111
Gatifloxacin	GAT	5		✓		9169
Gentamicin	CN	10		✓	✓	9026
Gentamicin	CN	30			✓	9125
Gentamicin	CN	120		✓		9124
Glucose-6-Phosphate	G6P	100				9216 ♦
Glucose-6-Phosphate	G6P	200				9217 ♦
Glucose-6-Phosphate	G6P	400				9218 ♦
Imipenem	IMI	10		✓	✓	9079
Imipenem + Cloxacillin	IMI+CL					9086
Imipenem + EDTA	IM+ED					9183
Imipenem + Phenylboronic acid	IMI+BO					9085

Antibiotic discs in cartridge					
Description		µg	CLSI ¹	EUCAST ^{3,4}	Ref.*
Kanamycin	K	30	✓		9027
Lefamulin	LMU	5			9249 ♦
Lefamulin	LMU	20	✓		9250 ♦
Levofloxacin	LEV	5	✓	✓	9102
Lincomycin	MY	2			9028
Lincomycin	MY	15			9116
Linezolid	LNZ	10		✓	9155
Linezolid	LNZ	30	✓		9136
Lomefloxacin	LOM	10	✓		9113
Mecillinam	MEC	10	✓	✓	9156
Meropenem	MRP	10	✓	✓	9068
Meropenem + Cloxacillin	MR+CL			✓ ⁹	9175
Meropenem + Dipicolinic acid	MR+DP			✓ ⁹	9177
Meropenem + EDTA	MR+ED			✓ ⁹	9178
Meropenem + Phenylboronic acid	MR+BO			✓ ⁹	9176
Methicillin	MET	5			9029
Metronidazole	MTZ	5			9076
Metronidazole	MTZ	10			9188
Metronidazole	MTZ	50			9119
Mezlocillin	MEZ	75			9062
Minocycline	MN	30	✓	✓	9030
Moxifloxacin	MXF	5	✓	✓	9103
Mupirocin	MUP	5			9189
Mupirocin	MUP	200	✓	✓	9157
Nafcillin	NAF	1	✓		9174
Nalidixic acid	NA	30	✓	✓	9001
Netilmicin	NET	10		✓	9170
Netilmicin	NET	30	✓		9033
Nitrofurantoin	F	50			9181
Nitrofurantoin	F	100		✓	9158
Nitrofurantoin	F	300	✓		9034
Nitroxoline	NI	30		✓	9209
Norfloxacin	NOR	10		✓	9035
Novobiocin	NO	5			9117
Novobiocin	NO	30			9063
Ofloxacin	OFX	5	✓	✓	9080
Oritavancin	ORI	5			9201
Oxacillin	OX	1	✓	✓	9036
Oxacillin	OX	5			9135
Oxolinic acid	OA	2			9002
Oxytetracycline	OT	30			9065
Pefloxacin	PEF	5	✓	✓	9091
Penicillin G	P	1 unit		✓	9130
Penicillin G	P	2 units			9127
Penicillin G	P	10 units	✓		9037
Penicillin G 1 + Clavulanic acid 1	PC	1 (1 unit/1)			9226 ♦
Penicillin G 1 + Clavulanic acid 2	PC	2 (1 unit/2)			9227 ♦
Penicillin G 1 + Clavulanic acid 3	PC	3 (1 unit/3)			9239 ♦
Penicillin G 1 + Clavulanic acid 5	PC	5 (1 unit/5)			9240 ♦
Penicillin G 1 + Clavulanic acid 10	PC	10 (1 unit/10)			9241 ♦
Penicillin G 1 + Clavulanic acid 20	PC	20 (1 unit/20)			9254 ♦
Phenoxymethylpenicillin	PV	10			9171
Phenylboronic acid	BO				9193
Pipemidic acid	PI	20			9003
Piperacillin	PRL	30		✓	9159
Piperacillin	PRL	100	✓		9038
Piperacillin-tazobactam	TZP	36 (30/6)		✓	9160
Piperacillin-tazobactam	TZP	110 (100/10)	✓		9100
Polymyxin B	PB	100 units			9066

Antibiotic discs in cartridge						
Description		µg		CLSI ¹	EUCAST ^{3,4}	Ref.*
Polymyxin B	PB	300	units	✓		9120
Rifampicin	RD	5		✓	✓	9118
Rifampicin	RD	30				9039
Roxithromycin	RXT	15				9060
Sisomicin	SIS	30				9046
Sodium Fusidate	FC	30				9131
Spectinomycin	SPC	100		✓		9067
Streptomycin	S	10		✓		9040
Streptomycin	S	300		✓	✓	9162
Sulbactam	SU	6				9212 ♦
Sulbactam	SU	10				9211 ♦
Sulbactam	SU	20				9129
Sulfadiazine	SUZ	300				9150
Sulfafurazole	SF	300				9041
Sulfamethoxazole	SMX	50				9084
Sulfamethoxazole	SMX	100				9187
Sulfaprim	SXT	50				9132
Sulfonamide	S3	300		✓		9126
Tedizolid	TZD	2			✓	9243
Teicoplanin	TEC	30		✓	✓	9050
Telithromycin	TEL	15		✓	✓	9172
Temocillin	TMO	30			✓ ⁹	9186
Tetracycline	TE	30		✓	✓	9043
Ticarcillin	TC	75		✓	✓	9070
Ticarcillin-clavulanic acid	TTC	85	(75/10)	✓	✓	9096
Tigecycline	TGC	15		✓	✓	9147
Tobramycin	TOB	10		✓	✓	9044
Tobramycin	TOB	30				9163
Trimethoprim	TM	2.5				9083
Trimethoprim	TM	5		✓	✓	9110
Trimethoprim-sulfamethoxazole	SXT	25	(1.25/23.75)	✓	✓	9042
Ulifloxacin	ULI	5		✓		9207
Vancomycin	VA	5			✓	9164
Vancomycin	VA	30		✓		9045

●, to be stored at -20°C ♦, not CE marked

Storage Unless otherwise indicated, discs may be stored at any temperature between -20°C and +8°C.

CE Marking Unless otherwise indicated, all products intended for clinical applications are CE marked according to the European Directive 98/79/EC for In Vitro Diagnostic Medical Devices.

*** Packaging** Each pack contains 5 cartridges of 50 discs (5x50 discs = 250 discs).
Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.
Example: ref. 9045/1 indicates Vancomycin 30 µg in one single cartridge of 50 discs.

Note that testing and reporting antimicrobial agents for which there are no interpretive criteria are the responsibility of the chief microbiologist and such decisions should be made with input from the infectious disease clinicians.

For full details on specific organism/agent combinations refer to current CLSI and EUCAST recommendations.
Also see the EUCAST guidance at http://www.eucast.org/clinical_breakpoints/when_there_are_no_breakpoints/

Antibiotic discs in cartridge for Veterinary use				
Description		µg	CLSI ⁶	Ref.*
Ceftiofur	FUR	30	✓	9251
Enrofloxacin	ENR	5	✓	9233
Florfenicol	FFC	30	✓	9234
Neomycin	N	30	✓	9032
Spiramycin	SP	100		9088
Tiamulin	T	30	✓	9094
Tylosin	TY	30	✓	9082

- * **Packaging** Each pack contains 5 cartridges of 50 discs (5x50 discs = 250 discs).
Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.
Example: ref. 9082/1 indicates Tylosin 30 µg in one single cartridge of 50 discs.

CE Marking All antifungal discs are CE marked according to the European Directive 98/79/EC for *In Vitro* Diagnostic Medical Devices.

1. CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 29th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2019.
2. CLSI. Performance Standards for Antimicrobial Disk Susceptibility Tests; 13th ed. CLSI standard M02. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
3. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 9.0, 2019. <http://www.eucast.org>
4. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST, Version 9.0, 2019. <http://www.eucast.org>
5. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing - Version 7.0 (January 2019).
6. CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals, 4th ed. CLSI supplement VET08. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
7. CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals; 5th ed. CLSI standard VET01. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
8. EUCAST technical guidance on the use of the combination disk test (CDT) for confirmation of ESBL in Enterobacterales. New disk potencies for combination disks containing cefotaxime and ceftazidime without and with clavulanic acid (February 2019).
9. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0 (July 2017).
10. DIN 58940-2 Medical microbiology - Susceptibility testing of microbial pathogens to antimicrobial agents - Part 2: Active substance carriers for the agar diffusion test; 2007-10.
11. FDA (1978) Codes of Fed.Rebs. 21.Part 460.
12. WHO (1977) Tech rep.Ser.n°610.

Antimycotic discs in cartridge




Description		µg	CLSI ¹³	Ref.*
Amphotericin B	AMB	10		9137
Amphotericin B	AMB	20		9071
Caspofungin	CAS	5	✓	9165
Clotrimazole	CLO	50		9097
Econazole	ECN	10		9072
Fluconazole	FLU	25	✓	9166
Fluconazole	FLU	100		9069
Flucytosine	AFY	1		9073
Flucytosine	AFY	10		9148
Griseofulvin	AGF	10		9074
Itraconazole	ITC	8		9139
Itraconazole	ITC	50		9107
Ketoconazole	KCA	10		9075
Ketoconazole	KCA	15		9140
Miconazole	MCL	10		9077
Nystatin	NY	100 IU		9078
Posaconazole	POS	5	✓	9167
Voriconazole	VO	1	✓	9168

* **Packaging** Each pack contains 5 cartridges of 50 discs (5x50 discs = 250 discs).
Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.
Example: ref. 9075/1 indicates Ketoconazole 10 µg in one single cartridge of 50 discs.


¹³ CLSI. Performance Standards for Antifungal Susceptibility Testing of Yeasts; 1st ed. CLSI Supplement M60. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.

Multodisc with eight antimicrobial agents



Description		µg	CLSI	EUCAST	Packaging	Ref.
ANAEROBES	Carbenicillin	CAR	100		100/pack	95200
	Cefoxitin	FOX	30			
	Clindamycin	CD	2			
	Chloramphenicol	C	30			
	Erythromycin	E	15			
	Metronidazole	MTZ	5			
	Ampicillin	AMP	10			
	Tetracycline	TE	30			
ENTEROCOCCI	Tigecycline	TGC	15		100/pack	95210
	Ampicillin	AMP	2			
	Linezolid	LNZ	10			
	Nitrofurantoin	F	100			
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Vancomycin	VA	5			
	Gentamicin	CN	30			
	Teicoplanin	TEC	30	✓		
ENTEROBACTERIA URINE	Levofloxacin	LEV	5	✓	100/pack	95230
	Ampicillin	AMP	10			
	Cefuroxime	CXM	30	✓		
	Ceftazidime	CAZ	10			
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Gentamicin	CN	10	✓		
	Nitrofurantoin	F	100	✓		
	Norfloxacin	NOR	10	✓		
ENTEROBACTERIA 1	Amikacin	AK	30	✓	100/pack	95220
	Aztreonam	ATM	30	✓		
	Piperacillin + Tazobactam	TZP	36			
	Chloramphenicol	C	30	✓		
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Netilmicin	NET	10			
	Cefotaxime	CTX	5			
	Nitrofurantoin	F	100			
ENTEROBACTERIA 2	Levofloxacin	LEV	5	✓	100/pack	95240
	Cefoxitin	FOX	30	✓		
	Piperacillin + Tazobactam	TZP	36			
	Cefuroxime	CXM	30	✓		
	Streptomycin	S	300			
	Netilmicin	NET	10			
	Ampicillin + Sulbactam	AMS	20	✓		
	Tobramycin	TOB	10	✓		
PSEUDOMONAS	Amikacin	AK	30	✓	100/pack	95250
	Ticarcillin	TC	75			
	Piperacillin + Tazobactam	TZP	36			
	Ceftazidime	CAZ	10			
	Gentamicin	CN	10	✓		
	Netilmicin	NET	10	✓		
	Levofloxacin	LEV	5	✓		
	Tobramycin	TOB	10	✓		

Multodisc with eight antimicrobial agents



Description		µg	CLSI	EUCAST	Packaging	Ref.
STAPH	Cefoxitin	FOX	30	✓	✓	100/pack 95260
	Levofloxacin	LEV	5	✓	✓	
	Netilmicin	NET	10	✓	✓	
	Tetracycline	TE	30	✓	✓	
	Erythromycin	ERY	15	✓	✓	
	Gentamicin	CN	10	✓	✓	
	Rifampicin	RD	5	✓	✓	
	Linezolid	LNZ	10	✓	✓	
STREPTO	Penicillin G	P	10 IU	✓	✓	100/pack 95290
	Levofloxacin	LEV	5	✓	✓	
	Teicoplanin	TEC	30	✓	✓	
	Vancomycin	VA	5	✓	✓	
	Erythromycin	E	15	✓	✓	
	Tetracycline	TE	30	✓	✓	
	Linezolid	LNZ	10	✓	✓	
Rifampicin	RD	5	✓	✓		
ACINETOBACTER	Doripenem	DOR	10	✓	✓	100/pack 95270
	Meropenem	MRP	10	✓	✓	
	Levofloxacin	LEV	5	✓	✓	
	Gentamicin	CN	10	✓	✓	
	Trimethoprim-Sulfamethoxazole	SXT	25	✓	✓	
	Tobramycin	TOB	10	✓	✓	
	Ciprofloxacin	CIP	5	✓	✓	
Netilmicin	NET	10	✓	✓		
YEASTS	Caspofungin	CAS	5	✓	✓	100/pack 95280
	Fluconazole	FLU	25	✓	✓	
	Posaconazole	POS	5	✓	✓	
	Voriconazole	VO	1	✓	✓	
	Amphotericin B	AMB	20	✓	✓	
	Ketoconazole	KCA	10	✓	✓	
	Flucytosine	AFY	1	✓	✓	
	Nystatin	NY	100 IU	✓	✓	

EUCAST ESBL disk kit			Packaging	Ref.
Description		µg		
Cefotaxime	CTX	5	50 Test	99002
Cefotaxime + Clavulanic acid	CTL	15 (5+10)		
Ceftazidime	CAZ	10		
Ceftazidime + Clavulanic acid	CAL	20 (10+10)		
Cefepime	FEP	30		
Cefepime + Clavulanic acid	FEL	40 (30+10)		


ESBL discs kit (according to EUCAST)			Packaging	Ref.
Description		µg		
Cefotaxime	CTX	30	50 Test	99004
Cefotaxime + Clavulanic acid	CTL	40 (30+10)		
Ceftazidime	CAZ	30		
Ceftazidime + Clavulanic acid	CAL	40 (30+10)		
Cefepime	FEP	30		
Cefepime + Clavulanic acid	FEL	40 (30+10)		

ESBL discs kit (according to CLSI)			Packaging	Ref.
Description		µg		
Cefotaxime	CTX	30	50 Test	99005
Cefotaxime + Clavulanic acid	CTL	40 (30+10)		
Ceftazidime	CAZ	30		
Ceftazidime + Clavulanic acid	CAL	40 (30+10)		

ESBL (chromosomal ind. AmpC) discs kit (according to EUCAST)			Packaging	Ref.
Description		µg		
Cefepime	FEP	30	100 Test	99006
Cefepime + Clavulanic acid	FEL	40 (30+10)		

KPC&MBL discs kit (according to EUCAST)			Packaging	Ref.
Description		µg		
Meropenem	MRP	10	50 Test	99003
Meropenem + Phenylboronic acid	MR+BO			
Meropenem + Cloxacillin	MR+CL			
Meropenem + EDTA	MR+EDTA			

KPC&MBL&OXA-48 discs kit			Packaging	Ref.
Description		µg		
Meropenem	MRP	10	50 Test	99007
Meropenem + Phenylboronic acid	MR+BO			
Meropenem + Cloxacillin	MR+CL			
Meropenem + EDTA	MR+EDTA			
Temocillin	TMO	30		



AmpC discs kit			
Description	µg	Packaging	Ref.
Cefotaxime	30	50 Test	99009
Cefotaxime + Cloxacillin			
Ceftazidime	30		
Ceftazidime + Cloxacillin			

ESBL+AmpC screen discs kit			
Description	µg	Packaging	Ref.
Cefotaxime	30	50 Test	99008
Cefotaxime + Clavulanic acid	40 (30+10)		
Cefotaxime + Cloxacillin			
Cefotaxime + Clavulanic acid + Cloxacillin			

Disc Dispenser		Ref.
Description		
		
Disc Dispenser for 6 cartridges		91203
Disc Dispenser for 8 cartridges		91200


Densitometer		Ref.
Description		
		
DEN-1B McFarland densitometer with built-in adapter for 18 mm tubes		96881
A-16 Adapter for 16 mm tubes		96882
CKG18 , calibration kit for DEN-1, for 18mm glass tubes		96883
CKG16 , calibration kit for DEN-1, for 16mm glass tubes		96884

Rotating Plate Inoculator		Ref.
Description		
AUTOMATIC ROTATION PLATER PRO for confluent microbial growth on the 90 mm petri dish		96796
		

Agar Dilution panels


according to CLSI and EUCAST standards for determining MIC through agar dilution.

Patent pending



Description	packaging	ref.
AD Fosfomycin 0.25-256 Comprises 12 wells, containing the antibiotic incorporated into an agar medium in different concentrations, i.e. 11 two-fold dilutions (0.25-256 µg/mL)	1 panel 6 panels	77001 77061
Blood AD Clindamycin 0.03-32 A 12-well panel containing Clindamycin incorporated into a blood agar medium in different concentrations, i.e. 11 two-fold dilutions (0.03-32 µg/mL)	1 panel 6 panels	77002 77062
Customized versions 12-well panel with a customized set of antibiotic concentrations.	minimum order 500 panels	to be defined

NP range for antimicrobial resistance testing



Description	packaging	ref.
Rapid ESBL NP test® A 8-test panel used for the identification of ESBL producers among Enterobacteriaceae. The dried up substrata in the wells are rehydrated with a lysed microbial suspension. After 20 minutes incubation at 36±2°C the result is read and interpreted by a clear color change.	2 panel 8 tests/panel	76036
RapidResa Polymyxin Acinetobacter NP test® A 8-test panel for Rapid Detection of Colistin Resistance in <i>Acinetobacter baumannii</i> .	2 panel 8 tests/panel	76046

ATCC® derived CultiControl
freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices

Description	WDCM	CultiControl Ref.	notes
<i>Acinetobacter baumannii</i>		ATCC® 19606 ^{TM*} 89174	
<i>Acinetobacter baumannii</i>		ATCC® BAA-747 ^{TM*} 89141	
<i>Actinomyces odontolyticus</i>		ATCC® 17929 ^{TM*} 89114	
<i>Aeromonas hydrophila</i>		ATCC® 35654 ^{TM*} 89169	
<i>Aeromonas hydrophila</i>	63	ATCC® 7966 ^{TM*} 89119	
<i>Aggregatibacter aphrophilus</i>		ATCC® 7901 ^{TM*} 89091	
<i>Aspergillus brasiliensis</i>	53	ATCC® 16404 ^{TM*} 89021	
<i>Bacillus cereus</i>		ATCC® 10876 ^{TM*} 89155	
<i>Bacillus cereus</i>	1	ATCC® 11778 ^{TM*} 89022	
<i>Bacillus subtilis subsp. spizizenii</i>	3	ATCC® 6633 ^{TM*} 89023	
<i>Bacteroides fragilis</i>		ATCC® 23745 ^{TM*} 89113	
<i>Bacteroides fragilis</i>		ATCC® 25285 ^{TM*} 89078	
<i>Bacteroides ovatus</i>		ATCC® 8483 ^{TM*} 89111	
<i>Bacteroides ovatus</i>		ATCC® BAA-1296 ^{TM*} 89193	
<i>Bacteroides thetaiotaomicron</i>		ATCC® 29741 ^{TM*} 89079	
<i>Bordetella bronchiseptica</i>		ATCC® 4617 ^{TM*} 89139	
<i>Burkholderia cepacia</i>		ATCC® 25416 ^{TM*} 89147	
<i>Burkholderia cepacia</i>		ATCC® 25608 ^{TM*} 89166	
<i>Campylobacter jejuni subsp. jejuni</i>	156	ATCC® 29428 ^{TM*} 89167	
<i>Campylobacter jejuni subsp. jejuni</i>	5	ATCC® 33291 ^{TM*} 89086	
<i>Campylobacter jejuni subsp. jejuni</i>		ATCC® 33560 ^{TM*} 89145	
<i>Candida albicans</i>	54	ATCC® 10231 ^{TM*} 89024	
<i>Candida albicans</i>		ATCC® 14053 ^{TM*} 89183	
<i>Candida albicans</i>		ATCC® 18804 ^{TM*} 89177	
<i>Candida albicans</i>		ATCC® 64124 ^{TM*} 89178	
<i>Candida albicans</i>		ATCC® 90028 ^{TM*} 89072	
<i>Candida krusei</i>		ATCC® 14243 ^{TM*} 89098	
<i>Candida parapsilosis</i>		ATCC® 22019 ^{TM*} 89071	
<i>Candida tropicalis</i>		ATCC® 750 ^{TM*} 89097	
<i>Citrobacter freundii</i>		ATCC® 43864 ^{TM*} 89146	
<i>Citrobacter freundii</i>		ATCC® 8090 ^{TM*} 89159	
<i>Clostridium difficile</i>		ATCC® 9689 ^{TM*} 89090	produces cytotoxin
<i>Clostridium histolyticum</i>		ATCC® 19401 ^{TM*} 89112	
<i>Clostridium perfringens</i>	7	ATCC® 13124 ^{TM*} 89053	
<i>Clostridium sordellii</i>		ATCC® 9714 ^{TM*} 89059	
<i>Clostridium sporogenes</i>	8	ATCC® 19404 ^{TM*} 89095	
<i>Cronobacter mutijensii</i>		ATCC® 51329 ^{TM*} 89158	
<i>Cronobacter sakazakii</i>	214	ATCC® 29544 ^{TM*} 89138	formerly Enterobacter sakazakii
<i>Eikenella corrodens</i>		ATCC® BAA-1152 ^{TM*} 89196	
<i>Enterobacter aerogenes</i>	175	ATCC® 13048 ^{TM*} 89156	
<i>Enterobacter cloacae subsp. cloacae</i>		ATCC® 49141 ^{TM*} 89200	
<i>Enterobacter cloacae subsp. cloacae</i>			control strain for the AmpC disk test; strong positive
<i>Enterococcus casseliflavus</i>		ATCC® BAA-1143 ^{TM*} 89065	
<i>Enterococcus faecalis</i>	9	ATCC® 700327 ^{TM*} 89195	
<i>Enterococcus faecalis</i>	87	ATCC® 19433 ^{TM*} 89025	
<i>Enterococcus faecalis</i>	210	ATCC® 29212 ^{TM*} 89026	
<i>Enterococcus faecalis</i>		ATCC® 33186 ^{TM*} 89115	
<i>Enterococcus faecalis</i>		ATCC® 49532 ^{TM*} 89066	high level Gentamicin-resistant and Streptomycin-sensitive
<i>Enterococcus faecalis</i>		ATCC® 49533 ^{TM*} 89067	high level Gentamicin-sensitive and Streptomycin-resistant
<i>Enterococcus faecalis</i>	85	ATCC® 51299 ^{TM*} 89173	Vancomycin resistant and high level aminoglycosides, vanB
<i>Enterococcus faecium</i>	10	ATCC® 19434 ^{TM*} 89171	
<i>Enterococcus faecium</i>		ATCC® 51559 ^{TM*} 89117	
<i>Enterococcus faecium</i>		ATCC® 6057 ^{TM*} 89152	
<i>Enterococcus faecium</i>		ATCC® BAA-2319 ^{TM*} 89172	vanA resistance
<i>Erysipelothrix rhusiopathiae</i>		ATCC® 19414 ^{TM*} 89187	
<i>Escherichia coli</i>		ATCC® 11303 ^{TM*} 89184	
<i>Escherichia coli</i>	13	ATCC® 25922 ^{TM*} 89027	
<i>Escherichia coli</i>		ATCC® 35218 ^{TM*} 89163	beta lactamase producer
<i>Escherichia coli</i>	12	ATCC® 8739 ^{TM*} 89028	
<i>Fluoribacter bozemanii</i>		ATCC® 33217 ^{TM*} 89157	

ATCC® derived CultiControl
freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices

Description	WDCM	CultiControl Ref.	notes
<i>Fusobacterium nucleatum subsp. nucleatum</i>		ATCC® 25586™*	89118
<i>Gardnerella vaginalis</i>		ATCC® 14018™*	89099
<i>Geobacillus stearothermophilus</i>		ATCC® 7953™*	89203
<i>Haemophilus haemolyticus</i>		ATCC® 33390™*	89123
<i>Haemophilus influenzae</i>		ATCC® 10211™*	89120
<i>Haemophilus influenzae</i>		ATCC® 19418™*	89160
<i>Haemophilus influenzae</i>		ATCC® 33391™*	89176
<i>Haemophilus influenzae</i>		ATCC® 33533™*	89124
<i>Haemophilus influenzae</i>		ATCC® 49247™*	89077
<i>Haemophilus influenzae</i>		ATCC® 49766™*	89076
<i>Haemophilus influenzae</i>		ATCC® 9007™*	89142
<i>Issatchenkia orientalis</i>		ATCC® 6258™*	89073
<i>Klebsiella pneumoniae</i>		ATCC® BAA-1144™*	89150
<i>Klebsiella pneumoniae</i>		ATCC® BAA-1705™*	89088
<i>Klebsiella pneumoniae</i>		ATCC® BAA-1706™*	89087
<i>Klebsiella pneumoniae</i>		ATCC® BAA-2146™*	89069
<i>Klebsiella pneumoniae subsp. pneumoniae</i>		ATCC® 13883™*	97 89089
<i>Klebsiella pneumoniae subsp. pneumoniae</i>		ATCC® 31488™*	89199
<i>Klebsiella pneumoniae subsp. pneumoniae</i>		ATCC® 4352™*	192 89192
<i>Klebsiella pneumoniae subsp. pneumoniae</i>		ATCC® 700603™*	89070
<i>Lactobacillus acidophilus</i>		ATCC® 4356™*	98 89080
<i>Lactobacillus fermentum</i>		ATCC® 9338™*	89100
<i>Lactobacillus leichmannii</i>		ATCC® 4797™*	89081
<i>Lactobacillus paracasei subsp. paracasei</i>		ATCC® BAA-52™*	89055
<i>Lactococcus lactis subsp. lactis</i>		ATCC® 19435™*	16 89082
<i>Legionella pneumophila subsp. fraseri</i>		ATCC® 33156™*	180 89151
<i>Legionella pneumophila subsp. pneumophila</i>		ATCC® 33152™*	107 89052
<i>Listeria grayi</i>		ATCC® 25401™*	89101
<i>Listeria innocua</i>		ATCC® 33090™*	17 89029
<i>Listeria ivanovii subsp. ivanovii</i>		ATCC® 19119™*	18 89030
<i>Listeria monocytogenes</i>		ATCC® 13932™*	21 89085
<i>Listeria monocytogenes</i>		ATCC® 15313™*	89188
<i>Listeria monocytogenes</i>		ATCC® 19111™*	20 89031
<i>Listeria monocytogenes</i>		ATCC® 19115™*	89051
<i>Listeria monocytogenes</i>		ATCC® 35152™*	109 89148
<i>Listeria monocytogenes</i>		ATCC® 7644™*	89060
<i>Listeria monocytogenes</i>		ATCC® BAA-751™*	89143
<i>Micrococcus luteus</i>		ATCC® 10240™*	89096
<i>Micrococcus luteus</i>		ATCC® 4698™*	111 89102
<i>Moraxella (Branhamella) catarrhalis</i>		ATCC® 25238™*	89103
<i>Neisseria gonorrhoeae</i>		ATCC® 19424™*	89074
<i>Neisseria gonorrhoeae</i>		ATCC® 31426™*	89075
<i>Neisseria gonorrhoeae</i>		ATCC® 49226™*	89104
<i>Neisseria gonorrhoeae</i>		ATCC® 49981™*	89122
<i>Neisseria meningitidis</i>		ATCC® 13090™*	89164
<i>Nocardia brasiliensis</i>		ATCC® 19296™*	89189
<i>Peptostreptococcus anaerobius</i>		ATCC® 27337™*	89165
<i>Plesiomonas shigelloides</i>		ATCC® 14029™*	89094
<i>Porphyromonas gingivalis</i>		ATCC® 33277™*	89162
<i>Prevotella melaninogenica</i>		ATCC® 25845™*	89134
<i>Propionibacterium acnes</i>		ATCC® 11827™*	89135
<i>Proteus hauseri</i>		ATCC® 13315™*	89190
<i>Proteus mirabilis</i>		ATCC® 12453™*	89049
<i>Proteus mirabilis</i>		ATCC® 25933™*	89032
<i>Proteus mirabilis</i>		ATCC® 29906™*	23 89083
<i>Proteus mirabilis</i>		ATCC® 35659™*	89105

ATCC® derived CultiControl
freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices

Description	WDCM	CultiControl Ref.	notes
<i>Proteus mirabilis</i>		ATCC® 43071™*	89106
<i>Proteus vulgaris</i>		ATCC® 6380™*	89107
<i>Providencia stuartii</i>		ATCC® 33672™*	89125
<i>Pseudomonas aeruginosa</i>	24	ATCC® 10145™*	89108
<i>Pseudomonas aeruginosa</i>		ATCC® 15442™*	89109
<i>Pseudomonas aeruginosa</i>	25	ATCC® 27853™*	89033
<i>Pseudomonas aeruginosa</i>	26	ATCC® 9027™*	89034
<i>Pseudomonas fluorescens</i>	115	ATCC® 13525™*	89110
			recommended for CAMP test for <i>Listeria</i> monocytogenes
<i>Rhodococcus equi</i>	28	ATCC® 6939™*	89035
<i>Saccharomyces cerevisiae</i>	58	ATCC® 9763™*	89036
<i>Salmonella enterica subsp. arizonae</i>		ATCC® 13314™*	89154
<i>Salmonella enterica subsp. enterica serovar Enteritidis</i>	30	ATCC® 13076™*	89084
<i>Salmonella enterica subsp. enterica serovar Hillingdon</i>		ATCC® 9184™*	89185
<i>Salmonella enterica subsp. enterica serovar Paratyphi</i>		ATCC® 9150™*	89161
<i>Salmonella enterica subsp. enterica serovar Typhimurium</i>	121	ATCC® 13311™*	89054
<i>Salmonella enterica subsp. enterica serovar Typhimurium</i>	31	ATCC® 14028™*	89037
<i>Salmonella enterica subsp. enterica serovar Typhimurium</i>		ATCC® 49416™*	89197
<i>Serratia marcescens</i>		ATCC® 14756™*	89191
<i>Serratia marcescens</i>		ATCC® 8100™*	89121
<i>Shigella boydii</i>		ATCC® 9207™*	89179
<i>Shigella flexneri</i>	126	ATCC® 12022™*	89038
<i>Shigella flexneri</i>		ATCC® 9199™*	89198
<i>Shigella sonnei</i>		ATCC® 25931™*	89058
<i>Shigella sonnei</i>		ATCC® 9290™*	89180
<i>Staphylococcus aureus</i>		ATCC® 33862™*	89042
<i>Staphylococcus aureus</i>	193	ATCC® 6538™*	89044
<i>Staphylococcus aureus subsp. aureus</i>		ATCC® 19095™*	89137
<i>Staphylococcus aureus subsp. aureus</i>	34	ATCC® 25923™*	89040
<i>Staphylococcus aureus subsp. aureus</i>	131	ATCC® 29213™*	89041
<i>Staphylococcus aureus subsp. aureus</i>		ATCC® 33591™*	89116
<i>Staphylococcus aureus subsp. aureus</i>	211	ATCC® 43300™*	89043
<i>Staphylococcus aureus subsp. aureus</i>		ATCC® 49476™*	89181
<i>Staphylococcus aureus subsp. aureus</i>		ATCC® 700698™*	89092
<i>Staphylococcus aureus subsp. aureus</i>		ATCC® 700699™*	89093
<i>Staphylococcus aureus subsp. aureus</i>	35	ATCC® 9144™*	89182
<i>Staphylococcus aureus subsp. aureus</i>		ATCC® BAA-44™*	89170
<i>Staphylococcus epidermidis</i>	36	ATCC® 12228™*	89045
<i>Staphylococcus epidermidis</i>	132	ATCC® 14990™*	89202
<i>Staphylococcus haemolyticus</i>		ATCC® 29970™*	89126
<i>Staphylococcus saprophyticus</i>	159	ATCC® 15305™*	89153
<i>Staphylococcus xylosum</i>		ATCC® 29971™*	89133
<i>Stenotrophomonas maltophilia</i>		ATCC® 13637™*	89149
<i>Stenotrophomonas maltophilia</i>		ATCC® 17666™*	89194
<i>Streptococcus agalactiae</i>		ATCC® 13813™*	89046
<i>Streptococcus anginosus</i>		ATCC® 33397™*	89127
<i>Streptococcus bovis</i>	133	ATCC® 33317™*	89061
			group B; non-hemolytic in absence of CAMP Factor group G; type 1

ATCC® derived CultiControl freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices

Description	WDCM	CultiControl Ref.	notes
<i>Streptococcus dysgalactiae</i> <i>subsp. equisimilis</i>	derived from	ATCC® 12388™*	89128 group C
<i>Streptococcus mitis</i>	derived from	ATCC® 6249™*	89129
<i>Streptococcus mutans</i>	derived from	ATCC® 25175™*	89062
<i>Streptococcus pneumoniae</i>	derived from	ATCC® 27336™*	89063
<i>Streptococcus pneumoniae</i>	derived from	ATCC® 49619™*	89047 low level penicillin resistance by oxacillin test
<i>Streptococcus pneumoniae</i>	derived from	ATCC® 700671™*	89175
<i>Streptococcus pyogenes</i>	derived from	ATCC® 19615™*	89048 group A
<i>Streptococcus pyogenes</i>	derived from	ATCC® 49399™*	89130 group A
<i>Streptococcus salivarius</i>	derived from	ATCC® 13419™*	89131
<i>Streptococcus salivarius subsp. thermophilus</i>	derived from	ATCC® 19258™*	134 89186
<i>Streptococcus sanguinis</i>	derived from	ATCC® 10556™*	89064
<i>Trichophyton mentagrophytes</i>	derived from	ATCC® 9533™*	89140
<i>Vibrio alginolyticus</i>	derived from	ATCC® 17749™*	89144
<i>Vibrio parahaemolyticus</i>	derived from	ATCC® 17802™*	37 89056
<i>Yersinia enterocolitica subsp. enterocolitica</i>	derived from	ATCC® 23715™*	160 89168 biotype 1; serotype 8
<i>Yersinia enterocolitica subsp. enterocolitica</i>	derived from	ATCC® 9610™*	38 89050 biovar 1; serogroup O:8



* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Liofilchem® s.r.l. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

Look for the ATCC Licensed Derivative® Emblem for products derived from ATCC® cultures.

The ATCC Licensed Derivative® program was created to ensure quality and address global responsibilities for handling biological materials. A royalty fee will be charged on behalf of ATCC® for products derived from ATCC® biological material. In order to receive these products, an authorized representative from your organization must agree to the terms of the End User Agreement. It is mandatory that your organization registers for this program before we are able to process your order for lyophilized microorganism preparations containing ATCC® biological material. If you have additional questions about the program contact Liofilchem Customer Service.

NCTC® derived CultiControl freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices

Description	WDCM	CultiControl Ref.	notes
<i>Acinetobacter baumannii</i>	derived from NCTC 12156	70174	
<i>Actinomyces odontolyticus</i>	derived from NCTC 9935	70114	
<i>Aeromonas hydrophila</i>	derived from NCTC 12902	70169	
<i>Aeromonas hydrophila</i>	derived from NCTC 8049	63	70119
<i>Aspergillus brasiliensis</i>	derived from NCPF 2275	53	70021
<i>Bacillus cereus</i>	derived from NCTC 10320	1	70022
<i>Bacillus cereus</i>	derived from NCTC 10320	1	
<i>Bacillus cereus</i>	derived from NCTC 7464		70155
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i>	derived from NCTC 10400	3	70023
<i>Bacteroides fragilis</i>	derived from NCTC 10581		70113
<i>Bacteroides fragilis</i>	derived from NCTC 9343		70078
<i>Bacteroides ovatus</i>	derived from NCTC 11153		70111
<i>Bordetella bronchiseptica</i>	derived from NCTC 8344		70139
<i>Burkholderia cepacia</i>	derived from NCTC 10743		70147
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	derived from NCTC 11322	156	70167
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	derived from NCTC 11351		70145
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	derived from NCTC 13367	5	70086
<i>Candida albicans</i>	derived from NCPF 3179	54	70024
<i>Candida albicans</i>	derived from NCPF 3939		70072
<i>Citrobacter freundii</i>	derived from NCTC 9750		70159
<i>Clostridium bifermentans</i>	derived from NCTC 506		70020
<i>Clostridium difficile</i>	derived from NCTC 11209		70090
<i>Clostridium histolyticum</i>	derived from NCTC 503		70112
<i>Clostridium perfringens</i>	derived from NCTC 8237	7	70053
<i>Clostridium sordellii</i>	derived from NCTC 13356		70059
<i>Clostridium sporogenes</i>	derived from NCTC 276		70032
<i>Clostridium sporogenes</i>	derived from NCTC 532	8	70095
<i>Cronobacter sakazakii</i>	derived from NCTC 11467	214	70138
<i>Enterobacter aerogenes</i>	derived from NCTC 10006	175	70156
<i>Enterococcus faecalis</i>	derived from NCTC 12697	87	70026
<i>Enterococcus faecalis</i>	derived from NCTC 13379	85	70173
<i>Enterococcus faecalis</i>	derived from NCTC 13763	210	70115
<i>Enterococcus faecalis</i>	derived from NCTC 775	9	70025
<i>Enterococcus faecium</i>	derived from NCTC 7171	10	70171
<i>Erysipelothrix rhusiopathiae</i>	derived from NCTC 8163		70187
<i>Escherichia coli</i>	derived from NCTC 11954		70163
<i>Escherichia coli</i>	derived from NCTC 12241	13	70027
<i>Escherichia coli</i>	derived from NCTC 12923	12	70028
<i>Fluoribacter bozemanae</i>	derived from NCTC 11368		70157
<i>Gardnerella vaginalis</i>	derived from NCTC 10915		70099
<i>GeoBacillus stearothermophilus</i>	derived from NCTC 10007		70203
<i>Haemophilus influenzae</i>	derived from NCTC 12699		70077
<i>Haemophilus influenzae</i>	derived from NCTC 12975		70076
<i>Haemophilus influenzae</i>	derived from NCTC 13377		70120
<i>Haemophilus influenzae</i>	derived from NCTC 4560		70160
<i>Haemophilus influenzae</i>	derived from NCTC 8143		70176
<i>Haemophilus influenzae</i>	derived from NCTC 8468		70136
<i>Haemophilus influenzae</i>	derived from NCTC 8469		70142
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	derived from NCTC 13368		70070
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	derived from NCTC 13635	192	70192
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	derived from NCTC 9633	97	70089
<i>LactoBacillus acidophilus</i>	derived from NCTC 12980	98	70080
<i>LactoBacillus delbrueckii</i> subsp. <i>bulgaricus</i>	derived from NCTC 12712	102	70209
<i>Lactococcus lactis</i> subsp. <i>lactis</i>	derived from NCTC 6681	16	70082
<i>Legionella pneumophila</i> subsp. <i>Fraseri</i>	derived from NCTC 11233	180	70151
<i>Legionella pneumophila</i> subsp. <i>pneumophila</i>	derived from NCTC 11192	107	70052
<i>Listeria grayi</i>	derived from NCTC 10812		70101
<i>Listeria innocua</i>	derived from NCTC 11288	17	70029
<i>Listeria ivanovi</i> subsp. <i>ivanovii</i>	derived from NCTC 11846	18	70030
<i>Listeria monocytogenes</i>	derived from NCTC 10357		70188
<i>Listeria monocytogenes</i>	derived from NCTC 10527	21	70085
<i>Listeria monocytogenes</i>	derived from NCTC 13372		70060
<i>Listeria monocytogenes</i>	derived from NCTC 13627	20	70031

NCTC® derived CultiControl freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices

Description	WDCM	CultiControl Ref.	notes
<i>Listeria monocytogenes</i>	derived from NCTC 7973	109	70148
<i>Micrococcus luteus</i>	derived from NCTC 2665	111	70102
<i>Micrococcus luteus</i>	derived from NCTC 7743		70096
<i>Moraxella (Branhamella) catarrhalis</i>	derived from NCTC 11020		70103
<i>Neisseria gonorrhoeae</i>	derived from NCTC 12700		70104
<i>Neisseria gonorrhoeae</i>	derived from NCTC 8375		70074
<i>Neisseria meningitidis</i>	derived from NCTC 10026		70164
<i>Nocardia brasiliensis</i>	derived from NCTC 11294		70189
<i>Peptostreptococcus anaerobius</i>	derived from NCTC 11460		70165
<i>Plesiomonas shigelloides</i>	derived from NCTC 10360		70094
<i>Porphyromonas gingivalis</i>	derived from NCTC 11834		70162
<i>Prevotella melaninogenica</i>	derived from NCTC 12963		70134
<i>Proteus hauseri</i>	derived from NCTC 4175		70190
<i>Proteus mirabilis</i>	derived from NCTC 11938	23	70083
<i>Pseudomonas aeruginosa</i>	derived from NCTC 10332	24	70108
<i>Pseudomonas aeruginosa</i>	derived from NCTC 12903	25	70033
<i>Pseudomonas aeruginosa</i>	derived from NCTC 12924	26	70034
<i>Pseudomonas aeruginosa</i>	derived from NCTC 13359		70109
<i>Pseudomonas fluorescens</i>	derived from NCTC 10038	115	70110
<i>Rhodococcus equi</i>	derived from NCTC 1621	28	70035
<i>Saccharomyces cerevisiae</i>	derived from NCTC 10716	58	70036
<i>Salmonella enterica subsp. arizonae</i>	derived from NCTC 8297		70154
<i>Salmonella enterica subsp. enterica serovar Abony</i>	derived from NCTC 6017	29	70132
<i>Salmonella enterica subsp. enterica serovar Enteritidis</i>	derived from NCTC 12694	30	70084
<i>Salmonella enterica subsp. enterica serovar Typhimurium</i>	derived from NCTC 12023	31	70037
<i>Salmonella enterica subsp. enterica serovar Typhimurium</i>	derived from NCTC 74	121	70054
<i>Serratia marcescens</i>	derived from NCTC 13382		70121
<i>Shigella flexneri</i>	derived from NCTC 12698		
<i>Staphylococcus aureus</i>	derived from NCTC 10788	193	70044
<i>Staphylococcus aureus</i>	derived from NCTC 12493	212	70039
<i>Staphylococcus aureus subsp. aureus</i>	derived from NCTC 10655		70137
<i>Staphylococcus aureus subsp. aureus</i>	derived from NCTC 12973	131	70041
<i>Staphylococcus aureus subsp. aureus</i>	derived from NCTC 12981	34	70040
<i>Staphylococcus aureus subsp. aureus</i>	derived from NCTC 13373	211	70043
<i>Staphylococcus aureus subsp. aureus</i>	derived from NCTC 6571	35	70182
<i>Staphylococcus epidermidis</i>	derived from NCTC 11047	132	70202
<i>Staphylococcus epidermidis</i>	derived from NCTC 13360	36	70045
<i>Staphylococcus haemolyticus</i>	derived from NCTC 11042		70126
<i>Staphylococcus saprophyticus</i>	derived from NCTC 7292	159	70153
<i>Staphylococcus xylosum</i>	derived from NCTC 11043		70133
<i>Stenotrophomonas maltophilia</i>	derived from NCTC 10257		70149
<i>Streptococcus agalactiae</i>	derived from NCTC 8181		70046
<i>Streptococcus anginosus</i>	derived from NCTC 10713		70127
<i>Streptococcus bovis</i>	derived from NCTC 8177	133	70061
<i>Streptococcus mutans</i>	derived from NCTC 10449		70062
<i>Streptococcus pneumoniae</i>	derived from NCTC 12977		70047
<i>Streptococcus pyogenes</i>	derived from NCTC 12696		70048
<i>Streptococcus salivarius subsp. thermophilus</i>	derived from NCTC 12958	134	70186
<i>Streptococcus sanguinis</i>	derived from NCTC 7863		70064
<i>Vibrio alginolyticus</i>	derived from NCTC 12160		70144
<i>Vibrio parahaemolyticus</i>	derived from NCTC 10903	37	70056
<i>Yersinia enterocolitica subsp. enterocolitica</i>	derived from NCTC 10598	160	70168
<i>Yersinia enterocolitica subsp. enterocolitica</i>	derived from NCTC 12982	38	70050



NCTC® is a registered trademark of the UK Health Security Agency (UKHSA). Liofilchem is licensed to use these trademarks and to sell products derived from NCTC® cultures.

ATCC® derived Quanti-CultiControl		
freeze-dried microorganisms		
Packaging: 1 vial containing 1 pellet + 1 rehydration fluid vial		
Quantitative: <100 CFU / 0.1 mL inoculum		
Applications: Growth Promotion Testing, Microbial Limits Testing, Microbial Enumeration Testing		
Description	Quanti-CultiControl Ref.	notes
<i>Aspergillus brasiliensis</i> derived from ATCC® 16404™*	89501	
<i>Bacillus cereus</i> derived from ATCC® 10876™*	89502	
<i>Bacillus subtilis</i> subsp. spizizenii derived from ATCC® 6633™*	89503	
<i>Bacteroides fragilis</i> derived from ATCC® 25285™*	89505	
<i>Bifidobacterium animalis</i> subsp. animalis derived from ATCC® 25527™*	89539	
<i>Brevundimonas diminuta</i> derived from ATCC® 19146™*	89506	
<i>Burkholderia cepacia</i> derived from ATCC® 25416™*	89507	
<i>Candida albicans</i> derived from ATCC® 10231™*	89508	
<i>Clostridium sporogenes</i> derived from ATCC® 11437™*	89513	
<i>Clostridium sporogenes</i> derived from ATCC® 19404™*	89514	
<i>Enterobacter aerogenes</i> derived from ATCC® 13048™*	89516	
<i>Enterococcus faecalis</i> derived from ATCC® 29212™*	89517	
<i>Escherichia coli</i> derived from ATCC® 8739™*	89519	
<i>Geobacillus stearothermophilus</i> derived from ATCC® 12980™*	89521	
<i>Geobacillus stearothermophilus</i> derived from ATCC® 7953™*	89522	
<i>Kocuria rhizophila</i> derived from ATCC® 9341™*	89523	
<i>Lactobacillus fermentum</i> derived from ATCC® 9338™*	89524	
<i>Listeria monocytogenes</i> derived from ATCC® 19115™*	89525	serotype 4b
<i>Micrococcus luteus</i> derived from ATCC® 4698™*	89526	
<i>Pseudomonas aeruginosa</i> derived from ATCC® 27853™*	89527	
<i>Pseudomonas aeruginosa</i> derived from ATCC® 9027™*	89528	
<i>Salmonella enterica</i> subsp. enterica serovar Typhimurium derived from ATCC® 14028™*	89531	
<i>Salmonella enterica</i> subsp. enterica serovar Typhimurium derived from ATCC® 13311™*	89530	
<i>Staphylococcus aureus</i> subsp. aureus derived from ATCC® 25923™*	89533	recommended for CAMP test
<i>Staphylococcus aureus</i> derived from ATCC® 6538™*	89535	
<i>Staphylococcus aureus</i> derived from ATCC® 6538P™*	89534	
<i>Staphylococcus epidermidis</i> derived from ATCC® 12228™*	89537	
<i>Streptococcus pyogenes</i> derived from ATCC® 19615™*	89538	group A



* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Liofilchem® s.r.l. is licensed to use these trademarks and to sell products derived from ATCC® cultures. Look for the ATCC Licensed Derivative® Emblem for products derived from ATCC® cultures. The ATCC Licensed Derivative® program was created to ensure quality and address global responsibilities for handling biological materials. A royalty fee will be charged on behalf of ATCC® for products derived from ATCC® biological material. In order to receive these products, an authorized representative from your organization must agree to the terms of the End User Agreement. It is mandatory that your organization registers for this program before we are able to process your order for lyophilized microorganism preparations containing ATCC® biological material. If you have additional questions about the program contact Liofilchem Customer Service.

Dip-Slide for the screening of Urinary Infections

Product with minimum order



Description	Packaging	Ref.
Uritest (MacConkey Agar / Cetrимide Agar / CLED Andrade Agar) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20/pack	500152
Uritest N (MacConkey Agar / Cetrимide Agar / CLED Agar) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20/pack 120/pack 500/pack	500232 51023 51123
Uritest C (MacConkey MUG Agar / Cetrимide Agar / CLED Andrade Agar) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20/pack 120/pack	500242 # 51024 #
Uritest EC (MacConkey Agar / E.Coli Agar / CLED Agar) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20/pack 120/pack	500412 51041
Uritest 2 (CLED Agar / MacConkey Agar) Total bacterial count, enterobacteria.	20/pack 120/pack	500302 51030
Uritest M (MacConkey Agar / Malt Selective Agar / CLED Andrade Agar) Enterobacteria, yeasts, total bacterial count.	20/pack 120/pack	500182 51018
Uritest Malto (MacConkey Agar / Malt Agar / CLED Agar) Enterobacteria, yeasts, total bacterial count.	20/pack 120/pack	500402 51040
Uritest EF (MacConkey Agar / Bile Aesculin Agar / CLED Agar) Enterobacteria, enterococci, total bacterial count.	20/pack 120/pack	500702 51070
Uritest Penta (CLED Agar / MacConkey Agar / Cetrимide Agar / Malt Agar / Bile Aesculin Agar) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20/pack 120/pack	500142 51014

Dip-Slide for the screening of Urogenital Infections



Description	Packaging	Ref.
Vagitest (Chocolate / Rogosa / Thayer Martin) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10/pack	50020

Dip-Slide for the screening of Dermatological Infections

Description	Packaging	Ref.
Dermatest (Malto / Corn Meal / DTM-Taplin) Pathogenic dermatophytes and others non pathogenic fungi.	10/pack	50021
DermaSlide (DTM modified / Sabouraud CAF Agar) Pathogenic dermatophytes and others non pathogenic fungi.	20/pack	500222

Flexible Contact Slide for microbial screening of surfaces (ISO 18593)



Description	Packaging	Ref. # Product with minimum order
Contact Slide 1 (PCA + TTC + Neutralizing / VRBG Agar + Neutralizing) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20/pack 120/pack	525262 53526
Contact Slide 2 (PCA + TTC + Neutralizing / Rose Bengal CAF Agar + Neutralizing) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20/pack 120/pack	525272 53527
Contact Slide 3 (PCA+TTC + Neutralizing / Vogel Johnson Agar) Total bacterial count with inactivation of disinfectants, staphylococci.	20/pack 120/pack	525302 53530
Contact Slide 4 (Cetrimide Agar + Neutralizing / Rose Bengal CAF Agar + Neutralizing) <i>Pseudomonas</i> spp., yeasts and moulds.	20/pack 120/pack	525282 53528
Contact Slide 5 (V.R.B.G. Agar + Neutralizing / Bile Aesculin Agar + Neutralizing) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20/pack 120/pack	525312 53531
Contact Slide 8 (PCA + TTC + Neutralizing / VRBL Agar + Neutralizing) Total bacterial count, coliforms, with inactivation of disinfectants.	20/pack 120/pack	525622 53562
Contact Slide 9 (PCA + TTC + Neutralizing / PCA + TTC + Neutralizing) Total bacterial count with inactivation of disinfectants.	20/pack 120/pack	525322 53532
Contact Slide 10 (Rose Bengal CAF Agar + Neutralizing / Vogel Johnson Agar) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20/pack 120/pack	525452 # 53545 #
Contact Slide Chrom 1 (Chromatic E.coli-coliform / Chromatic Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20/pack 120/pack	525292 53529
Contact Slide Chrom 2 (Chromatic E.coli-coliform / PCA+TTC + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20/pack	525392
Contact Slide Chrom 3 (Chromatic E.coli-coliform / Baird Parker Agar + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20/pack	525462 #
Contact Slide TSA + Neutralizing / TSA + Neutralizing Total bacterial count with inactivation of disinfectants.	20/pack	525472
Contact Slide TSA + Neutralizing / Rose Bengal CAF Agar + Neutralizing Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20/pack 120/pack	525482 53548


Systems for microbial identification directly from food and surfaces

Description	Packaging	Ref.
<p>Food System System for detection and identification of pathogenic microorganisms from food. VALIDATED to ISO 16140 standard Microbiology of food and animal feeding stuffs - Protocol for the validation of alternative methods</p> 	<p>40 test 8 test</p>	<p>71680 79680</p>
<p>HACCP System Plus System for enumeration and identification of pathogenic microorganisms from surfaces.</p> 	<p>40 test 8 test</p>	<p>71655 79655</p>

Coliforms in water testing

Description	Packaging	Ref.
 <p>Water Detect Coliforms Chromogenic medium for detecting <i>E. coli</i> and others coliforms in water samples.</p>	<p>20/pack 12 bottles x 100 mL</p>	<p>630600 630601</p>

Pesticide residue detection

Description	Packaging	Ref.
 <p>MAIA Pesticide MultiTest Colorimetric system for the detection of organophosphate, organochloride and carbamate pesticides residues in water, food and drink.</p>	<p>96 test</p>	<p>79700</p>

Antibiotic residue detection		
Description	Packaging	Ref.
MeRA Test Detection of antibiotic and sulphamide residues in meat by <i>Geobacillus stearothermophilus</i> .	50 vials	80356
MiRA Test Detection of antibiotic and sulphamide residues in milk by <i>Geobacillus stearothermophilus</i> .	50 vials	80355
PAR TEST Agar Medium with spores of <i>Geobacillus stearothermophilus</i> subsp. <i>calidolactis</i> for residual antibiotic detection in biological liquids (milk, blood and urine).	20 plates 100 plates	11200 11200*
BLANK DISCS 6 mm Paper discs.	250 discs	9999
Antibiotic Test Antibiotic residue detection in urine and other biological fluids by <i>Bacillus subtilis</i> .	50 discs	80350
Rapid Antibiotic Test Antibiotic residue detection in foodstuffs by <i>Geobacillus stearothermophilus</i> var. <i>calidolactis</i> .	50 tests	80351
Suspension Spore BAS E8 Sterile Suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁸ CFU/mL).	10 mL 2 mL	82012 82212
Suspension Spore BAS E7 Sterile Suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁷ CFU/mL).	10 mL	82013
Suspension Spore BAS E6 Sterile Suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁶ CFU/mL).	10 mL	82014
Suspension Spore BAS E5 Sterile suspension of <i>Bacillus subtilis</i> spores in distilled water for antibiotic residue detection in milk and other biological fluids (1-5x10 ⁵ CFU/mL).	10 mL	82015

Solutions for sample collection		
Description	Packaging	Ref.
SWAB Solution SKELETON Solution to hydrate samples collection swabs (2001/471/CE).	20 vials	20157
SWAB Solution SURFACES Solution to hydrate samples collection swabs before surfaces checking.	20 vials	20160

Swabs w/ and w/o transport medium		
Description	Packaging	Ref.
BIO – CULT AMIES + charcoal Swabs for collection and transport of clinical, environmental samples and microorganisms.	150 Swabs	85307
BIO – CULT AMIES Swabs for collection and transport of clinical, environmental samples and microorganisms.	100 Swabs	85317
BIO CULT CARY BLAIR Swabs for collection and transport of clinical, environmental samples and microorganisms.	150 Swabs	85310
BIO – CULT Stuart + charcoal Swabs for collection and transport of clinical, environmental samples and microorganisms.	150 Swabs	85306
BIO – CULT Stuart Swabs for collection and transport of clinical, environmental samples and microorganisms.	100 Swabs	85308
Sterile Swab in Envelope Dry swabs for clinical and environmental samples collection.	100 Swabs	85270
Sterile Swab without Envelope Swabs for clinical and environmental samples collection.	100 Swabs	85275

Transport medium for Viruses, Chlamydia, Mycoplasma and Ureaplasma		
Description	Packaging	Ref.
VTM 3 mL	100/pack	26490
VTM 200 mL	6/pack	412450



Organic residue detection on surfaces

Description	Packaging	Ref.
 <p>Clean Test Swab for detecting protein, fat and sugar residues for the hygiene control of surfaces.</p>	50/pack 3/pack	86001 86001/3

Surface swabs for microbial detection

Description	Packaging	Ref.
		
Contam Swab Listeria	30/pack	86102
Contam Swab E.coli-coliforms	30/pack	86101
Contam Swab Salmonella	30/pack	86103
Contam Swab Contamination	30/pack	86100
Contam Swab MRSA	30/pack	86104 #
Contam Swab Yeasts & Moulds	30/pack	86105 #

ESC (Easy Surface Checking) swab (ISO 18593)

Description	Packaging	Ref.
  		
ESC Swab Neutralizing Rinse Solution 10 mL	100/pack	85601
ESC Swab Neutralizing Rinse Solution 5 mL	100/pack	85602
ESC Swab Buffered Peptone Water 10 mL	100/pack	85603
ESC Swab Buffered Peptone Water 5 mL	100/pack	85604
ESC Swab D/E Neutralizing Broth 10 mL	100/pack	85605
ESC Swab D/E Neutralizing Broth 5 mL	100/pack	85606
ESC Swab Lethen Broth 10 mL	100/pack	85607
ESC Swab Lethen Broth 5 mL	100/pack	85608
ESC Swab Lethen Broth 1 mL	100/pack	85621
ESC Swab Maximum Recovery Diluent 10 mL	100/pack	85609
ESC Swab Maximum Recovery Diluent 5 mL	100/pack	85610
ESC swab Mycoplasma Transport Broth	100/pack	85612
ESC swab Demi Fraser Broth 10 mL	100/pack	85614
Sampling Template 10 cm x 10 cm sterile	60/pack	96762

Self-contained for steam sterilization		# Product with minimum request
		
Description	Packaging	Ref.
STERILtest GST E6 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/disc).	100/pack	91100
STERILtest GST E5 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/disc).	100/pack	91101 #

Ampoules for steam sterilization		
		
Description	Packaging	Ref.
Steril Control GST E6 ampoules Ampoules with nutrient medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/ampoule).	20 Ampoules 4 mL	91210
	50 Ampoules 4 mL	91050
	50 Ampoules 1 mL + 10 negative controls	91040
Steril Control GST E5 ampoules Ampoules with nutrient medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/ampoule).	20 Ampoules 4 mL	91251
	50 Ampoules 4 mL	91052
Steril Control GST E4 ampoules Ampoules with nutrient medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 ⁴ CFU/ampoule).	20 Ampoules 4 mL	91253
	50 Ampoules 4 mL	91054
		

Strips for steam sterilization		# Product with minimum request
		
Description	Packaging	Ref.
Strip Control GST E6 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips	91055
	20 tubes	
Strip Control GST E5 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	100 strips	91155
	20 strips	91056 #
Strip Control GST E4 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips	91057 #
	20 tubes	
Strip Control GST E3 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips	91058 #
	20 tubes	

Spores suspensions for steam sterilization		
Description	Packaging	Ref.
Suspension Spore GST E8 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82002
Suspension Spore GST E7 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82003
Suspension Spore GST E6 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82004
Suspension Spore GST E5 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82005

Steel coupon for steam sterilization		
Description	Packaging	Ref.
Steam Control GST E6 SS Steel coupons with spores of <i>Geobacillus stearothermophilus</i> (1-5x10 ⁶ CFU/mL) and tubes of nutrient broth with pH indicator.	20 coupons + 20 tubes	91087
	50 coupons	91097

D value calculation		
Description		Ref.
D Value & Count Test for steam sterilization		S0003

Strips for irradiation sterilization		
Description	# Product with minimum request	
	Packaging	Ref.
Strip Control PUM E7 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁷ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91071 #
Strip Control PUM E6 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91072
Strip Control PUM E5 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91073 #
Strip Control PUM E4 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91074 #
Strip Control PUM E3 Strips with spores of <i>B. pumilus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91075 #

Strips for dry heat or ethylene oxide sterilization		# Product with minimum request	
Description	Packaging	Ref.	
Strip Control BAT E8 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁸ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91059 #	
Strip Control BAT E7 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁷ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91062 #	
Strip Control BAT E6 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91063	
	100 strips	91067	
Strip Control BAT E5 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91064 #	
Strip Control BAT E4 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91065 #	
Strip Control BAT E3 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91066 #	

Spores suspension for dry heat or ethylene oxide sterilization



Description	Packaging	Ref.
Suspension Spore BAT E10 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ¹⁰ CFU/mL).	10 mL vial	82006
Suspension Spore BAT E9 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁹ CFU/mL).	10 mL vial	82007
Suspension Spore BAT E8 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82008
Suspension Spore BAT E7 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82009
Suspension Spore BAT E6 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82010
Suspension Spore BAT E5 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82011

Self-contained for ethylene oxide sterilization

Self-contained for ethylene oxide sterilization		# Product with minimum request	
Description	Packaging	Ref.	
STERILtest BAT E6 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/disc).	100/pack	91150	
STERILtest BAT E5 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/disc).	100/pack	91151 #	
Crusher Steel pincer for the easy break of self-contained vials.	1 pincer	96790	

Steel coupons for hydrogen peroxide sterilization		# Product with minimum request	
Description	Packaging	Ref.	
OXI Control E6 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91082	
OXI Control E6 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon).	100 coupons	91182	
OXI Control E5 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91083 #	
OXI Control E4 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁴ CFU/coupon).	100 coupons	91181	
OXI Control E4 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ³ CFU/coupon).	100 coupons	91180	

Poliflex coupons for hydrogen peroxide sterilization		# Product with minimum request	
Description	Packaging	Ref.	
OXI Control E6 PLX Poliflex coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91085	
OXI Control E5 PLX Poliflex coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91086 #	

Strips for hydrogen peroxide sterilization		# Product with minimum request	
Description	Packaging	Ref.	
OXI Control E6 strip Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91088	
OXI Control E5 strip Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91089 #	

Spore suspensions for hydrogen peroxide sterilization		Packaging	Ref.
Suspension Spore GST VHP E8 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82030	
Suspension Spore GST VHP E7 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82031	
Suspension Spore GST VHP E6 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82032	
Suspension Spore GST VHP E5 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82033	

PVC coupon for hydrogen peroxide sterilization		Packaging	Ref.
OXI Control E6 PVC PVC coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91092	

Aluminum coupons for hydrogen peroxide sterilization		Packaging	Ref.
OXI Control E6 Aluminum Aluminum coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91098	

Glass coupons for hydrogen peroxide sterilization		Packaging	Ref.
OXI Control E6 Glass Glass coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91096	

Ready culture media in 55 mm contact plates		# Product with minimum order
Description	Packaging	Ref.
Baird Parker Agar + Neutralizing Selective medium for <i>Staphylococcus aureus</i> isolation with inactivation of disinfectants.	20/pack	15331
Cetrimide Agar Selective medium for <i>Pseudomonas aeruginosa</i> isolation.	20/pack	15332
Cetrimide Agar + Neutralizing Selective medium for <i>Pseudomonas aeruginosa</i> isolation with inactivation of disinfectants.	20/pack	15386 #
Legionella BCYE Agar Non-selective medium for <i>Legionella</i> spp. isolation.	20/pack	15334
Legionella GVPC Agar Selective medium for <i>Legionella</i> spp. isolation.	20/pack	15376
MacConkey Agar Selective medium for gram-negatives isolation.	20/pack	15329
Mannitol Salt Agar Selective medium for staphylococci isolation.	20/pack	15328
m-Faecal Coliform Agar Selective medium for fecal coliforms isolation and enumeration.	20/pack	15368 #
Plate Count Agar Medium for surfaces and air total bacterial count.	20/pack	15325
Plate Count Agar + Neutralizing Medium for surfaces total bacterial count with inactivation of disinfectants.	20/pack	15336
Plate Count Agar + Neutralizing (Irradiated) Medium for surfaces total bacterial count with inactivation of disinfectants.	20/pack	15336S
Plate Count Agar + TTC + Neutralizing Medium for surfaces total bacterial count with inactivation of disinfectants.	20/pack	15360
Rose Bengal Agar + CAF Selective Medium for yeasts and moulds isolation and enumeration.	20/pack	15374
Rose Bengal Agar + CAF + Neutralizing Selective Medium for yeasts and moulds isolation and enumeration.	20/pack	15385 #
R2 Agar Medium for heterotrophic microorganisms count in drinking water.	20/pack	15354 #
Sabouraud Dextrose Agar Medium for yeasts and moulds isolation.	20/pack	15327
Sabouraud Dextrose Agar + Lactamator 500 IU + Neutralizing (Irradiated) Medium for yeasts and moulds isolation with inactivation of β -lactam antibiotics and disinfectants.	20/pack	15345S#
Sabouraud Dextrose Agar + 0.1% Penase + Neutralizing (Irradiated) Medium for yeasts and moulds isolation with inactivation of disinfectants and penicillins.	20/pack	15333S#
Sabouraud Agar + Neutralizing (Irradiated) Medium for yeasts and moulds isolation with inactivation of disinfectants.	20/pack 400/pack	15323S 15323KS
Sabouraud CAF Agar (Irradiated) Selective medium for yeasts and moulds isolation.	20/pack	15380S
Sabouraud CAF Agar + Neutralizing Selective medium for yeasts and moulds isolation with inactivation of disinfectants.	20/pack 20/pack irradiated	15365 15365S
Slanetz Bartley Agar + Neutralizing Selective medium for fecal streptococci enumeration and isolation with inactivation of disinfectants.	20/pack	15362
Tryptic Soy Agar (Irradiated) Medium for surfaces and air total bacterial count.	20/pack	15340S
Tryptic Soy Agar + Neutralizing (Irradiated) Medium for surfaces and air total bacterial count with inactivation of disinfectants.	20 plates 400/pack	15339S 15339KS
		
Tryptic Soy Agar + 0.1% Penase + Neutralizing Medium for surfaces and air total bacterial count with inactivation of disinfectants and penicillins.	20/pack	15330 #
Tryptic Soy Agar + 0.1% Penase + Neutralizing (Irradiated) Medium for surfaces and air total bacterial count with inactivation of disinfectants and penicillins.	20/pack	15330S #

Ready culture media in 55 mm contact plates		
Description	Packaging	# Product with minimum order Ref.
Tryptic Soy Agar + Lactamator 500 UI + Neutralizing (Irradiated) Irradiated medium for non fastidious bacteria growth with inactivation of β -lactam antibiotics and disinfectants.	20/pack	15341S
Violet Red Bile Glucose Agar Selective Medium for enterobacteria isolation and enumeration.	20/pack	15375
Violet Red Bile Lactose Agar Selective Medium for coliforms isolation and enumeration.	20/pack	15326
Violet Red Bile Glucose Agar + Neutralizing Selective Medium for <i>Enterobacteriaceae</i> isolation and enumeration with inactivation of disinfectants.	20/pack	15382
Violet Red Bile Lactose Agar + Neutralizing Selective Medium for coliforms isolation and enumeration with inactivation of disinfectants.	20/pack	15361

Lockable contact plates


- gamma irradiated
- triple wrapped
- ventilated
- 270 days shelf life
- storage at 10-25°C
- data matrix barcode on each plate



Description	packaging	Ref.
Tryptic Soy Agar + Neutralizing	20/pack	15203S
Sabouraud Dextrose Agar + Neutralizing	20/pack	15204S
DG 18 Agar + Neutralizing (not irradiated)	20/pack	15205

Ready culture media in 60 mm petri dishes



Product with minimum order



Description	Packaging	Ref.
MEMBRANE DISC FILTERS Mixed cellulose ester membrane for the isolation and enumeration of bacteria by Membrane Filter Technique.		
47 mm - 0.45 µm - white	100 filters	97000
47 mm - 0.45 µm - black	100 filters	97013
47 mm - 0.22 µm - white	100 filters	97014
Aeromonas Agar Selective medium for <i>Aeromonas</i> spp. isolation.	20/pack	163622 #
Baird Parker Agar Selective medium for <i>Staphylococcus aureus</i> isolation.	ISO 6888-1 20/pack	163512
Bile Aesculin Azide Agar Selective medium for group D streptococci isolation.	ISO 7899-2 20/pack	163572
C-EC MF Agar Chromogenic and fluorogenic medium for isolating and detecting <i>Escherichia coli</i> and other coliforms in water.	20/pack	163862
Chromatic™ EC X-Gluc Agar Chromogenic medium for the detection of <i>E. coli</i> in water and foodstuffs.	20/pack	163722 #
Iron Sulfite Agar Medium for thermophilic anaerobes determination.	20/pack	163372 #
Legionella Agar BCYE Agar Non-selective medium for <i>Legionella</i> spp. isolation.	ISO 11731 20/pack	163362
Legionella Agar (GVPC) Selective medium for <i>Legionella</i> spp. isolation.	ISO 11731 20/pack	163392
MacConkey Agar Selective medium for enterobacteria isolation.	E.P. 20/pack	163542
m-CP Agar Selective medium for <i>Clostridium perfringens</i> isolation from water.	98/83/EC 20/pack	163612
m-ENDO Agar LES Medium for coliforms determination in water.	20/pack	163432
m-Faecal Coliform Agar Selective medium for fecal coliforms isolation and enumeration.	APHA 20/pack	163442
M-Green Yeast & Mould Agar + CAF Selective medium for detection and count of yeasts and moulds in beverages.	20/pack	163812 #
M-Green Yeast And Mold Agar pH modified Selective medium for detection and count of yeasts and moulds in beverages.	20/pack	163842L
M-Green Yeast & Mould Agar Selective medium for detection and count of yeasts and moulds in beverages.	20/pack	163642L
Mycosel Agar Selective medium for pathogenic fungi isolation.	20/pack	163412 #
Orange Serum Agar Medium for the isolation and enumeration of microorganisms associated with spoilage of citrus fruits and derivatives.	20/pack	163732 #
Plate Count Agar Medium for bacterial count in water and foodstuffs determination.	ISO 4833 20/pack	163452
Cetrimide Agar Selective medium for <i>Pseudomonas aeruginosa</i> isolation.	E.P. 20/pack	163472
Pseudomonas CN Agar Selective medium for <i>Pseudomonas</i> spp. isolation.	ISO 12780 20/pack	163592
R2 Agar Medium for heterotrophic microorganisms count in drinking water.	E.P. 20/pack	163672
Sabouraud Dextrose Agar Medium for yeasts and moulds isolation.	E.P. 20/pack	163402
Slanetz-Bartley Agar (m-Enterococcus A.) Selective medium for fecal streptococci enumeration and isolation.	ISO 7899-2 20/pack	163462
SPS Agar Medium for <i>C. perfringens</i> and sulfite reducing bacteria isolation and count.	APHA 20/pack	163522

Ready culture media in 60 mm petri dishes

Product with minimum order



Description

Packaging

Ref.

Streptococcal KF+TTC Agar

Selective medium for fecal streptococci isolation and enumeration.

APHA

20/pack

163482

TBX Agar
Selective chromogenic medium for detecting and enumerating *E. coli* in food.

ISO 16649

20/pack

163652

Tergitol Agar

Selective medium for coliforms enumeration and identification.

ISO 9308-1

20/pack

163562

Tergitol + TTC Agar

Selective medium for coliforms enumeration and identification.

ISO 9308-1

20/pack

163492

Tryptone Glucose Extract Agar + TTC

Medium for cultivating and enumerating microorganisms in water and dairy products.

20/pack

163742L

Tryptose Sulfite Cycloserine Agar

Selective medium for the isolation and enumeration of *Clostridium perfringens* in water, according to ISO 14189.

20/pack

163872

Wort Agar

Medium for the cultivation, isolation and enumeration or enrichment of fungi, especially of yeasts.

20/pack

163802

Yeast Extract Agar


Medium for microorganisms enumeration in water.

ISO 6222

20/pack

163582 #

Dry media in 60 mm petri dishes
100 pads, packed in 10 piles of 10



Description	Packaging	Ref.
Easy Dry Chromatic Coliform ISO	100/pack	87500
Easy Dry Tryptic Soy	100/pack	87512
Easy Dry Yeast Extract	100/pack	87513
Easy Dry MacConkey	100/pack	87514
Easy Dry Chromatic Salmonella	100/pack	87515
Easy Dry C-EC MF	100/pack	87516
Easy Dry Tergitol	100/pack	87517
Easy Dry Tergitol + TTC	100/pack	87518
Easy Dry Chromatic Coli Coliform	100/pack	87519
Easy Dry M-Green Yeast & Mould	100/pack	87511
Easy Dry Sabouraud Dextrose	100/pack	87510
Easy Dry Streptococcal KF + TTC	100/pack	87509
Easy Dry Plate Count + TTC	100/pack	87501
Easy Dry Violet Red Bile Glucose	100/pack	87502
Easy Dry TBX	100/pack	87503
Easy Dry m-Endo LES	100/pack	87504
Easy Dry Pseudomonas CN	100/pack	87505
Easy Dry m-Faecal Coliform	100/pack	87506
Easy Dry Slanetz Bartley	100/pack	87507
Easy Dry Bile Aesculin Azide	100/pack	87508
Easy Dry Orange Serum	100/pack	87520

Ready BAGS

Ready to use culture media in 3 L or 5 L bags.






Description	Packaging	Ref.
Buffered NaCl Peptone Solution pH 7.0 3 L bag	3/pack	499070
Buffered Peptone Water 3 L bag	3/pack	499030
Buffered Peptone Water 5 L bag	3/pack	499035
Listeria Demi Fraser Broth 3 L bag	3/pack	499020
Listeria Demi Fraser Broth 5 L bag	3/pack	499025
Listeria UVM Mod. Broth 3 L bag	3/pack	499010
Listeria UVM Mod. Broth 5 L bag	3/pack	499015
Maximum Recovery Diluent 3 L bag	3/pack	499040
Maximum Recovery Diluent 5 L bag	3/pack	499045
Tryptic Soy Broth 3 L bag	3/pack	499060





90 mm ready-to-use agar plates



Ref.	Description	Packaging	Intended use	Min. order qty (packs)
10406	AEROMONAS Agar	20/pack		
11030	ANAEROBIC Agar	20/pack		20
10447	Antibiotic Agar No. 1	20/pack		40
11041	AZIDE Agar (Sheep Blood 5%)	20/pack	IVD	
10027	BACILLUS CEREUS Agar (Mossel)	20/pack	IVD	
10007	BACILLUS CEREUS Agar (PEMBA)	20/pack	IVD	
10020	BAIRD PARKER Agar	20/pack	IVD	
				
10521	BAIRD PARKER Agar + RPF	20/pack		
10473	BAT Agar	20/pack		40
10458	Bifidobacterium Agar	20/pack		
10021	BIGGY (NICKERSON) Agar	20/pack	IVD	
10057	BILE AESCULIN Agar	20/pack	IVD	
10411	BILE ESCULIN AZIDE Agar w VANCOMYCIN	20/pack		10
10047	BISMUTH SULFITE Agar	20/pack		40
10446	Bismuth Sulphite Agar ISO 6579-1	20/pack	IVD	
10142	Blood Agar (Sheep Blood 7%)	20/pack	IVD	
10353	Bordet Gengou Agar (Sheep Blood 15%)	20/pack	IVD	
10060	Brain Heart Infusion Agar	20/pack	IVD	
10022	BRILLIANT GREEN Agar	20/pack	IVD	
10245	Brucella Blood Agar w Hemin and Vitamin K1	20/pack	IVD	
10460	Burkholderia Cepacia Selective Agar USP	20/pack		
11506	BURKHOLDERIA CEPACIA SELECTIVE Agar(BCSA)	20/pack	IVD	
10148	Campylobacter Agar (Sheep Blood 10%)	20/pack	IVD	
10050	Campylobacter Agar (Sheep Blood 5%)	20/pack	IVD	
10409	Campylobacter CCDA Agar	20/pack	IVD	
10145	Campylobacter Karmali Agar	20/pack	IVD	
10146	Campylobacter Preston Agar	20/pack		20
10602	Campylobacter Skirrow Agar	20/pack	IVD	
10079	CASITONE Agar	20/pack		40
10033	Cetrimide Agar	20/pack	IVD	
10023	Chocolate Agar	20/pack	IVD	
11023	Chocolate BACITRACIN Agar	20/pack	IVD	
11628	Chromatic BACILLUS CEREUS	20/pack		
				
11612	Chromatic CANDIDA	20/pack	IVD	
				

90 mm ready-to-use agar plates

Ref.	Description	Packaging	Intended use	Min. order qty (packs)
11632	Chromatic Clostridium difficile	20/pack	IVD	
				
11613	Chromatic Coli Coliform	20/pack		
				
11630	Chromatic Coliform Agar ISO	20/pack		
				
11640	Chromatic Colistin	20/pack	IVD	
				
11619	Chromatic CRE	20/pack	IVD	
				
11636	Chromatic Cronobacter Isolation Agar	20/pack		
				
11611	Chromatic Detection	20/pack	IVD	
				

90 mm ready-to-use agar plates

Ref.	Description	Packaging	Intended use	Min. order qty (packs)
11610	Chromatic E.COLI O157	20/pack	IVD	
				
11622	Chromatic ESBL	20/pack	IVD	
				
11629	Chromatic ESBL + AmpC	20/pack		20
11639	Chromatic GBS	20/pack	IVD	
				
11618	Chromatic MH	20/pack	IVD	
				
10599	Chromatic MRSA	20/pack	IVD	
				
11631	Chromatic OXA-48	20/pack	IVD	
				

90 mm ready-to-use agar plates				
Ref.	Description	Packaging	Intended use	Min. order qty (packs)
11635	Chromatic Pseudomonas	20/pack		30
11614	Chromatic Salmonella	20/pack	IVD	
				
11616	Chromatic STAPH AUREUS	20/pack	IVD	
				
11641	Chromatic Super CAZ/AVI	20/pack		
				
11633	Chromatic Vibrio	20/pack	IVD	
11621	Chromatic VRE	20/pack	IVD	
				
10026	CLED Agar	20/pack	IVD	
10004	CLED Andrade Agar	20/pack		30
11060	Clostridium Agar (Sheep Blood 5%)	20/pack	IVD	
10025	Columbia Agar (Horse Blood 5%)	20/pack	IVD	
11025	Columbia Agar (Sheep Blood 5%)	20/pack	IVD	
11517	Columbia Agar(Sheep Blood 5%) + Vancomycin	20/pack		20
11024	Columbia CNA Agar (Sheep Blood 5%)	20/pack	IVD	
11124	Columbia CNA Modified Agar (Sheep blood 5%)	20/pack	IVD	
11507	CORN MEAL Agar	20/pack	IVD	
10017	Czapek Dox Agar	20/pack	IVD	
10457	D/E Neutralizing Agar	20/pack		40
11052	DERMATOPHYTE (D.T.M.) Agar	20/pack	IVD	
11027	Desoxycholate Agar	20/pack		40
10098	Dichloran Glycerol (DG18) Agar	20/pack		
10097	DICHLORAN Rose Bengal CAF Agar	20/pack		
10013	DNase Test Agar	20/pack	IVD	
10018	Drigalski Lactose Agar	20/pack		20
10048	E.M.B. Levine Agar	20/pack	IVD	
10010	Endo Agar	20/pack		40
11623	Enterobacter (Cronobacter) sakazakii Isolation Agar	20/pack		
11057	Enterococco Agar	20/pack	IVD	
11501	ENTEROCOCCUS Agar w Vancomycin	20/pack	IVD	
10062	Fastidious Anaerobe Agar w/ Horse Blood (FAA-HB)	20/pack	IVD	
11054	GARDNERELLA Agar (Sheep Blood 5%)	20/pack	IVD	
10442	Gelatin Peptone Agar	20/pack		

90 mm ready-to-use agar plates

Ref.	Description	Packaging	Intended use	Min. order qty (packs)
10080	Haemophilus Test Agar	20/pack	IVD	
10043	HEKTOEN ENTERIC Agar	20/pack	IVD	
10082	Helicobacter pylori Agar	20/pack	IVD	
10605	Helicobacter pylori EGG YOLK EMULSION Agar	20/pack		20
10028	ISOSENSITEST Agar	20/pack	IVD	
11072	King's B Medium	20/pack		
10128	Legionella Agar (GVPC)	20/pack		
10448	Legionella BCYE + AB Agar	20/pack	IVD	
10051	Legionella BCYE Agar	20/pack	IVD	
10424	Legionella BCYE Agar w Vancomycin + Colistin	20/pack		15
10412	Legionella BCYE Agar w/o Cysteine	20/pack	IVD	
10127	Legionella MWY Agar	20/pack		
11230	Listeria Aesculin Agar	20/pack		
10041	Listeria PALCAM Agar	20/pack	IVD	
10091	M 17 Agar	20/pack		10
11053	m-ENDO Agar LES	20/pack		
11210	m-FAECAL COLIFORM Agar	20/pack		20
10090	M.R.S. Agar	20/pack	IVD	
10049	M.S.R.V. Agar	20/pack		
11503	MacConkey Agar w/o CRYSTAL VIOLET	20/pack		20
11508	MacConkey Agar w/o NaCl	20/pack		20
10603	MacConkey II Agar	20/pack	IVD	
10129	MacConkey MUG Agar	20/pack		10
				
11514	MacConkey S-CT Agar E.coli O157	20/pack	IVD	
10005	MacConkey Sorbitol Agar	20/pack	IVD	
10029	MacConkey Agar	20/pack	IVD	
10438	MacConkey Agar No. 2	20/pack		20
10418	Malt Agar	20/pack		20
10066	Malt Extract Agar	20/pack		50
10030	Mannitol Salt Agar	20/pack	IVD	
10426	Marine Salt Agar	20/pack		40
10416	Middlebrook 7H11 Agar	20/pack	IVD	
10433	Milk Plate Count Agar	20/pack		40
11518	Mueller Hinton Agar + Cloxacillin	20/pack		10
10335	Mueller Hinton Chocolate Agar	20/pack	IVD	
10132	Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L beta-NAD)	20/pack	IVD	
10031	Mueller Hinton II Agar	20/pack	IVD	
10131	Mueller Hinton II Agar (Sheep Blood 5%)	20/pack	IVD	
11206	Mueller Hinton II Agar + 2% NaCl	20/pack	IVD	
11205	Mycoplasma Agar	20/pack	IVD	
11070	Mycosel Agar	20/pack	IVD	
10334	Neomycin Blood Agar (Sheep Blood 5%)	20/pack		30
10044	Nutrient Agar ISO 16266	20/pack	IVD	
11512	Nutrient Agar ISO 21528	20/pack	IVD	
11513	Nutrient Agar ISO 6579	20/pack	IVD	
10620	O.A. Listeria Agar	20/pack	IVD	
				

90 mm ready-to-use agar plates


Ref.	Description	Packaging	Intended use	Min. order qty (packs)
11055	O.G.Y.E. Agar	20/pack		
10003	O.G.Y.E. Agar + Gentamicin	20/pack		20
10089	ORANGE SERUM Agar	20/pack		10
10600	Oxacillin Resistance Staphylococcus Agar	20/pack		40
11200	PAR Test Agar	20/pack	IVD	
10032	Plate Count Agar	20/pack		
11199	POTATO DEXTROSE Agar	20/pack		
11034	Pseudomonas CN Agar	20/pack		
11033	Pseudomonas Isolation Agar	20/pack		20
10061	Pseudomonas Selective Agar	20/pack		
10014	Purple Lactose Agar	20/pack	IVD	
11509	RPMI Agar	20/pack	IVD	
10019	R2A Agar	20/pack		
10456	Rapid Y/M Agar	20/pack		40
				
10039	Rogosa Agar	20/pack	IVD	
10034	Rose Bengal CAF Agar	20/pack		
10036	S.S. Agar	20/pack	IVD	
11335	Sabouraud Agar + Gentamicin	20/pack	IVD	
11135	Sabouraud Agar Modified	20/pack	IVD	
11035	Sabouraud CAF Agar	20/pack	IVD	
11236	Sabouraud CAF Agar + Actidione	20/pack	IVD	
10235	Sabouraud CAF Agar + Gentamicin	20/pack	IVD	
10437	Sabouraud CAF Agar + Neutralizing	20/pack		40
10441	Sabouraud CAF Agar 50 mg	20/pack		40
10035	Sabouraud Dextrose Agar	20/pack	IVD	
10425	SCEDOSPORIUM SELECTIVE Agar	20/pack	IVD	
10405	Schaedler CNA Agar (Sheep Blood 5%)	20/pack	IVD	
11065	Schaedler K Agar (Sheep Blood 5%)	20/pack	IVD	
10065	Scheduler KKV Agar(Sheep Blood 5%)	20/pack	IVD	
10067	Schaedler KVN Agar (Sheep Blood 5%)	20/pack		30
10046	Serum Tellurite Agar	20/pack		40
11058	SLANETZ BARTLEY Agar(m-ENTEROCOCCUS)	20/pack	IVD	
11196	SPS Agar	20/pack	IVD	
10038	STREPTOCOCCAL KF+TTC Agar	20/pack		
10419	STROGYLOIDES STERCORALIS Agar	20/pack		30
10415	SUGAR FREE Agar	20/pack		40
10015	T.S.A. + ACTIDIONE	20/pack		20
10522	TBX Agar	20/pack		
				
11195	TCBS Agar	20/pack	IVD	
10417	TERGITOL Agar	20/pack		40
10428	Tergitol TTC Agar (ISO 9308-1)	20/pack		
11040	THAYER MARTIN Agar	20/pack	IVD	
11250	TINSDALE Agar	20/pack		20
10037	Tryptic Soy Agar	20/pack	IVD	
11038	Tryptic Soy Agar (Horse Blood 5%)	20/pack		40
11037	Tryptic Soy Agar (Sheep Blood 5%)	20/pack	IVD	

90 mm ready-to-use agar plates

Ref.	Description	Packaging	Intended use	Min. order qty (packs)
10006	Tryptic Soy Agar + 0,6% YEAST EXTRACT	20/pack	IVD	
10085	Tryptic Soy Agar + 1% Penase + Neutralizing	20/pack		20
10085	Tryptic Soy Agar + 1% Penase + Neutralizing	20/pack		20
11137	Tryptic Soy Agar + 5% Penase + Neutralizing	20/pack		20
10009	Tryptic Soy Agar + NaCl 3%	20/pack		20
10074	Tryptic Soy Agar+NEUTRALIZING	20/pack		
10432	Tryptone Soy Yeast Extract Agar	20/pack		
10449	Tryptose Sulfite Cycloserine Agar	20/pack		
10407	VANCOMYCIN SCREEN Agar	20/pack	IVD	
11184	VIOLET RED BILE GLUCOSE Agar	20/pack		
11183	VIOLET RED BILE LACTOSE Agar	20/pack		
10012	VIOLET RED BILE LACTOSE Agar + MUG	20/pack		20
11185	VOGEL JOHNSON Agar	20/pack		20
10053	WILKINS CHALGREN Agar	20/pack		20
10408	WILKINS CHALGREN Agar +5% SHEEP BLOOD	20/pack	IVD	
10054	WURTZ LACTOSE Agar	20/pack	IVD	
10056	X.L.D. Agar	20/pack	IVD	
10413	XLD Agar EP, USP, JP Formulation	20/pack	IVD	
10069	XLT 4 Agar	20/pack	IVD	
10011	YEAST GLUCOSE CHLORAMPHENICOL Agar	20/pack		
10400	YEAST MOLD Agar	20/pack		40
10052	YERSINIA SELECTIVE Agar	20/pack	IVD	
10430	YSG Agar	20/pack		40

Irradiated 90 mm ready-to-use agar plates

Ref.	Description	Packaging	Min. order qty (packs)
10075S	Sabouraud Agar + Neutralizing (Irradiated)	20/pack	
10109S	Sabouraud CAF Agar + Penase + Cephalosporinase + Neutralizing (Irradiated)	20/pack	40
11035S	Sabouraud CAF Agar (Irradiated)	20/pack	20
10035S	Sabouraud Dextrose Agar (Irradiated)	20/pack	
10114S	Sabouraud Dextrose Agar (Irradiated) 30 mL	20/pack	
10088S	Sabouraud Dextrose Agar + 0.1% Penase + Neutralizing (Irradiated)	20/pack	
10108S	Sabouraud Dextrose Agar + Lactamator 500 IU + Neutralizing (Irradiated)	20/pack	
10086S	Tryptic Soy Agar + 0.1% Penase + Neutralizing (Irradiated)	20/pack	30
10111S	Tryptic Soy Agar + 0.1% Penase + Neutralizing (Irradiated) 30 mL	20/pack	20
10085S	Tryptic Soy Agar + 1% Penase + Neutralizing (Irradiated)	20/pack	
10103S	Tryptic Soy Agar + Lactamator 1000 IU (Irradiated)	20/pack	
10105S	Tryptic Soy Agar + Lactamator 1000 IU + Neutralizing (Irradiated)	20/pack	30
10112S	Tryptic Soy Agar + Lactamator 1000 IU + Neutralizing (Irradiated) 30 mL	20/pack	20
10102S	Tryptic Soy Agar + Lactamator 500 IU (Irradiated)	20/pack	30
10116S	Tryptic Soy Agar + Lactamator 500 IU (Irradiated) 30 mL	20/pack	20
10104S	Tryptic Soy Agar + Lactamator 500 IU + Neutralizing (Irradiated)	20/pack	
10110S	Tryptic Soy Agar + Lactamator 1000 IU + Neutralizing (Irradiated) 30 mL	20/pack	20
10106S	Tryptic Soy Agar + Penase + Cephalosporinase + Neutralizing (Irradiated)	20/pack	30
10037S	Tryptic Soy Agar (Irradiated)	20/pack	
10058S	Tryptic Soy Agar 30 mL (Irradiated)	20/pack	
10076S	Tryptic Soy Agar 18 mL (Irradiated)	20/pack	30
10059S	Tryptic Soy Agar + Neutralizing 30 mL (Irradiated)	20/pack	
10074S	Tryptic Soy Agar + Neutralizing (Irradiated)	20/pack	

Ready culture media in two sectors petri dishes		# Product with minimum order	
Description	Packaging	Ref.	
			
Baird Parker / MacConkey Selective media for <i>Staphylococcus aureus</i> and gram-negatives isolation.	20/pack	18500 #	
Baird Parker / Sabouraud CAF Selective media for <i>Staphylococcus aureus</i> , yeasts and moulds isolation.	20/pack	18390 #	
Brilliant Green / S.S. Selective media for <i>Salmonella</i> spp. isolation.	20/pack	18012 #	
CLED / MacConkey Media for total bacterial count and gram-negatives isolation in urine.	20/pack	18502 #	
Columbia C.N.A. / Chocolate Selective medium for staphylococci, streptococci, pneumococci and non selective medium for fastidious bacteria isolation.	20/pack	18507 #	
Columbia C.N.A. / Gardnerella V. Selective media for streptococci, staphylococci, pneumococci and <i>Gardnerella vaginalis</i> isolation.	20/pack	18422 #	
Columbia C.N.A. / MacConkey Selective media for streptococci, staphylococci, pneumococci and gram-negative isolation.	20/pack	18327 #	
D.T.M. / Sabouraud Selective media for dermatophytes, yeasts, moulds isolation.	20/pack	18595	
Gardnerella V. / Thayer Martin Selective media for <i>Gardnerella vaginalis</i> and pathogenic neisseriae isolation.	20/pack	18379 #	
Hektoen Enteric / Yersinia Selective media for enteric pathogens and <i>Yersinia enterocolitica</i> isolation.	20/pack	18391 #	
Hektoen Enteric / S.S. Selective media for enteric pathogens isolation.	20/pack	18503 #	
MacConkey / Sabouraud CAF Selective media for coliforms and pathogenic fungi detection.	20/pack	18018 #	
MacConkey / S.S. Selective media for enteric pathogens isolation.	20/pack	18505 #	
MacConkey / T.S.A. Blood Media for Gram-negative and fastidious bacteria isolation.	20/pack	18380 #	
Schaedler K / Schaedler KKV Bi-plate for detection of anaerobic organisms and selective isolation of Gram-negative anaerobic bacilli from clinical specimens.	20/pack	18025	
T.S.A. Blood / Chromatic Detection Medium for Gram-negative fastidious bacteria isolation and urinary bacterial identification.	20/pack	18008 #	

140 mm ready-to-use agar plates				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
10224	BAIRD PARKER Agar	10/pack	IVD	
10250	Brilliant Green Agar ISO 6579 : 2017	10/pack		10
10248	Chromatic Bacillus Cereus	10/pack		10
10246	Chromatic MH	10/pack	IVD	10
11132	Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L beta-NAD) (140mm)	10/pack	IVD	10
10231	Mueller Hinton II Agar	10/pack	IVD	
12031	Mueller Hinton II Agar 120 mm x 120 mm	5/pack	IVD	
12032	Mueller Hinton II Agar (Sheep Blood 5%) 120 mm x 120 mm	5/pack	IVD	
10249	Purple Lactose Agar	10/pack	IVD	
10233	RPMI Agar	10/pack	IVD	
10234	Rose Bengal CAF Agar	10/pack		20
10226	Tryptic Soy Agar + Neutralizing	10/pack		10
10223	O.A. Listeria Agar	10/pack		

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
24154	ACETAMIDE Broth 5 mL	20/pack		
414060	Albumin Solution 200 mL	6/pack		10
401990	ALKALINE PEPTONE WATER 100 mL	6/pack	IVD	
24100	ALKALINE PEPTONE WATER	20/pack	IVD	
451404	Alkaline Peptone Water (screw cap) 225 mL	25/pack	IVD	5
414130	Alkaline Saline Peptone Water ISO 225 mL	6/pack		25
26486	Alkaline Saline Peptone Water ISO 10 mL	100/pack		
26477	Amos Kont Jones Medium 10 mL	100/pack		
20114	AZIDE DEXTROSE Broth	10/pack		
24114	AZIDE DEXTROSE Broth	20/pack		5
26114	AZIDE DEXTROSE Broth	100/pack		4
21114	AZIDE DEXTROSE Broth (D.C.)	10/pack		
24414	AZIDE DEXTROSE Broth (D.C.)	20/pack		
26414	AZIDE DEXTROSE Broth (D.C.)	100/pack		4
402710	BACILLUS CEREUS Agar Base 100 mL	6/pack		6
420110	Baird Parker Agar Base 100 mL	6/pack		
452680	Baird Parker Agar Base 200 mL	25/pack		
420010	Baird Parker Agar+RPF Suppl. 90 mL	6/pack		
402670	BAT Broth 90 mL	6/pack		20
30091	BIGGY (NICKERSON) Agar	10/pack	IVD	20
442350	BIGGY (NICKERSON) Agar 100 mL	6/pack	IVD	
33085	BILE AESCULIN Agar	10/pack	IVD	20
442320	BILE AESCULIN Agar 100 mL	6/pack	IVD	
24120	BILE AESCULIN Broth	20/pack	IVD	
403800	Blaurock semi-liquid modified hepatic Medium 25 mL (capacity 125 mL)	6/pack		40
470340	Bolton Broth 500 mL	6/pack		
400260	Bolton Broth 1000 mL	6/pack		6
414150	Bolton Broth 90 mL	6/pack		
414080	Bolton Broth (screw cap) 225 mL (capacity 250 mL)	6/pack		
451405	Bolton Broth (screw cap) 225 mL (capacity 250 mL)	25/pack		
30080	BORDET GENGOU Agar (Sheep Blood 15%)	10/pack	IVD	20
30084	Brain Heart Infusion Agar	10/pack	IVD	
402380	Brain Heart Infusion Agar 100 mL	6/pack	IVD	12
24480	Brain Heart Infusion Broth 9 mL	20/pack	IVD	
24104	Brain Heart Infusion Broth	20/pack	IVD	
24141	Brain Heart Infusion Broth 2 mL	20/pack	IVD	10
412010	Brain Heart Infusion Broth 200 mL	6/pack	IVD	6
27502	Brain Heart Infusion Broth 5 mL	50/pack	IVD	4
26104	Brain Heart Infusion Broth 10 mL	100/pack	IVD	
402330	BRILLIANT GREEN Agar 100 mL	6/pack	IVD	6
24102	BRILLIANT GREEN BILE Broth 2%	20/pack		
402560	BRILLIANT GREEN BILE Broth 2% 100 mL	6/pack		10
24418	BRUCELLA Broth	20/pack		10
401890	BUFFER SOLUTION pH 7 100 mL	6/pack		20
26470	Buffered NaCl Peptone Solution modified pH 7.0 + NT	100/pack		
452650	Buffered NaCl Peptone Solution modified pH 7.0 + NT 90 mL	25/pack		3
402410	Buffered NaCl Peptone Solution pH 7.0 100 mL	6/pack		
403410	Buffered NaCl Peptone Solution pH 7.0 300 mL (capacity 500 mL)	6/pack		20
402730	Buffered NaCl Peptone Solution pH 7.0 90 mL (capacity 125 mL)	6/pack		
402620	Buffered NaCl Peptone Solution pH 7.0 + N 100 mL	6/pack		
402740	Buffered NaCl Peptone Solution pH 7.0 + N 90 mL (capacity 125 mL)	6/pack		20
400310	Buffered NaCl Peptone Solution pH 7.0 + N (flip-off cap) 1000 mL	6/pack		10
402630	Buffered NaCl Peptone Solution pH 7.0 + NT 90 mL	6/pack		20
24099	Buffered Peptone Water	20/pack		
24199	Buffered Peptone Water 9 mL	20/pack		
414030	Buffered Peptone Water 90 mL	6/pack		
412090	Buffered Peptone Water 200 mL	6/pack		
414020	Buffered Peptone Water 225 mL	6/pack		
454030	Buffered Peptone Water 90 mL	25/pack		5
451402	Buffered Peptone Water 225 mL	25/pack		
24463	Buffered Peptone Water (Double Concentration)	20/pack		

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
414170	Buffered Peptone Water double concentration 225 mL	6/pack		
400040	Buffered Peptone Water EP, USP 100 mL	6/pack		6
400140	Buffered Peptone Water EP, USP 300 mL (capacity 500 mL)	6/pack		10
400240	Buffered Peptone Water EP, USP 1000 mL	6/pack		
413030	Campylobacter Agar 200 mL	6/pack	IVD	6
24404	Campylobacter Broth	20/pack	IVD	
470290	CARY BLAIR TRANSPORT MEDIUM 500 mL	6/pack	IVD	
24459	CASO Broth (Double Concentration)	20/pack		10
402270	CETRIMIDE Agar 100 mL	6/pack	IVD	
412270	CETRIMIDE Agar 200 mL	6/pack	IVD	
29000	CHECK-SET Broth Irradiated	20/pack	IVD	
30099	Chocolate Agar	10/pack	IVD	
442220	Chocolate Agar 100 mL	6/pack	IVD	
470120	Chocolate Agar 500 mL	6/pack	IVD	4
30081	CHRISTENSEN UREA Agar 10 mL	10/pack	IVD	
412100	CHRISTENSEN UREA Agar 200 mL	6/pack	IVD	
481110	Chromatic Candida 100 mL	6/pack	IVD	
481120	Chromatic Coli Coliform 100 mL	6/pack		
453100	Chromatic Coli Coliform 200 mL	25/pack		
481190	Chromatic Coliform Agar ISO 100 mL	6/pack		
453110	Chromatic Coliform Agar ISO 200 mL	25/pack		
481130	Chromatic Detection 100 mL	6/pack	IVD	
482190	Chromatic E.COLI O157 200 mL	6/pack		6
481140	Chromatic Salmonella 100 mL	6/pack	IVD	
481160	Chromatic STAPH AUREUS 100 mL	6/pack		3
402180	CLED Agar 100 mL	6/pack	IVD	
412180	CLED Agar 200 mL	6/pack	IVD	
470110	CLED Agar 500 mL	6/pack	IVD	
452210	Columbia Agar BASE 200 mL	6/pack	IVD	
24071	COOKED MEAT MEDIUM	20/pack		20
24483	Cronobacter Selective Broth	20/pack		
403210	D/E Neutralizing Broth 100 mL	6/pack		34
412460	D/E Neutralizing Broth 200 mL	6/pack		17
424010	DEMI FRASER Broth 90 mL	6/pack		6
414000	DEMI FRASER Broth 225 mL	6/pack		
452401	DEMI FRASER Broth 90 mL	25/pack		5
451400	DEMI FRASER Broth 225 mL	25/pack		
33086	DERMATHOPHYTE (D.T.M.) Agar	10/pack	IVD	
402200	DERMATHOPHYTE (D.T.M.) Agar 100 mL	6/pack	IVD	
402340	DESOXYCHOLATE Agar 100 mL	6/pack	IVD	
26487	Dextrose Tryptone Broth 9 mL	100/pack		
20076	DEY ENGLELY NEUTRALIZING B	20/pack		
412430	Dichloran Glycerol (DG18) Agar 200 mL	6/pack		
452430	Dichloran Glycerol (DG18) Agar 200 mL	25/pack		
412410	Dichloran Rose Bengal CAF Agar 200 mL	6/pack		
452410	Dichloran Rose Bengal CAF Agar 200 mL	25/pack		
21097	Distilled Water 9 mL	10/pack		
20097	Distilled Water 10 mL	20/pack		
473030	Distilled Water 500 mL	6/pack		
402220	Drigalski Lactose Agar 100 mL	6/pack		20
20122	E.C. Broth	10/pack		
24122	E.C. Broth	20/pack		
26122	E.C. Broth	100/pack		
402350	E.M.B. Levine Agar 100 mL	6/pack	IVD	
24101	E.V.A. Broth (ETHYL VIOLET AZIDE)	20/pack		15
24096	EE Broth-Mossel	25/pack		
453080	EE Broth-Mossel (flip-off cap) 100 mL	25/pack		3
402480	EE Broth-Mossel (screw cap) 100 mL	6/pack		
452480	EE Broth-Mossel (screw cap) 100 mL	25/pack		
403010	Endo Agar 100 mL	6/pack		
26492	Eugon LT SUP 9 mL	100/pack		
402750	Eugon LT SUP 90 mL (capacity 125 mL)	6/pack		

Ready-to-use culture media in tubes and bottles

Ref.	Description	Packaging	IVD	Min. order qty (packs)
495030	Fluid A 100 mL	6/pack		5
400210	Fluid A 1000 mL	6/pack		10
459503	Fluid A 100 mL	25/pack		5
495170	Fluid A (Perforable) 500 mL (capacity 1000 mL)	6/pack		10
495150	Fluid A (Perforable) 1000 mL	6/pack		10
400210S	Fluid A Irradiated 1000 mL	6/pack		10
400110	Fluid A (flip-off cap) 300 mL (capacity 500 mL)	6/pack		
453010	Fluid A (flip-off cap) 100 mL	25/pack		3
495040	Fluid D 100 mL	6/pack		5
459504	Fluid D 100 mL	25/pack		
495190	Fluid D (Perforable) 500 mL	6/pack		20
495200	Fluid D (Perforable) 1000 mL	6/pack		10
400170	Fluid D (flip-off cap) 300 mL (capacity 500 mL)	6/pack		10
453070	Fluid D (flip-off cap) 100 mL	25/pack		3
495160	Fluid K (Perforable) 900 mL (capacity 1000 mL)	6/pack		10
21241	Fluid Thioglycollate Medium	10/pack	IVD	
24124	Fluid Thioglycollate Medium	20/pack	IVD	
24241	Fluid Thioglycollate Medium 20 mL	20/pack		10
26124	Fluid Thioglycollate Medium	100/pack		
495020	Fluid Thioglycollate Medium (crimp cap) 100 mL	6/pack		10
400020	Fluid Thioglycollate Medium (flip-off cap) 100 mL	6/pack		10
453020	Fluid Thioglycollate Medium (flip-off cap) 100 mL	25/pack		
493000	Fluid Thioglycollate Medium (perforable cap) 100 mL	6/pack		10
452060	Fluid Thioglycollate Medium (screw cap) 100 mL	6/pack		
463100	Fluid Thioglycollate Medium (screw cap) 900 mL (capacity 1000 mL)	6/pack		6
453060	Fluid Thioglycollate Medium (screw cap) 100 mL	25/pack		
470300	Fluid Thioglycollate Medium (wide neck) 500 mL	6/pack		
402500	Fluid Thioglycollate Medium + 1% Tween 80 100 mL	6/pack		10
452500	Fluid Thioglycollate Medium + 1% Tween 80 100 mL	25/pack		3
453280	Fluid Thioglycollate Medium + 1% Tween 80 (flip-off cap) 100 mL	25/pack		
34094	Gelatin Peptone Agar	10/pack		
401960	Gelisato Agar 150 mL	6/pack		6
20105	Glucose Broth	10/pack	IVD	
24105	Glucose Broth	20/pack		5
26105	Glucose Broth	100/pack	IVD	
24482	Glucose OF Medium ISO 21528	20/pack		
26473	Glucose Yeast Extract Agar 20 mL	100/pack		
495130	Glycine 3% (Perforable) 1000 mL	6/pack		
414070	GN (Hajna) Broth 225 mL	6/pack	IVD	
24119	GN HAJNA Broth	20/pack	IVD	
24091	Haemophilus Test Broth	20/pack	IVD	
402230	HEKTOEN ENTERIC Agar 100 mL	6/pack	IVD	
412230	HEKTOEN ENTERIC Agar 200 mL	6/pack	IVD	
20090	Helicobacter pylori Test	20/pack	IVD	
490010	HEMO-AEROBIC CULTURING 80 mL	6/pack	IVD	
490050	HEMO-AEROBIC CULTURING Neonatal 9 mL	6/pack	IVD	
490030	HEMO-AEROBIC CULTURING-Pediatric 40 mL	6/pack	IVD	
490020	HEMO-ANAEROBIC CULTURING 80 mL	6/pack	IVD	
490060	HEMO-ANAEROBIC CULTURING Neonatal 9 mL	6/pack	IVD	
490040	HEMO-ANAEROBIC CULTURING Pediatric 40 mL	6/pack	IVD	
				
33118	I.U.T.M. MEDIUM	10/pack		80
20121	INOCULUM Broth 7 ML	20/pack	IVD	
20081	INOCULUM SOLUTION 5 ML	20/pack	IVD	
403180	Iron Sulphite Agar 100 mL	6/pack		

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
412350	Kanamycin Aesculin Azide Agar 200 mL	6/pack		
31203	King's A Medium 5 mL	20/pack		
31202	King's B Medium 5 mL	20/pack		
25314	KIT STREPTO FECALI NELLE ACQUE	5 tubes x 10 mL 1 bottle x 50 mL		20
30087	KLIGLER Iron Agar	10/pack	IVD	
30089	KLIGLER Iron Agar (Butt=5cmm/Slope=6cm)	10/pack		20
412390	L.B. Medium 200 mL	6/pack		6
20106	Lactose Broth	10/pack		
402090	Lactose Broth 100 mL	6/pack		
21106	Lactose Broth (Double Concentration)	10/pack		
24468	Lactose Sulfite Broth	20/pack		
24481	Lauryl Sulphate Tryptose Broth 9 mL	20/pack		
21453	Lauryl Sulphate Tryptose Broth (LST)	10/pack		
26453	Lauryl Sulphate Tryptose Broth (LST)	100/pack		
21454	Lauryl Sulphate Tryptose Broth (LST) DOUBLE	10/pack		
26454	Lauryl Sulphate Tryptose Broth (LST) DOUBLE	100/pack		4
24457	Lauryl Sulphate Tryptose Broth (LST) w/ MUG 9 mL	20/pack		10
26457	Lauryl Sulphate Tryptose Broth (LST) w/ MUG 9 mL	100/pack		4
21455	Lauryl Sulphate Tryptose Broth (LST) w/ MUG	10/pack		40
26455	Lauryl Sulphate Tryptose Broth (LST) w/ MUG	100/pack		4
402600	Legionella BCYE Agar Base 90 mL	6/pack		
24170	Lothian Broth 9 mL	20/pack		15
414120	Lethen Modified Broth 225 ml	6/pack		
85701	LIM Broth 5 mL	100/pack		
414090	Listeria Enrichment Broth 225 mL	6/pack		
20131/9	Listeria Fraser Broth 9 mL	10/pack		30
24131	Listeria Fraser Broth 10 mL	20/pack		
413980	Listeria Fraser Broth 225 ml	6/pack		
24471	Listeria Motility Medium	20/pack	IVD	30
33087	Listeria Palcam Agar	10/pack		10
24435	Listeria U.V.M. Broth	20/pack		5
461030	Listeria UVM 1 Broth 225 ml	6/pack		
461020	Listeria UVM 2 Broth 225 ml	6/pack		6
30116	Loeffler Medium	10/pack	IVD	
34127/1	Lowenstein Jensen + Amikacin 40ug/ml	10/pack	IVD	
34127	Lowenstein Jensen + Amikacin 5ug/ml	10/pack	IVD	
34138/1	Lowenstein Jensen + Capreomycin 10ug/ml	10/pack	IVD	
34138/3	Lowenstein Jensen + Capreomycin 20ug/ml	10/pack	IVD	
34138/4	Lowenstein Jensen + Capreomycin 30ug/ml	10/pack	IVD	
35090	Lowenstein Jensen + Capreomycin 30ug/ml	10/pack	IVD	
34138/2	Lowenstein Jensen + Capreomycin 40ug/ml	10/pack	IVD	
34131/2	Lowenstein Jensen + Clarithromycin 32ug/ml CE	10/pack	IVD	
34131/1	Lowenstein Jensen + Clarithromycin 4ug/ml	10/pack	IVD	
34139/2	Lowenstein Jensen + Clofazimine 10ug/ml	10/pack	IVD	
34139/1	Lowenstein Jensen + Clofazimine 5ug/ml	10/pack	IVD	
34137/2	Lowenstein Jensen + Cycloserine 10ug/ml	10/pack	IVD	
34137/3	Lowenstein Jensen + Cycloserine 20ug/ml	10/pack	IVD	
34137/1	Lowenstein Jensen + Cycloserine 30ug/ml	10/pack	IVD	
34137/4	Lowenstein Jensen + Cycloserine 40ug/ml	10/pack	IVD	
34137/5	Lowenstein Jensen + Cycloserine 50ug/ml	10/pack	IVD	
34126/6	Lowenstein Jensen + Ethambutol 10ug/ml	10/pack	IVD	
34126/4	Lowenstein Jensen + Ethambutol 1ug/ml	10/pack	IVD	
34126/1	Lowenstein Jensen + Ethambutol 2ug/ml	10/pack	IVD	
35030	Lowenstein Jensen + Ethambutol 2ug/ml	10/pack	IVD	
34126/5	Lowenstein Jensen + Ethambutol 3ug/ml	10/pack	IVD	
34126/2	Lowenstein Jensen + Ethambutol 4ug/ml	10/pack	IVD	
34126/3	Lowenstein Jensen + Ethambutol 5ug/ml	10/pack	IVD	
34132/1	Lowenstein Jensen + Ethionamide 10ug/ml	10/pack	IVD	
34132/2	Lowenstein Jensen + Ethionamide 20ug/ml	10/pack	IVD	
35040	Lowenstein Jensen + Ethionamide 20ug/ml	10/pack	IVD	
34132/3	Lowenstein Jensen + Ethionamide 30ug/ml	10/pack	IVD	

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
35041	Lowenstein Jensen + Ethionamide 30ug/ml	10/pack	IVD	
34132/4	Lowenstein Jensen + Ethionamide 40ug/ml	10/pack	IVD	
34123	Lowenstein Jensen + Isoniazid 0.1ug/ml	10/pack	IVD	
35001	Lowenstein Jensen + Isoniazid 0.20ug/ml	10/pack	IVD	
34123/1	Lowenstein Jensen + Isoniazid 0.2ug/ml	10/pack	IVD	
35002	Lowenstein Jensen + Isoniazid 1 ug/ml	10/pack	IVD	
34123/4	Lowenstein Jensen + Isoniazid 10ug/ml	10/pack	IVD	
34123/2	Lowenstein Jensen + Isoniazid 1ug/ml	10/pack	IVD	
34123/3	Lowenstein Jensen + Isoniazid 5ug/ml	10/pack	IVD	
34143/1	Lowenstein Jensen + Kanamycin 10ug/ml	10/pack	IVD	
34143/2	Lowenstein Jensen + Kanamycin 20ug/ml	10/pack	IVD	
35060	Lowenstein Jensen + Kanamycin 20ug/ml	10/pack	IVD	
34143/3	Lowenstein Jensen + Kanamycin 30ug/ml	10/pack	IVD	
35061	Lowenstein Jensen + Kanamycin 30ug/ml	10/pack	IVD	
34146/1	Lowenstein Jensen + Levofloxacin 2ug/ml	10/pack	IVD	
34135/1	Lowenstein Jensen + Nicotinamide 10ug/ml	10/pack	IVD	
34135/2	Lowenstein Jensen + Nicotinamide 20 ug/ml	10/pack	IVD	
34135/3	Lowenstein Jensen + NICOTINAMIDE 30ug/ml	10/pack	IVD	
34128/2	Lowenstein Jensen + Ofloxacin 10ug/ml	10/pack	IVD	
34128/5	Lowenstein Jensen + Ofloxacin 20ug/ml	10/pack	IVD	
34128/3	Lowenstein Jensen + Ofloxacin 25ug/ml	10/pack	IVD	
34128/4	Lowenstein Jensen + Ofloxacin 2ug/ml	10/pack	IVD	
35080	Lowenstein Jensen + Ofloxacin 2ug/ml	10/pack	IVD	
34128/1	Lowenstein Jensen + Ofloxacin 5ug/ml	10/pack	IVD	
34145	Lowenstein Jensen + PACT	10/pack	IVD	
34129/4	Lowenstein Jensen + PAS 0.1ug/ml	10/pack	IVD	
34129/3	Lowenstein Jensen + PAS 0.5ug/ml	10/pack	IVD	
34129/2	Lowenstein Jensen + PAS 10ug/ml	10/pack	IVD	
34129/1	Lowenstein Jensen + PAS 1ug/ml	10/pack	IVD	
35070	Lowenstein Jensen + PAS 1ug/ml	10/pack	IVD	
34129/5	Lowenstein Jensen + PAS 5ug/ml	10/pack	IVD	
34136	Lowenstein Jensen + Pefloxacin 2ug/ml	10/pack	IVD	
35147	Lowenstein Jensen + PNB 500 ug/ml	10/pack	IVD	
34124/2	Lowenstein Jensen + Pyrazinamide 15ug/ml	10/pack	IVD	
35050	Lowenstein Jensen + Pyrazinamide 1ug/ml	10/pack	IVD	
34124/4	Lowenstein Jensen + Pyrazinamide 200ug/ml	10/pack	IVD	
34124/3	Lowenstein Jensen + Pyrazinamide 20ug/ml	10/pack	IVD	
34124/1	Lowenstein Jensen + Pyrazinamide 5ug/ml	10/pack	IVD	
34144	Lowenstein Jensen + Pyruvate 0.2%	10/pack	IVD	
34130/1	Lowenstein Jensen + Rifabutin 10ug/ml	10/pack	IVD	
34130/2	Lowenstein Jensen + Rifabutin 30ug/ml	10/pack	IVD	
34130/3	Lowenstein Jensen + Rifabutin 50ug/ml	10/pack	IVD	
34121/2	Lowenstein Jensen + Rifampicin 10ug/ml	10/pack	IVD	
34121	Lowenstein Jensen + Rifampicin 15ug/ml	10/pack	IVD	
34121/6	Lowenstein Jensen + Rifampicin 20ug/ml	10/pack	IVD	
34121/3	Lowenstein Jensen + Rifampicin 25ug/ml	10/pack	IVD	
34121/5	Lowenstein Jensen + Rifampicin 40ug/ml	10/pack	IVD	
35010	Lowenstein Jensen + Rifampicin 40ug/ml	10/pack	IVD	
34121/4	Lowenstein Jensen + Rifampicin 50ug/ml	10/pack	IVD	
34121/1	Lowenstein Jensen + Rifampicin 5ug/ml	10/pack	IVD	
34122	Lowenstein Jensen + Rifapentin 9ug/ml	10/pack	IVD	
34125/2	Lowenstein Jensen + Streptomycin 10ug/ml	10/pack	IVD	
35021	Lowenstein Jensen + Streptomycin 10ug/ml	10/pack	IVD	
34125/3	Lowenstein Jensen + Streptomycin 25ug/ml	10/pack	IVD	
34125/4	Lowenstein Jensen + Streptomycin 2ug/ml	10/pack	IVD	
34125/1	Lowenstein Jensen + Streptomycin 4ug/ml	10/pack	IVD	
35020	Lowenstein Jensen + Streptomycin 4ug/ml	10/pack	IVD	
34125/5	Lowenstein Jensen + Streptomycin 50ug/ml	10/pack	IVD	
35148	Lowenstein Jensen + TCH 2 ug/ml	10/pack	IVD	
30118	Lowenstein Jensen Medium 9 mL	10/pack	IVD	
35000	Lowenstein Jensen Medium 8 mL	10/pack	IVD	
31118	Lowenstein Jensen Medium 9 mL	20/pack	IVD	

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
30119	Lowenstein Jensen Medium w/o Glycerol	10/pack	IVD	
403600	LPT Diluent Broth 90 mL (capacity 125 mL)	6/pack		
24139	LYSINE DECARBOXYLASE Broth	20/pack	IVD	
30098	LYSINE IRON Agar	10/pack	IVD	
412040	LYSINE IRON Agar 200 mL	6/pack		6
30040	M 17 Agar	10/pack		40
402110	M 17 Agar 100 mL	6/pack		
452340	M 17 Agar 200 mL	25/pack		
403020	m-FC Agar + ROSOLIC ACID 100 mL	6/pack		
402120	M.R.S. Agar 100 mL	6/pack	IVD	
412310	M.S.R.V. Agar 200 mL	6/pack		
452240	MacConkey Agar 200 mL	25/pack	IVD	
402240	MacConkey Agar 100 mL	6/pack	IVD	
412240	MacConkey Agar 200 mL	6/pack	IVD	
470090	MacConkey Agar (Harm EP) 500 mL	6/pack	IVD	
20126	MacConkey Broth	6/pack		
402490	MacConkey Broth 100 mL	6/pack		
24126	MacConkey Broth	20/pack		
452490	MacConkey Broth 100 mL	25/pack		
494000	MacConkey Broth (Perforable) 100 mL	6/pack		
453090	MacConkey Broth (flip-off cap) 100 mL	25/pack		3
24226	MacConkey Broth singled strength 5 mL	20/pack		10
402470	Malt Extract Agar 100 mL	6/pack		
402290	Mannitol Salt Agar 100 mL	6/pack	IVD	
412290	Mannitol Salt Agar 200 mL	6/pack	IVD	
470080	Mannitol Salt Agar 500 mL	6/pack	IVD	
20071	Maximum Recovery Diluent 9 mL	20/pack		
402590	Maximum Recovery Diluent 100 mL	6/pack		
402660	Maximum Recovery Diluent 90 mL	6/pack		
412400	Maximum Recovery Diluent 200 mL	6/pack		
412420	Maximum Recovery Diluent 225 mL	6/pack		
26071	Maximum Recovery Diluent 9 mL	100/pack		
452420	Maximum Recovery Diluent 225 mL (capacity 250 mL)	25/pack		
80356	MeRA Test	50 tests		
30368	Middlebrook 7H10 Agar	10/pack	IVD	10
37000	Middlebrook 7H11	10/pack	IVD	
37001	Middlebrook 7H11 + AMIKACIN 2ug/ml	10/pack	IVD	30
37002	Middlebrook 7H11 + AMIKACIN 4ug/ml	10/pack	IVD	30
37056	Middlebrook 7H11 + CYCLOSERINE 30ug/ml	10/pack	IVD	30
37006	Middlebrook 7H11 + ETHAMBUTOL 7.5ug/ml	10/pack	IVD	30
37011	Middlebrook 7H11 + ETHIONAMIDE 10ug/ml	10/pack	IVD	30
37016	Middlebrook 7H11 + ISONIAZIDE 0.2ug/ml	10/pack	IVD	30
37017	Middlebrook 7H11 + ISONIAZIDE 1ug/ml	10/pack	IVD	30
37051	Middlebrook 7H11 + OFLOXACIN 2ug/ml	10/pack	IVD	30
37026	Middlebrook 7H11 + PAS 8ug/ml	10/pack	IVD	30
37031	Middlebrook 7H11 + PYRAZINAMIDE 25ug/ml	10/pack	IVD	30
37036	Middlebrook 7H11 + RIFABUTIN 1ug/ml	10/pack	IVD	30
37041	Middlebrook 7H11 + RIFAMPICIN 1ug/ml	10/pack	IVD	30
37046	Middlebrook 7H11 + STREPTOMYCIN 2ug/ml	10/pack	IVD	30
24436	Middlebrook 7H9 Broth	20/pack	IVD	
401950	MILK Agar 150 mL	6/pack		
463120	Milk Plate Count Agar 500 mL	6/pack		
452660	Milk Plate Count Agar 200 mL	25/pack		
24466	Mineral Modified Glutamate Medium	20/pack		
26466	Mineral Modified Glutamate Medium	100/pack		4
24467	Mineral Modified Glutamate Medium (Double Concentration)	20/pack		
26467	Mineral Modified Glutamate Medium (Double Concentration)	100/pack		4
80355	MiRA Test	50 tests		
31204	MIU Agar	20/pack	IVD	
24343	MIU Semisolid Agar	20/pack		10
24452	MODIFIED LAURYL SULFATE Tryptose Broth	20/pack		
24489	MOPS Buffered Listeria Enrichment Broth	20/pack		

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
24413	MOSSEL AND MARTIN w MANNITOL	20/pack		15
24472	Motility Nitrate CP Medium	20/pack		
24342	MOTILITY TEST MEDIUM	20/pack	IVD	
26342	Motility Test Medium	100/pack	IVD	
403700	MPC-1 Medium 25 mL (capacity 125 mL)	6/pack		40
24149	MR-VP MEDIUM	20/pack	IVD	
495140	MRD + Tween + Lecithin modified (Perforable) 1000 mL	6/pack		
412320	MRS Agar 200 mL	6/pack	IVD	
452320	MRS Agar 200 mL	25/pack	IVD	
402160	MRS Agar pH 5 100 mL	6/pack		
470380	MSRV Agar (Semisolid) 500 mL	6/pack		
26478	MSRV Agar (semisolid) 20 mL	100/pack		
432090	MTSB with Novobiocin 225 mL	6/pack		6
402250	Mueller Hinton II Agar 100 mL	6/pack	IVD	
412250	Mueller Hinton II Agar 200 mL	6/pack	IVD	
470070	Mueller Hinton II Agar 500 mL	6/pack	IVD	
24107	Mueller Hinton II Broth	20/pack	IVD	
402020	Mueller Hinton II Broth 100 mL	6/pack	IVD	
463240	Mueller Hinton II Broth 500 mL	6/pack		10
23002	Mueller Hinton II Broth w/ Horse Blood 11 mL	10/pack	IVD	30
402030	Muller Kauffmann Broth 100 mL	6/pack		6
24108	Muller Kaufmann Broth	20/pack	IVD	
20072	Muller Kauffmann TETRATHIONATE NOV.B_MKTTn	20/pack		
26072	Muller Kauffmann TETRATHIONATE NOV.B_MKTTn	100/pack		
403110	Muller Kauffmann TETRATHIONATE NOV.B_MKTTn 100 mL	6/pack		
20162	Mycoplasma Selective Broth 2 mL	20/pack	IVD	
20158	Mycoplasma TRANSPORT Broth 6 mL	20/pack	IVD	
80100	NEUTRALIZING SOLUTION 100 mL	6/pack		
413050	Neutralizing w/ Saponin 200 mL	6/pack		6
463170	Neutralizing w/ Saponin 250 mL (capacity 500 mL)	6/pack		6
463210	Neutralizing w/o Histidine 250 mL (capacity 500 mL)	6/pack		6
30183	Nutrient Agar 10 ml	10/pack	IVD	
30083	Nutrient Agar ISO 16266	10/pack		30
402190	Nutrient Agar ISO 16266 100 mL	6/pack		
412190	Nutrient Agar ISO 16266 200 mL	6/pack		
470060	Nutrient Agar ISO 16266 500 mL	6/pack		6
30103	Nutrient Agar ISO 6579 9 mL	20/pack		
413080	Nutrient Agar ISO 6579 200 mL	6/pack		
402000	Nutrient Broth 100 mL	6/pack	IVD	
24103	Nutrient Broth 10 mL	20/pack	IVD	
27503	Nutrient Broth 5 mL	50/pack	IVD	
26103	Nutrient Broth 10 mL	100/pack	IVD	
24153	Nutrient GELATIN	20/pack		10
24345	OF Medium w/ Glucose	20/pack	IVD	
403170	OGYE Agar 100 mL	6/pack		
452690	Orange Serum Agar 200 mL	25/pack		3
401870	Peptone NEUTRALIZING SOLUTION 45 mL	6/pack		100
463160	Peptone Neutralizing Solution 250 mL (capacity 500 mL)	6/pack		20
402130	Peptone Water 100 mL	6/pack	IVD	
24098	Peptone Water	20/pack	IVD	
452640	Peptone Water 90 mL	25/pack		
470320	Peptone Water 500 mL	6/pack		10
402530	Peptone Water (Perforable) 100 mL	6/pack		
414010	Peptone Water pH 8.4 + NaCl 1% 225 mL	6/pack	IVD	
30117	PERGOLA MEDIUM	10/pack	IVD	
24446	PHENOL RED Broth	20/pack		
30085	PHENYLALANINE Agar	10/pack		50
412170	PHENYLALANINE Agar 200 mL	6/pack	IVD	
462040	PHOSPHATE BUFFER NON STERILE 250 mL	6/pack		6
463190	Phosphate Buffer Solution pH 7.2 1000 mL	6/pack		
20095	Physiological Solution 7 mL	20/pack	IVD	
20196	Physiological Solution 9 mL	20/pack	IVD	

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
20197	Physiological Solution 10 mL	20/pack	IVD	
20079	Physiological Solution 4.5 mL	20/pack	IVD	
471120	Physiological Solution 240 mL	6/pack	IVD	
473000	Physiological Solution 500 mL	6/pack	IVD	
24142	Physiological Solution 3 mL	20/pack	IVD	
470400	Physiological Solution 250 mL (capacity 500 mL)	6/pack	IVD	
463200	Physiological Solution 1000 mL	6/pack	IVD	
26196	Physiological Solution 9 mL	100/pack	IVD	
463250	Physiological solution (screw cap) 350 mL (capacity 500 mL)	6/pack		50
495220	Physiological Solution (Perforable) 500 mL	6/pack		5
495230	Physiological Solution (Perforable) 1000 mL	6/pack		5
471130	Physiological Solution 38 mL	6/pack		15
34073	Plate Count Agar 22 mL	10/pack		
33070	Plate Count Agar	10/pack		20
402260	Plate Count Agar 100 mL	6/pack		
401940	Plate Count Agar 150 mL	6/pack		
412260	Plate Count Agar 200 mL	6/pack		
470180	Plate Count Agar 500 mL	6/pack		
452260	Plate Count Agar 100 mL	25/pack		
462260	Plate Count Agar 200 mL	25/pack		
413100	Potato Dextrose Agar 200 mL	6/pack		
414160	Preston Broth 90 mL	6/pack		
24484	Preston Broth ISO 10272 10 mL	20/pack		15
412130	Pseudomonas Agar BASE 200 mL	6/pack	IVD	
24151	Purple Glucose Agar	20/pack		
31151	Purple Glucose Agar 15 mL	20/pack		
20140	Purple Lactose Broth	10/pack		20
24140	Purple Lactose Broth	20/pack		5
402580	R2A Agar 100 mL	6/pack		
472580	R2A Agar 500 mL	6/pack		
463110	R2A Agar 500 mL screw cap	6/pack		6
452700	R2A Agar 200 mL	25/pack		
80351	RAPID ANTIBIOTIC TEST	50 tests		
412470	Rapid Y/M Agar 200 mL	6/pack		
24450	RAPPAPORT Broth w/o SOY	20/pack	IVD	
26400	Rappaport Vassiliadis Soy (RVS) Broth	100/pack	IVD	
24400	Rappaport Vassiliadis Soy (RVS) Broth	20/pack	IVD	
402550	Rappaport Vassiliadis Soy (RVS) Broth 100 mL	6/pack	IVD	
24202	RESTORING Broth	20/pack		
24491	Ringer 1/4 Solution	20/pack		50
442480	Rose Bengal CAF Agar 100 mL	6/pack		
24461	RPMI Broth	20/pack	IVD	
24462	RPMI Broth Double strength	20/pack		40
24411	S.F. Broth	20/pack		20
403050	S.I.M. Medium 100 mL	6/pack	IVD	
402300	S.S. Agar 100 mL	6/pack	IVD	
412300	S.S. Agar 200 mL	6/pack	IVD	
24109	Sabouraud Broth	20/pack	IVD	
471070	Sabouraud Broth 500 mL	6/pack	IVD	4
30024	Sabouraud CAF + Actidione Agar	10/pack	IVD	
31024	Sabouraud CAF + Actidione Agar	20/pack	IVD	
402370	Sabouraud CAF Agar 100 mL	6/pack	IVD	
412370	Sabouraud CAF Agar 200 mL	6/pack	IVD	
31023	Sabouraud CAF Agar	20/pack	IVD	30
30093	Sabouraud Dextrose Agar	10/pack	IVD	
402280	Sabouraud dextrose Agar 100 mL	6/pack	IVD	
412280	Sabouraud Dextrose Agar 200 mL	6/pack	IVD	
470040	Sabouraud Dextrose Agar 500 mL	6/pack	IVD	
452280	Sabouraud Dextrose Agar 200 mL	25/pack	IVD	
402040	Sabouraud Dextrose Broth 100 mL	6/pack	IVD	10
452040	Sabouraud Dextrose Broth (screw cap) 100 mL	25/pack	IVD	
442280	Sabouraud Agar Modified 100 mL	6/pack		6

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
24135	Salmonella Differential Broth	20/pack		30
24430	Schaedler Broth	20/pack		10
24143	Selenite Broth 5 mL	20/pack		
24110	Selenite Broth	20/pack	IVD	
402050	Selenite Broth 100 mL	6/pack	IVD	
412050	Selenite Broth 200 mL	6/pack	IVD	5
470020	Selenite Broth 500 mL	6/pack	IVD	4
26110	Selenite Broth	100/pack	IVD	
463130	Selenite Broth 1000 mL	6/pack	IVD	
402390	Selenite Cystine Broth 100 mL	6/pack		
24510	Selenite Cystine Broth	20/pack		
24485	Selenite Cystine Broth ISO 6579 10 mL	20/pack		15
24479	SIM Medium 9 mL	20/pack	IVD	
26095	Sim Medium	100/pack	IVD	
30011	SIMMONS CITRATE Agar	10/pack	IVD	
412030	SIMMONS CITRATE Agar 200 mL	6/pack	IVD	
33065	SPS Agar	10/pack	IVD	20
442490	SPS Agar 100 mL	6/pack	IVD	
401930	SPS Agar 150 mL	6/pack	IVD	
20199	STERI TEST MEDIUM	20/pack		5
20200	STERI TEST VHP MEDIUM (1% Pyr) 4 mL	20/pack		
495210	Sterile Acetate Buffer 20 mM, pH 5.0 (Perforable) 100 mL (capacity 300 mL)	6/pack		100
400300	Sterile Distilled Water 100 mL	6/pack		
400290	Sterile Distilled Water 1000 mL	6/pack		
414190	Sterile Distilled Water 200 mL	6/pack		20
402150	STREPTOCOCCAL KF Agar 100 mL	6/pack		
24412	STREPTOCOCCUS Broth	20/pack	IVD	
455210	Sugar Free Agar 100 mL	25/pack		4
24440	SULFATE REDUCERS Broth	20/pack		50
20089	SUSPENSION Broth 5 mL	20/pack		50
20157	SWAB SOLUTION SKELETON	20/pack		
20160	SWAB SOLUTION SURFACES 5 mL	20/pack		
26488	Swarm Agar 10 mL	100/pack		
30096	T.S.I. Agar	10/pack	IVD	
442100	TAT Broth 90 mL (capacity 200 mL)	6/pack		
442160	TAT Broth 190mL (capacity 200 mL)	6/pack		3
463070	TAT Broth 490 mL (capacity 1000 mL)	6/pack		
463020	TAT Broth 990 mL (capacity 1000 mL)	6/pack		
495180	TAT Broth (Perforable) 500 mL (capacity 1000 mL)	6/pack		
481170	TBX Agar 100 mL	6/pack		
451170	TBX Agar 200 mL	25/pack		
30022	TCBS Agar	10/pack		10
403140	TCBS Agar 100 mL	6/pack	IVD	
24451	TETRATHIONATE Broth	20/pack	IVD	
33040	THAYER MARTIN Agar	10/pack		10
20171	Thioglycollate Medium w Vit.K1 & Hemin	10/pack	IVD	
24145	Todd Hewitt Broth w/ Colistin/Nalid.a.	20/pack	IVD	
24111	Todd Hewitt Broth 10 mL	20/pack	IVD	
412060	Todd Hewitt Broth 200 mL	6/pack	IVD	6
27501	Todd Hewitt Broth 5 mL	50 x 5 mL	IVD	4
25004	TOTAL COLIFORMS IN WATER	5 tubes x 10 mL 1 bottle x 50 mL		30
402520	TRIBUTYRIN Agar 100 mL	6/pack		10
24115	Trichomonas Broth 5 mL	20/pack	IVD	
30196	Triple Sugar Iron Agar 15 mL	10/pack		
30100	TRIPLE SUGAR IRON Agar	10/pack		
413120	TRIPLE SUGAR IRON Agar 200 mL	6/pack		6
30082	Tryptic Soy Agar	10/pack	IVD	
442290	Tryptic Soy Agar 100 mL	6/pack	IVD	
414110	Tryptic Soy Agar 225 mL	6/pack		25
432290	Tryptic Soy Agar 200 mL	6/pack	IVD	

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
470010	Tryptic Soy Agar 500 mL	6/pack	IVD	
26475	Tryptic Soy Agar 20 mL	100/pack x 20 mL	IVD	
452290	Tryptic Soy Agar 200 mL	25/pack	IVD	
24513	Tryptic Soy Broth 10 mL	20/pack	IVD	
24469	Tryptic Soy Broth 9 mL	20/pack	IVD	
27500	Tryptic Soy Broth 5 mL	50 x 5 mL	IVD	
463090	Tryptic Soy Broth 900 mL (capacity 1000 mL)	6/pack		6
453208	Tryptic Soy Broth 200 mL	25/pack		5
26513	Tryptic Soy Broth 10 mL	100/pack	IVD	
463270	Tryptic Soy Broth 1000 mL	6/pack		50
492000	Tryptic Soy Broth (Perforable) 100 mL	6/pack		10
495010	Tryptic Soy Broth (crimp cap) 100 mL	6/pack		10
442080	Tryptic Soy Broth (screw cap) 200 mL	6/pack	IVD	
432080	Tryptic Soy Broth (screw cap) 225 mL	6/pack	IVD	
470370	Tryptic Soy Broth (screw cap) 500 mL	6/pack		6
455208	Tryptic Soy Broth (screw cap) 100 mL	25/pack	IVD	
453030	Tryptic Soy Broth (flip-off cap) 100 mL	25/pack	IVD	
452080	Tryptic Soy Broth (screw cap) 100 mL	6/pack	IVD	
452080S	Tryptic Soy Broth (triple wrapped and gamma-irradiated) 100 mL	6/pack	IVD	
24460	Tryptic Soy Broth + 0.3% Lecithin + 2% Tween 80	20/pack		20
400270	Tryptic Soy Broth + 1% Tween 80 100 mL	6/pack		6
412110	Tryptic Soy Broth + 1% Tween 80 200 mL	6/pack		6
453270	Tryptic Soy Broth + 1% Tween 80 (flip-off cap) 100 mL	25/pack		3
452090	Tryptic Soy Broth + 1% Tween 80 (screw cap) 100 mL	25/pack		
27002	Tryptic Soy Broth + 15% Glycerol	20/pack		5
26002	Tryptic Soy Broth + 15% Glycerol	100/pack		
442090	Tryptic Soy Broth + Capitol 4 90 mL (capacity 200 mL)	6/pack		
442150	Tryptic Soy Broth + Capitol 4 190mL (capacity 200 mL)	6/pack		
463080	Tryptic Soy Broth + Capitol 4 490 mL (capacity 1000 mL)	6/pack		
463010	Tryptic Soy Broth + Capitol 4 990 mL (capacity 1000 mL)	6/pack		
403500	Tryptic Soy Broth + Lecithin + Tween 80 90 mL	6/pack		25
400230	Tryptic Soy Broth EP, USP 1000 mL	6/pack		5
400030	Tryptic Soy Broth EP, USP (flip-off cap) 100 mL	6/pack	IVD	
24445	Tryptic Soy Broth w 15% GLYCEROL	20/pack		10
21115	Tryptic Soy Broth w PHENOL RED	10/pack		10
26474	Tryptone Glucose Broth 10 mL	100/pack		10
412440	Tryptone Soy Yeast Extract Agar 200 mL	6/pack		
24464	Tryptone Soy Yeast Extract Broth	20/pack		
31201	Tryptone Sulfite Neomycin Agar 10 mL	20/pack		15
401980	Tryptone WATER 100 mL	6/pack	IVD	
24136	Tryptone WATER	20/pack	IVD	
24073	TRYPTOPHAN Broth	20/pack		
26073	TRYPTOPHAN Broth	100/pack		
30097	Tryptose Agar	10/pack	IVD	
402320	Tryptose Agar 100 mL	6/pack		4
24112	Tryptose Broth	20/pack	IVD	
455209	Tryptose Phosphate Broth (screw cap) 100 mL	25/pack	IVD	
402700	Tryptose SULFITE CYCLOSERINE A.B. 90 mL	6/pack		
402720	Tryptose Sulfite Cycloserine Agar Base 100 mL	6/pack		
412720	Tryptose Sulfite Cycloserine Agar Base 200 mL	6/pack		
24456	Tryptose SULPHITE Agar	20/pack	IVD	
463180	Tween 80 2% 600 mL (capacity 1000 mL)	6/pack		
80110	Urea 40% SUPPLEMENT 100 mL	6/pack	IVD	
24416	Urea Broth	20/pack	IVD	
403060	Urea INDOLE Broth 100 mL	6/pack	IVD	
20340	VAGITUBE	20/pack	IVD	
402460	Violet Red Bile Lactose Agar 100 mL	6/pack		
402540	Violet Red Bile Glucose Agar 100 mL	6/pack		
470031	Violet Red Bile Glucose Agar 500 mL	6/pack		5
450254	Violet Red Bile Glucose Agar 100 mL	25/pack		
452600	Violet Red Bile Glucose Agar 200 mL	25/pack		
34076	Violet Red Bile Lactose Agar	10/pack		

Ready-to-use culture media in tubes and bottles				
Ref.	Description	Packaging	IVD	Min. order qty (packs)
453140	Violet Red Bile Lactose Agar 200 mL	25/pack		
432300	Wort Agar 200 mL	6/pack		10
442300	Wurtz Lactose Agar 100 mL	6/pack	IVD	
402570	XLD Agar 100 mL	6/pack	IVD	
31102	Yeast Extract Agar	20/pack		10
403120	Yeast Extract Agar 100 mL	6/pack		
34074	Yeast Extract Agar	20/pack		
412120	Yeast Extract Agar 200 mL	6/pack		
26074	Yeast Extract Agar	100/pack		4
403090	Yeast Glucose Chloramphenicol Agar 100 mL	6/pack		
413090	Yeast Glucose Chloramphenicol Agar 200 mL	6/pack		
453091	Yeast Glucose Chloramphenicol Agar 200 mL	25/pack		
24432	Yersinia Broth	20/pack	IVD	
24431	Yersinia ITC Broth	20/pack		25
402610	YSG Broth 90 mL	6/pack		

GASPAK: Anaerobiosis, Microaerophilia, Capnophilia
Containers, Generators, Indicators, Disposable pouches.



Description	Packaging	Ref.
GasPak EZ Standard Container	1/pack	96720
GasPak EZ Standard Rack	1/pack	96721
GasPak EZ Large Container	1/pack	96722
GasPak EZ Large Rack	1/pack	96723
GasPak EZ Anaerobe Container System Sachet	20/pack	96724
Indicator of anaerobiosis	100/pack	96725
GasPak EZ Anaerobe Container System Sachet with indicator	20/pack	96726
GasPak EZ Campy Container System Sachet	20/pack	96727
GasPak EZ CO₂ Container System Sachet	20/pack	96728
GasPak CO₂ Indicator	50/pack	96729
GasPak EZ Anaerobe Pouch System	20/pack	96730
GasPak EZ CO₂ Pouch System	20/pack	96731
GasPak EZ Campy Pouch System	20/pack	96732

Dehydrated Culture Media			
Description		Packaging	Ref.
Acetamide Agar Medium for differentiation of nonfermentative, Gram-negative bacteria, especially <i>Pseudomonas aeruginosa</i> , on the basis of acetamide utilization.		500 g	610312
Acetamide Broth Broth for the differentiation of nonfermentative, Gram-negative bacteria, especially <i>Pseudomonas aeruginosa</i> , on the basis of acetamide utilization.		500 g	610313
Acetamide Broth ISO 16266 Liquid medium for detection of <i>P. aeruginosa</i> in water samples, according to ISO 16266.		500 g	610392
Aeromonas Agar Base Basal medium for selective isolation of <i>Aeromonas spp.</i> Ampicillin supplement		500 g 100 g 10 vials	610048 620048 81001
Alkaline Peptone Water Liquid enrichment medium for <i>Vibrio spp.</i> isolation.	APHA	500 g 100 g	610098 620098
Amies Transport Medium (with charcoal) Semi-solid medium for transport of clinical, environmental specimens and of microorganisms.		500 g 100 g 5 kg	610152 620152 6101525
Amies Transport Medium (w/o charcoal) Semi-solid medium for transport of clinical, environmental specimens and of microorganisms.		500 g 100 g	610191 620191
Amos Kont Jones Medium Liquid medium for detection of thermophilic microorganisms in food.		500 g	610386
Anaerobic Agar (Brewer) Medium for cultivating anaerobic microorganisms.		500 g	610320
Andrade Lactose Peptone Water Liquid medium for coliforms detection with andrade's indicator.		500 g	610118
Andrade Peptone Water Liquid enrichment medium with andrade's indicator.		500 g 100 g	610119 620119
Antibiotic Agar No.1 Surface medium for the antibiotic assay by Agar-diffusion method.	E.P.	500 g	610314
Antibiotic Broth No.3 Broth for turbidimetric assay of antibiotics.	U.S.P.	500 g	610316
Antibiotic Agar No.11 Medium for the microbiological assay of neomycin, used for the base and seed layers.	E.P.	500 g	610315
APT Agar Medium for counting and cultivating heterofermentative lactic acid bacteria. Tween 80 supplement	APHA	500 g 2x50 mL	610300 80031
Arginine Decarboxylase Broth Broth for the detection of the production of arginine decarboxylase by Enterobacteriaceae.		500 g	610306
Asparagine Enrichment Broth Liquid medium for <i>Pseudomonas aeruginosa</i> detection with MPN method in swimming pools and waters designed for bathing.	APHA	500 g 100 g	610138 620138
Azide Blood Agar Base Selective basal medium for streptococci isolation. Sheep Blood Defibrinated		500 g 100 g 50 mL	610153 620153 83296
Azide Dextrose Broth Selective liquid medium for fecal streptococci detection in water and food.	APHA	500 g 100 g	610003 620003
Azide Maltose Agar KF Selective medium for fecal streptococci isolation and enumeration. TTC 1% supplement	APHA	500 g 100 g 5x10 mL	610154 620154 80300
Bacillus Cereus Agar Base (MOSEL, MYP) Basal medium for selective isolation and enumeration of <i>Bacillus cereus</i> . Egg Yolk emulsion Supplement for lecithinase producers bacteria (<i>Staphylococcus</i> , <i>Bacillus</i> and <i>Clostridium</i>) isolation. <i>Bacillus cereus</i> supplement (Polymyxin B)	ISO 7932	500 g 100 g 4x50 mL 10 vials	610114 620114 80219 81016
Bacillus Cereus Agar Base (PEMBA) Basal medium for selective isolation and enumeration of <i>Bacillus cereus</i> . Egg Yolk emulsion Supplement for lecithinase producers bacteria (<i>Staphylococcus</i> , <i>Bacillus</i> and <i>Clostridium</i>) isolation. <i>Bacillus cereus</i> supplement (Polymyxin B)	ISO 21871	500 g 100 g 4x50 mL 10 vials	610136 620136 80219 81016
Baird Parker Agar Base Selective basal medium for <i>Staphylococcus aureus</i> isolation. Egg Yolk Tellurite emulsion Supplement for coagulase positive staphylococci isolation in food. R.P.F. supplement Selective and diagnostic supplement for isolation, enumeration and confirmation of <i>Staphylococcus aureus</i> from food and other specimens.	ISO 6888-1	500 g 100 g 4x50 mL 8 vials	610004 620004 80122 81057

Dehydrated Culture Media			
Description		Packaging	Ref.
Baird Parker Agar Base acc. to E.P. Selective basal medium for <i>Staphylococcus aureus</i> isolation.	E.P.	500 g	610338
Egg Yolk Tellurite emulsion Supplement for coagulase positive staphylococci isolation in food.	ISO 6888-1	4x50 mL	80122
Barnes Agar Base Selective basal medium for fecal streptococci isolation and identification.		500 g	610156
TTC 1% supplement		100 g	620156
		5x10 mL	80300
BAT Agar Medium for detection of <i>Alicyclobacillus spp.</i> in fruit juices and other beverages.		500 g	610361
Beer Universal Agar Medium for isolation of microorganisms contaminating beer.		500 g	610126
		100 g	620126
Biggy (Nickerson) Agar Selective medium for <i>Candida spp.</i> isolation and presumptive identification.		500 g	610135
		100 g	620135
Bile Aesculin Agar Medium for group D streptococci isolation.		500 g	610210
		100 g	620210
Bile Aesculin Broth Selective medium for group D streptococci isolation.	ISO 7899-2	500 g	610151
		100 g	620151
Bile Aesculin Azide Agar Selective medium for group D streptococci isolation.	ISO 7899-2	500 g	610001
		100 g	620001
Bismuth Sulfite Agar Selective medium for isolation and differentiation of <i>Salmonella typhi</i> and other salmonellae from clinical specimens and other materials.	U.S.P.	500 g	610301
Bismuth Sulphite Agar Selective medium for detection of <i>Salmonella spp</i> in food and and other materials.	ISO 6579-1	500 g	610387
Blood Agar Base Basal medium for fastidious microorganisms isolation and for haemolysis test.		500 g	610005
Sheep Blood Defibrinated		100 g	620005
		50 mL	83296
Blood Agar Base No.2 Basal medium for fastidious microorganisms isolation and for haemolysis test.		500 g	610188
Sheep Blood Defibrinated		100 g	620188
		50 mL	83296
Bolton Broth Basal Broth for enrichment of <i>Campylobacter spp.</i>	ISO 10272-1	500 g	610216
Campylobacter C.T.V.A. supplement (Cefoperazone, Trimethoprim, Vancomycin, Amphotericin B)	ISO 10272-1	100 g	620216
Selective supplement for enrichment of <i>Campylobacter spp.</i>		10 vials	81077
Bordet Gengou Agar Base Basal medium for selective isolation of <i>Bordetella spp.</i>		500 g	610006
Bordetella supplement (Cephalexin)		10 vials	81013
Sheep Blood Defibrinated		50 mL	83296
Brain Heart Infusion Agar Medium for cultivation of fastidious aerobic and anaerobic bacteria.		500 g	610007
		100 g	620007
Brain Heart Infusion Broth Medium for cultivation of fastidious aerobic and anaerobic bacteria.		500 g	610008
		100 g	620008
		5 kg	6100085
Brilliant Green Agar Selective medium for <i>Salmonella spp.</i> isolation other than <i>S. typhi</i> .	U.S.P.	500 g	610009
		100 g	620009
Brilliant Green Agar w/ Sulfadiazine Selective medium for the isolation of <i>Salmonella spp.</i>		500 g	610373
Brilliant Green Bile Broth 2% Selective medium for detection and confirmation of coliforms.	ISO 4831	500 g	610010
	ISO 4832	100 g	620010
		5 kg	6100105
Bromocresol Purple Azide Broth Broth for the confirmation of the enterococci presence, particularly in the bacteriological analysis of water.		500 g	610321
Glycerol supplement (Glycerol)		4x50 mL	80021
Brucella Agar Base Basal medium for selective isolation of <i>Brucella spp.</i> and other bacteria.	APHA	500 g	610079
Brucella supplement (Polymyxin B, Bacitracin, Cycloheximide, Nalidixic Acid, Nystatin, Vancomycin)		100 g	620079
		10 vials	81003
Brucella Broth Medium for the cultivation of <i>Brucella spp.</i> and other fastidious microorganisms.		500 g	610356
Bryant and Burkey Medium Medium for the selective enrichment of lactate-fermenting clostridia (<i>C. tyrobutyricum</i>) which are responsible for "late blowing" in brine-salted semi-hard and hard cheeses and other dairy products.		500 g	610142
Buffered Peptone Water Medium for <i>Salmonella</i> pre-enrichment in food.	ISO 6579	500 g	611014
	ISO 6785	100 g	621014
		5 Kg	6110145

Dehydrated Culture Media			
Description		Packaging	Ref.
Buffered Sodium Chloride Peptone Solution acc. to E.P. Medium for the microbial examination of non-sterile pharmaceutical products, e.g. solution for sample preparation or dilution or washing solution.	ISO 21149 E.P.	500 g	610335
Campylobacter Agar Basal medium for selective isolation of <i>Campylobacter spp.</i>	AOAC	500 g 100 g	611007 621007
Sheep Blood Defibrinated		50 mL	83296
Lysed horse blood		100 mL	83397
Campylobacter growth supplement (Sodium pyruvate, Sodium metabisulphate, Ferrous sulphate)		10 vials	81050
Campylobacter Blaser Wang supplement (Polymyxin B, Trimethoprim, Amphotericin B, Cephalotin, Vancomycin)		10 vials	81051
Campylobacter Preston supplement (Polymyxin B, Trimethoprim, Cyclohexinide, Rifampicin)		10 vials	81004
Campylobacter Skirrow supplement (Polymyxin B, Trimethoprim, Vancomycin)		10 vials	81055
Campylobacter Blood Free Medium Base Basal medium for selective isolation of <i>Campylobacter jejuni</i> , coli, landis from clinical samples.	ISO 10272-1	500 g 100 g	610130 620130
Campylobacter C.C.D.A. supplement (Cefoperazone, Amphotericin B).	ISO 10272-1	10 vials	81037
Campylobacter Karmali Agar Base Basal medium for selective isolation of <i>Campylobacter spp.</i>	AOAC	500 g 100 g	610200 620200
Campylobacter Karmali supplement (Sodium pyruvate, Cefoperazone, Vancomycin, Cyclohexinide)		10 vials	81036
Cary Blair Transport Medium Semi-solid medium for transport of clinical specimens and microorganisms.		500 g 100 g	611402 621402
Casein Peptone Lecithin Polysorbate Broth Base Broth for diluting samples of pharmaceutical, cosmetic and other raw materials or final products when determining microbial counts.	U.S.P.	500 g	610302
Tween 20 supplement		2x50 mL	80032
Cetrimide Agar Selective basal medium for <i>Pseudomonas aeruginosa</i> isolation.	E.P.	500 g 100 g 5 kg	610041 620041 6100415
Glycerol supplement		4x50 mL	80021
China-Blue Lactose Agar Non-selective medium for enumeration and differentiation of bacteria in dairy products.		500 g	610362
Chromatic™ EC X-GLUC Agar Chromogenic medium for detection and enumeration of E. coli in water and foodstuffs.		500 g 100 g	610602 620602
CLED Agar Differential medium for total bacterial count of urine.		500 g 100 g 5 kg	610012 620012 6100125
CLED Andrade Agar Medium for total bacterial count of urine.		500 g 100 g	610112 620112
Clostridium Broth Semisolid medium for clostridia and other anaerobes cultivation in food and clinical samples.	ISO 6461-1 E.P.	500 g 100 g	610056 620056
Clostridium difficile Agar Base Basal medium for selective isolation of <i>Clostridium difficile</i> .		500 g 100 g	610115 620115
Sheep Blood Defibrinated		50 mL	83296
Clostridium difficile supplement (D-Cyclozenine, Cefoxitin)		10 vials	81007
Clostridium Perfringens Agar Base Basal medium for selective isolation of <i>Clostridium perfringens</i> .	ISO 7937 ISO 14189	500 g 100 g	610207 620207
Egg Yolk emulsion	ISO 7932	4x50 mL	80219
Enrichment supplement for lecithinase producing bacteria (<i>Staphylococcus</i> , <i>Bacillus</i> and <i>Clostridium</i>) isolation.			
Clostridium perfringens (T.S.C.) supplement (D-Cycloserine).	ISO 7937	10 vials	81011
Kanamycin / Polymyxin B supplement		10 vials	81031
Columbia Agar Base Basal medium for fastidious microorganisms isolation.		500 g 100 g 5 kg	610013 620013 6100135
Sheep Blood Defibrinated		50 mL	83296
CNA (Staf/Strep) supplement (Colistin, Nalidixic acid)		10 vials	81048
<i>Gardnerella vaginalis</i> supplement (Gentamicin, Nalidixic acid, Amphotericin B)		10 vials	81040
Columbia CNA Agar Base Selective medium for isolation of streptococci and staphylococci.		500 g	610113
Sheep Blood Defibrinated		50 mL	83296
Cooked Meat Medium Medium for cultivating anaerobic bacteria.		500 g	610372
Corn Meal Agar Medium for the cultivation of fungi and chlamydospore production.		500 g 100 g	610123 620123

Dehydrated Culture Media			
Description		Packaging	Ref.
C.P.L.M. Selective with CAF Selective medium for <i>Trichomonas vaginalis</i> isolation and cultivation.		500 g	610159
Horse Serum		100 mL	83810
Cronobacter Selective Broth Liquid medium for the selective enrichment of <i>Cronobacter spp</i> (formerly <i>Enterobacter sakazakii</i>).	ISO 22964	500 g	610389
Vancomycin (5 mg) Supplement		10 vials	81064
Cystine Heart Agar Medium for cultivating <i>Francisella tularensis</i> and gram-negative cocci.		500 g	610506
Czapek Dox Agar Medium for cultivation of fungi and bacteria capable of using sodium nitrate as the sole nitrogen source.	ISO 6461-1	500 g 100 g	610095 620095
Czapek Dox Broth Liquid medium or cultivation of fungi and bacteria capable of using sodium nitrate as the sole nitrogen source.		500 g 100 g	610072 620072
DNase Test Agar Medium for detection of deoxyribonucleasic activity (DNase) of microorganisms.		500 g 100 g	610205 620205
Demi Fraser Broth Basal Broth for <i>Listeria</i> enrichment.	ISO 11290-1	500 g 100 g 5 kg	610091 620091 6100915
<i>Listeria</i> Fraser supplement (Ferric ammonium citrate)		10 vials	81028
<i>Listeria</i> Fraser supplement (liquid)		10x5 mL	80304
Dermatophyte (D.T.M.) Agar Selective medium for selective isolation of dermatophytic fungi.		500 g 100 g	610160 620160
Dermatophyte supplement (Chlortetracycline)		10 vials	81025
Desoxycholate Agar Selective medium for enterobacteria isolation.	APHA	500 g 100 g 5 kg	610014 620014 6100145
Desoxycholate Citrate Agar Selective medium for pathogenic enterobacteria isolation.		500 g 100 g	610015 620015
Dextrose Agar Medium for cultivation and maintaining of microorganisms.		500 g	610002
Dextrose Broth Broth for cultivation and maintaining of microorganisms.		500 g 100 g	610161 620161
Dextrose Tryptone Agar Medium for the enumeration of spores of mesophilic and thermophilic <i>Bacillus</i> in sugar, sweet desserts, aromatic preparation and other foodstuffs.		500 g	610198
Dextrose Tryptone Broth Broth for the enumeration of spores of mesophilic and thermophilic bacillus in sugar, sweet desserts, aromatic preparation and other foodstuffs.		500 g	610376
Dichloran-Glycerol (DG18) Agar Base Medium for counting yeasts and moulds in foodstuffs.		500 g	610238
Dichloran Rose-Bengal Chloramphenicol Agar Selective medium for counting yeasts and moulds.		500 g	610237
Drigalski Lactose Agar Selective medium for pathogenic enterobacteria isolation.		500 g 100 g	610016 620016
D/E Neutralizing Broth Broth for the neutralization and testing of antiseptics and disinfectants according to the procedure of Engley and Dey.		500 g	610088
Tween 80 supplement		2x50 mL	80031
E.C. Broth Selective Broth for fecal coliforms detection and for <i>E. coli</i> confirmation.	ISO 7251	500 g 100 g	610063 620063
E.C. Broth MUG Selective broth for fecal coliforms and <i>E. coli</i> detection in water.		500 g	610332
E.E. Broth Broth for enterobacteria selective enrichment.	E.P.	500 g 100 g	610017 620017
EMB Levine Agar Selective medium for gram-negative enteric bacteria isolation.	U.S.P.	500 g 100 g	610019 620019
Endo Agar Medium for coliforms detection, enumeration and confirmation.	APHA	500 g 100 g	610020 620020
Enterobacter (Cronobacter) sakazakii Isolation Agar Selective medium for <i>Cronobacter sakazakii</i> isolation.	ISO 22964	500 g	610623
Eugon Broth General purpose medium for cultivating a wide variety of microorganisms from clinical and non-clinical specimens.		500 g	610341
E.V.A. Broth (Ethyl Violet Azide) Selective broth for enterococci detection in water and other samples.	APHA	500 g 100 g	611365 621365
Fermentation Broth Base Liquid medium for the carbohydrate utilization test.	ISO 11290	500 g	610030

Dehydrated Culture Media			
Description		Packaging	Ref.
Fluid Thioglycollate Medium Medium for aerobic and anaerobic bacteria for sterility tests.	ISO 7937 E.P. U.S.P.	500 g 100 g 5 kg	610050 620050 6100505
Fraser Broth Base Enrichment medium for detection and enumeration of <i>Listeria monocytogenes</i> and <i>Listeria</i> spp. Half Fraser Supplement Fraser Supplement	ISO 11290	500 g 10 vials 10 vials	610375 81043 81046
Gassner Agar Selective agar for the detection and isolation of pathogenic Enterobacteriaceae in food and other material according to Gassner.		500 g	610366 #
G.C. Medium Basal medium for <i>Haemophilus</i> spp. and <i>Neisseria</i> spp. isolation. Horse Blood Defibrinated VITALEX growth supplement L.C.A.T. supplement (Lincomycin, Colistin, Amphotericin B, Trimethoprim) V.C.A.T. supplement (Vancomycin, Colistin, Amphotericin B, Trimethoprim) V.C.N. supplement (Vancomycin, Colistin, Nystatin) V.C.N.T. supplement (Vancomycin, Colistin, Nystatin, Trimethoprim) HAEMOPHILUS supplement (Bacitracin, Vancomycin, Clindamycin)		500 g 100 g 100 mL 10 vials 10 vials 10 vials 10 vials 10 vials 10 vials	610022 620022 83396 81023 81012 81041 81022 81024 81014
Gelatin Peptone Agar (Agar gelisato) Medium for total bacterial count of aerobes in ice creams.	AOAC	500 g 5 kg	610054 6100545
Giolitti Cantoni Broth Base Basal broth for <i>Staphylococcus aureus</i> enrichment. Potassium Tellurite 3.5% supplement	ISO 6888-3	500 g 100 g 5x10 mL	610100 620100 80291
Glucose OF Medium ISO 21528 Semisolid medium for confirmation of Enterobacteriaceae.	ISO 21528	500 g	610388
Glucose Salt Teepol Broth Medium for selective enrichment and enumeration (MPN method) of <i>Vibrio parahemolyticus</i> from foodstuffs.		500 g 100 g	610359 610359
G.N. Hajna Broth Broth for gram-negatives enrichment.		500 g	610163
Hektoen Enteric Agar Selective medium for pathogenic enterobacteria isolation.	ISO 21567	500 g 100 g 5 kg	610021 620021 6100215
Herellea Agar Selective medium for isolation from urethral and vaginal specimens.		500 g 5 kg	610164 6101645
HHD Broth Medium for differentiation of heterofermentative and homofermentative lactic bacteria.		500 g 100 g	610064 # 620064 #
Iron Sulfite Agar Medium for thermophilic anaerobes determination.		500 g 100 g	611401 621401
Isosensitest Agar Medium for susceptibility test.		500 g 100 g	611265 621265
Kanamycin Aesculin Azide Agar Selective medium for isolation of group D streptococci in food.	APHA	500 g 100 g	610190 620190
Kligler Iron Agar Differential medium for enterobacteria identification.	ISO 10273	500 g 100 g	610023 620023
Lactose Broth Broth for coliforms detection.	ISO 9308-2	500 g 100 g	611202 621202
Lactose Gelatin Medium Dehydrated medium for the identification of <i>C. perfringens</i> by lactose and gelatin test.	ISO 7937	500 g	610360
Lactose Sulfite Broth Selective medium for the detection and enumeration of spores of <i>Clostridium perfringens</i> in food products.	ISO 7937	500 g	610358
Lauryl Pepto Broth Broth for selective enrichment of coliforms in water and dairy products.	ISO 7251	500 g 100 g	610189 620189
Lauryl Sulfate Tryptose Broth Modified (m-LST) Selective medium for the enrichment of <i>Cronobacter sakazakii</i> .	ISO 22964	500 g	610247
Lauryl Tryptose Broth Broth for selective enrichment of coliforms in water and dairy products.	ISO 7251	500 g 100 g	610085 620085
Koser Citrate Medium Medium for enterobacteria differentiation on the basis of their capability of utilizing citrate as the sole carbon source.	AOAC	500 g 100 g	610165 620165
LB Agar Medium for cultivation of recombinant <i>Escherichia coli</i> .		500 g	610245
LB Medium (Luria Bertani) Medium for cultivation of recombinant <i>Escherichia coli</i> .		500 g 100 g	610084 620084

Dehydrated Culture Media			
Description		Packaging	Ref.
Legionella CYE Agar Base (w/o ACES and α-Chetoglutaric Acid) Basal medium for <i>Legionella spp.</i> isolation.	ISO 11731	500 g 100 g	610125 620125
Legionella (BCYE) growth supplement (Aces, L-cysteine, Ferric pyrophosphate, α-ketoglutarate).	ISO 11731	10 vials	81056
Legionella (BCYE) growth supplement w/o L-Cysteine		10 vials	81091
Legionella (BMPA) supplement (Cefamandole, Polymyxin B, Anisomycin)		10 vials	81002
Legionella (GVPC) supplement (Glycine, Vancomycin, Polymyxin B, Cycloheximide).	ISO 11731	10 vials	81008
Legionella (MWY) supplement (Anisomycin, Polymyxin B, Vancomycin, Glycine, Bromothynol Blue, Bromocrystal Purple).		10 vials	81019
Lethen Agar Base Basal medium for coefficient phenol detection of quaternary ammonium disinfectants.		500 g 100 g	610076 620076
Lethen supplement (Lecithin)		4x50 mL	80008
Tween 80 supplement		2x50 mL	80031
Lethen Broth Base Basal medium for coefficient phenol detection of quaternary ammonium disinfectants.	AOAC	500 g 100 g	610208 620208
Lethen supplement (Lecithin)		4x50 mL	80008
Tween 80 supplement		2x50 mL	80031
Lethen Broth Modified Medium for the evaluation of bactericidal activity of quaternary ammonium base.	ISO 21149	500 g	610328
Tween 80 supplement		2x50 mL	80031
Listeria Enrichment Broth Base (Lovett) Broth for enrichment of <i>Listeria spp.</i> according to Lovett formulation.	APHA	500 g	610326
Listeria Supplement for Enrichment (FDA/IDF)		10 vials	81045
Selective supplement for <i>Listeria</i> isolation.			
Listeria Fraser Broth Basal medium for <i>Listeria monocytogenes</i> selective enrichment.	ISO 11290-1	500 g 100 g	610166 620166
Listeria Fraser supplement (Ferric ammonium citrate)		10 vials	81028
Listeria Fraser supplement (liquid)		10x5 mL	80304
Listeria Oxford Agar Base Basal medium for selective isolation and differentiation of <i>L. monocytogenes</i> .	ISO 11290-1	500 g 100 g	610167 620167
Listeria Oxford supplement (Cefotetan, Fosfomycin)		10 vials	81027
Listeria Oxford Agar Base ISO 11290 Basal medium w/o antibiotics and acriflavine for the isolation of <i>L. monocytogenes</i> .	ISO 11290-1	500 g	610231
Listeria Oxford with CEF/FOS/COL/CYC/ACR supplement (Cefotetan, Fosfomycin, Colistin, Cycloheximide, Acriflavine)		10 vials	81629
Listeria Palcam Agar Base Basal medium for selective isolation and differentiation of <i>L. monocytogenes</i> .	ISO 11290-1	500 g 100 g	610168 620168
Listeria Palcam supplement (Polymyxin B, Ceftazidime, Acriflavina HCl)		10 vials	81026
Listeria U.V.M. 1 Broth Broth for <i>Listeria monocytogenes</i> selective enrichment.	USDA-FSIS	500 g 100 g	610031 620031
Listeria U.V.M. 2 Broth Broth for <i>Listeria monocytogenes</i> selective enrichment.	USDA-FSIS	500 g	610045
Litmus Milk Medium for the maintenance of lactic acid bacteria and for the determination of bacterial action on milk.		500 g	610318
Lowenstein Jensen Medium Selective medium for mycobacteria isolation.		500 g 100 g	610026 620026
Glycerol supplement (Glycerol)		4x50 mL	80021
Lysine Decarboxylase Broth Broth for the detection of the production of Lysine decarboxylase by <i>Enterobacteriaceae</i> .	APHA	500 g	610303
Lysine Iron Agar Differential medium for enterobacteria isolation.		500 g 100 g	610027 620027
M17 Agar Medium for lactic streptococci cultivation.		500 g 100 g	610192 620192
M.R.S. Agar Basal medium for <i>Lactobacillus</i> isolation.	ISO 15214	500 g 100 g	610024 620024
Tween 80 supplement		2x50 mL	80031
M.R.S. Broth Basal medium for <i>Lactobacillus</i> enrichment and isolation.	ISO 15214	500 g 100 g	610025 620025
Tween 80 supplement		2x50 mL	80031
M.S.R.V. (Modified Semisolid Rappaport Vassiliadis) Medium Base Semisolid basal medium for selective isolation of motile <i>Salmonella spp.</i>		500 g 100 g	610018 620018
Novobiocin supplement	ISO 6579	10 vials	81021
MacConkey Agar Selective and differential medium for coliforms and intestinal pathogens isolation.	E.P.	500 g 100 g 5 kg	610028 620028 6100285

Dehydrated Culture Media			
Description		Packaging	Ref.
MUG supplement		10 vials	81020
MacConkey Agar w/o Crystal Violet Differential moderately selective medium for enterobacteria isolation.		500 g 100 g	610195 620195
MacConkey Agar w/o Bile Salt Selective and differential medium for coliforms isolation with <i>Proteus spp.</i> swarming inhibition.		500 g	610128
MacConkey Broth Selective and differential broth for the detection of coliforms in water, milk and other materials of sanitary importance.	E.P.	500 g	610337
MacConkey Broth (Purple) Selective and differential broth for coliforms isolation (with BCP).		500 g 100 g	610171 620171
MacConkey Broth Modified (w/o NaCl) Selective and differential broth for coliforms isolation (with BCP).		500 g 100 g	610149 620149
MacConkey MUG Agar Selective and differential medium for <i>Escherichia coli</i> isolation.		500 g 5 kg	610170 6101705
MacConkey Sorbitol Agar Selective and differential medium for <i>E. coli</i> O157 isolation.		500 g 100 g	610108 620108
Cefixime Tellurite Supplement		10 vials	81082
Malonate Broth Broth for enterobacteria differentiation.		500 g 100 g	610172 620172
Malt Agar Medium for yeasts and moulds isolation and enumeration.	AOAC	500 g 100 g	610101 620101
Malt Extract Agar Medium for yeasts and moulds isolation and enumeration.		500 g 100 g	610173 620173
Malt Extract Broth Broth for yeasts and moulds culture.	AOAC	500 g 100 g	610069 620069
Mannitol Motility Test Medium Medium for differentiation of Enterobacteriaceae.		500 g 100 g	610235 620235
Mannitol Salt Agar Selective medium for pathogenic staphylococci isolation.	EP, USP, JP	500 g 100 g 5 kg	610029 620029 6100295
Maximum Recovery Diluent Isotonic diluent for microorganisms optimal recovery.	ISO 6887-1	500 g 100 g	610077 620077
m-CP Agar Base Basal medium for selective isolation of <i>C. perfringens</i> in waters.	98/83/EC	500 g 100 g	610199 620199
m-CP supplement		10 vials	81034
m-CP Supp.A Phenolphthalein diphosphate		10 vials	81070
m-CP Supp.B Ferric chloride		10 vials	81071
m-CP Supp.C Indoxyl B-D-Glucoside		10 vials	81072
m-ENDO Agar LES Medium for coliforms enumeration in water by membrane filtration.	APHA	500 g 100 g	611011 621011
m-ENDO Broth Broth for coliforms counting in water by membrane filtration.	APHA	500 g	610220
m-FC Agar Base Selective basal medium for fecal coliforms isolation and enumeration by membrane filtration method.	APHA	500 g 100 g 5 kg	611012 621012 6110125
Rosolic acid supplement		10 vials	81029
m-FC Broth Base Selective basal Broth for fecal coliforms isolation and enumeration by membrane filtration method.	APHA	500 g 100 g	610162 620162
Rosolic acid supplement		10 vials	81029
Microbial Content Test Agar Base Basal medium for bacterial count determination with inactivation of disinfectants.		500 g	610094
Lecithin supplement		4x50 mL	80007
Tween 80 supplement		2x50 mL	80031
Middlebrook 7 H 10 Agar Base Basal medium for mycobacteria cultivation.		500 g	611022
Middlebrook 7 H 10 supplement (Oleic acid, Glucose, NaCl)		4x50 mL	81035
Glycerol supplement (Glycerol)		4x50 mL	80021
Milk Plate Count Agar Medium for bacterial count in milk and dairy products.		500 g 100 g	610073 620073
Motility Indole Urea Agar Medium for Enterobacteriaceae identification.		500 g	610236
Motility Nitrate CP Medium Selective medium for motile nitrate-utilizing microorganisms (<i>P. aeruginosa</i> , <i>C. perfringens</i>).		500 g	610371

Dehydrated Culture Media			
Description		Packaging	Ref.
Motility Test Agar Medium for the detection of motility of gram-negative enteric bacilli.		500 g 100 g	610132 620132
MRS Broth w/o Glucose Liquid medium for cultivating lactobacilli from clinical and nonclinical specimens. Tween 80 supplement		500 g 100 g 2x50 mL	610144 620144 80031
MRS-IM + Maltose Agar Medium for counting probiotic bacteria. Tween 80 supplement		500 g 2x50 mL	610067 80031
MRS-IM + Glucose Agar Medium for counting probiotic bacteria. Tween 80 supplement Dicloxacillin Supplement		500 g 2x50 mL 10 vials	610068 80031 81094
MR-VP Broth Broth for carrying out Methyl Red (MR) and Voges-Proskauer (VP) tests.	ISO 6785 ISO 6579	500 g 100 g	610032 620032
Mueller Hinton Agar II Medium for susceptibility test (Kirby Bauer method), low concentration of thymine-thymidine with Ca ⁺⁺ and Mg ⁺⁺ .	CLSI	500 g 100 g 5 kg	610627 620627 6106275
Mueller Hinton Broth II Broth for susceptibility test by dilution method, low concentration of thymine-thymidine with Ca ⁺⁺ and Mg ⁺⁺ .	CLSI	500 g 100 g	610218 620218
Muller Kauffmann Broth Base Basal broth for <i>Salmonella spp.</i> selective enrichment. Muller Kauffmann supplement (Iodine, Potassium iodide, Brilliant green)	ISO 6785	500 g 100 g 2x50 mL Iodine 50 mL B.G. 0.1%	610035 620035 80047
Muller Kauffmann Tetrathionate Broth Base Basal medium for <i>Salmonella spp.</i> selective enrichment from water, foodstuffs and fecal samples. Iodine MKTT Solution Novobiocin MKTT Supplement	ISO 6579	500 g 100 g 10 tubes x 10 mL 10 vials	610239 620239 80009 81073
Nitrate Broth Broth for the identification of aerobic and facultative anaerobic gram-negative microorganisms by means of the nitrate reduction test.		500 g	610322
Nutrient Agar Medium for non fastidious bacteria growth.	ISO 16266	500 g 100 g 5 kg	610036 620036 6100365
Nutrient Agar Modified Medium for the cultivation of non-fastidious bacteria from clinical and non clinical specimens.		500 g	610385
Nutrient Broth Broth for non fastidious bacteria growth.	APHA	500 g 100 g 5 kg	610037 620037 6100375
Nutrient Broth no.2 ISO 10272-1 For the cultivation and enrichment of less fastidious bacteria for the establishment of permanent cultures, as well as a base in the preparation of special media, according to ISO 10272-1	ISO 10272-1	500 g	610393
O.F. Basal Medium Medium for the detection of oxidative and fermentative carbohydrate degradation, particularly for the differentiation of gram-negative intestinal bacteria.		500 g	610304
O.G.Y.E. Agar Base Basal medium for selective enumeration of yeasts and moulds. O.G.Y.E. supplement (Oxytetracycline)	ISO 21527-1 ISO 21527-2	500 g 100 g 5 kg 10 vials	610202 620202 6102025 81018
Orange Serum Agar Medium for isolation and count of bacteria and yeasts in citrus fruit and by-products.	APHA	500 g 100 g 5 kg	610078 620078 6100785
Orange Serum Broth Broth for cultivation and enumeration of microorganisms associated with spoilage of citrus products.		500 g	610307
Ornithine Decarboxylase Broth Broth for the detection of the production of ornithine decarboxylase by <i>Enterobacteriaceae</i> .		500 g	610305
P.P.L.O. Broth Broth for mycoplasma cultivation.	APHA	500 g	610109
Peptone Water Medium for indole production detection.	ISO 7251	500 g 100 g	610038 620038
Pfizer Selective Enterococcus Agar Medium for selective isolation, cultivation and enumeration of <i>Enterococcus spp.</i>		500 g	610319

Dehydrated Culture Media			
Description		Packaging	Ref.
Phenol Red Agar Base Base medium used with the addition of carbohydrates for differentiating pure culture of bacteria based on sugar fermentation reactions.		500 g	610308
Phenol Red Broth Base Broth for carbohydrates fermentation studies.		500 g 100 g	610174 620174
Phenylalanine Agar Medium for enterobacteria differentiation. Phenylalanine Reagent		500 g 100 g 2x25 mL	610039 620039 80272
Plate Count Agar Medium for total bacterial count detection in foodstuffs and water.	ISO 4833	500 g 100 g 5 kg	610040 620040 6100405
Potato Dextrose Agar Medium for yeasts and moulds enumeration in food and dairy products.	E.P.	500 g 100 g 5 kg	610102 620102 6101025
Potato Dextrose Broth Broth for yeasts and moulds enumeration in food and dairy products.	APHA	500 g 100 g 5 kg	610106 620106 6101065
Potato Infusion Agar Medium for <i>Brucella spp.</i> isolation.		500 g 100 g	610011 620011
Proteose Peptone Special peptone obtained by enzymatic hydrolysis of animal tissue.		500 g	612201
Pseudomonas Agar Base Basal medium for selective isolation of <i>Pseudomonas spp.</i> Glycerol supplement (Glycerol) C.F.C. (Pseudomonas) supplement (Cetrimide, Fucidin, Cephalordine) C.N. (Pseudomonas) supplement (Cetrimide, Nalidix acid) Pseudomonas PP Supplement	ISO 12780	500 g 100 g 4x50 mL 10 vials 10 vials 10 vials	610071 620071 80021 81049 81006 81093
Pseudomonas Agar F Medium for the enhancement of fluoescin production by <i>Pseudomonas spp.</i> Glycerol supplement	ISO 12780	500 g 4x50 mL	610309 80021
Pseudomonas Agar P Medium for the enhancement of pyocyanin production by <i>Pseudomonas spp.</i>	U.S.P.	500 g	610310
Purple Glucose Agar Medium for enterobacteria isolation and differentiation.	ISO 21528-2	500 g 100 g	610194 620194
Purple Lactose Agar Medium for coliforms isolation and differentiation.		500 g 100 g	610044 620044
R2 Agar Medium for heterotrophic microorganisms count in the examination of potable water.	E.P.	500 g 100 g	610129 620129
Rappaport Vassiliadis Broth Selective enrichment Broth for Salmonella isolation.	ISO 6785 ISO 6579	500 g 100 g	610175 620175
Reinforced Clostridial Agar Medium for isolation and enumeration of anaerobes, specially <i>Clostridium spp.</i>	ISO 6461-1	500 g 100 g	610096 620096
Rogosa Agar Selective basal medium for lactobacilli isolation and enumeration. Tween 80 supplement		500 g 100 g 2x50 mL	610176 620176 80031
Rogosa Broth Selective basal broth for lactobacilli culture. Tween 80 supplement		500 g 100 g 2x50 mL	610177 620177 80031
Rose Bengal Agar Base Basal Broth for selective isolation and enumeration of yeasts and moulds. Chloramphenicol supplement	APHA	500 g 100 g 5 kg 10 vials	610178 620178 6101785 81017
Rose Bengal CAF Agar Selective medium for yeasts and moulds isolation and enumeration.	APHA	500 g 100 g	610090 620090
Sabouraud Dextrose Agar Medium for yeasts and moulds isolation. Chloramphenicol supplement Gentamicin supplement	harmonized EP	500 g 100 g 5 kg 10 vials 10 vials	610103 620103 6101035 81017 81033
Sabouraud Dextrose Broth Broth for yeasts, moulds and acidophilic bacteria cultivation.	harmonized EP	500 g 100 g	610104 620104
Sabouraud CAF Agar Selective medium for yeasts and moulds isolation.	APHA	500 g 100 g	610203 620203
Sabouraud CAF (50 mg) Agar Selective medium for the isolation, identification and cultivation of pathogenic and nonpathogenic fungi.		500 g	610625
Sabouraud CAF Agar + Actidione Selective medium for dermatophyte isolation.		500 g 100 g	610179 620179

Dehydrated Culture Media			
Description		Packaging	Ref.
Sabouraud Maltose Agar Acid medium for the isolation of dermatophytes, other fungi and yeasts.	U.S.P.	500 g 100 g	610146 620146
Schaedler Agar Base Basal medium for anaerobic bacteria isolation.		500 g 100 g	610043 620043
Sheep Blood Defibrinated		50 mL	83296
Vitamin K 1% supplement		5x5 mL	80053
Schaedler supplement (Kanamycin, Vancomycin)		10 vials	81054
Schaedler Broth Broth for anaerobic bacteria cultivation.		500 g 100 g	610137 620137
Schwarz Differential Medium Medium used in the brewing industry for the differentiation of brewing yeasts from wild yeasts.		500 g	610127
Selenite Broth Broth for selective enrichment of <i>Salmonella spp.</i>	APHA	500 g 100 g 5 kg	610145 620145 6101455
Selenite Cystine Broth Medium for selective enrichment of <i>Salmonella spp.</i>	ISO 6785	500 g 100 g	610150 620150
S.F. Broth Broth medium for fecal streptococci isolation.	AOAC	500 g 100 g	610180 620180
Sierra Lipolytic Agar Medium for lipolytic bacteria isolation.		500 g	611015
SIM Medium Differential medium for indole, motility and H ₂ S tests.		500 g 100 g	610181 620181
Simmons Citrate Agar Differential medium for enterobacteria identification based on the utilization of citrate as the sole source of carbon.	ISO 10273	500 g 100 g	610046 620046
Slanetz-Bartley Agar Base Selective basal medium for fecal streptococci isolation and enumeration.	ISO 7899-2	500 g 100 g	610134 620134
TTC 1% supplement		5x10 mL	80300
Slanetz-Bartley Agar + TTC Selective medium for fecal streptococci isolation and enumeration.	ISO 7899-2	500 g 100 g 5 kg	610147 620147 6101475
SPS Agar Selective medium for <i>Clostridium perfringens</i> and sulphite-reducing bacteria isolation and enumeration.	APHA	500 g 100 g	610148 620148
S.S. Agar Modified Selective medium for <i>Salmonella</i> and <i>Shigella</i> isolation.		500 g 100 g 5 kg	610042 620042 6100425
Staphylococcus 110 Agar Selective medium for staphylococci isolation.		500 g	611366
Stuart Transport Medium Semisolid medium for transport of clinical specimens and microorganisms.		500 g 100 g	610182 620182
TBX Agar Selective chromogenic medium for <i>E. coli</i> detection in food.	ISO 16649-2	500 g 100 g	610224 620224
T.C.B.S. Agar Selective medium for <i>Vibrio spp.</i> isolation.		500 g 100 g	611010 621010
Tergitol Agar Base Basal medium for coliforms enumeration and identification.	ISO 9308-1	500 g 100 g	610184 620184
TTC 1% supplement		5x10 mL	80300
T.S.I. Agar Differential medium for enterobacteria identification.	U.S.P.	500 g 100 g	610055 620055
Tetrathionate Broth Base Basal broth for <i>Salmonella spp.</i> enrichment	APHA U.S.P.	500 g 100 g	610183 620183
Iodine solution (Iodine, Potassium iodide)		5x10 mL	81009
Todd Hewitt Broth Broth for haemolytic streptococci culture prior to serological identification.		500 g 100 g 5 kg	610051 620051 6100515
TOS Propionate Agar Base Basal medium for detection of bifidobacteria in milk products.	ISO 29981	500 g	610378
MUP Selective Supplement		10 vials	81101
Trichomonas Broth Selective base Broth for <i>Trichomonas spp.</i> isolation.		500 g 100 g	610061 620061
Horse Serum		100 mL	83810
Tryptic (CTA) Medium Semisolid medium for motility test and carbohydrates fermentation.		500 g 100 g	610185 620185
Tryptic Soy Agar Medium for the growth of a wide variety of microorganisms.	E.P.	500 g 100 g	610052 620052

Dehydrated Culture Media			
Description		Packaging	Ref.
Tryptic Soy Broth Broth for the growth of a wide variety of microorganisms.	E.P.	500 g 100 g 5 kg	610053 620053 6100535
TSYE Agar (Tryptic Soy Yeast Extract Agar) General use medium for bacteria cultivation.	ISO 11290	500 g 100 g	610349 620349
Tryptone Glucose Extract Agar Medium for cultivation and count of microorganisms in food and water.	APHA	500 g 100 g	610089 620089
Tryptone Soya Yeast Extract Broth Medium for the confirmation of <i>Listeria monocytogenes</i> .	ISO 11290-1	500 g 100 g	610241 620241
Tryptone Sulfite Neomycin Agar Medium for <i>Clostridium perfringens</i> detection in food.		500 g 100 g	610074 620074
Tryptone Water Liquid medium for <i>Escherichia coli</i> detection in food and water.		500 g 100 g	610206 620206
Tryptophan Broth Broth of subcultivation for coliform differentiation, and for indole testing in the bacteriological examination of water.	ISO 9308-1 ISO 6579 ISO 6785	500 g 100 g	610197 620197
Tryptose Agar Medium for the growth of a wide variety of microorganisms.	APHA	500 g 100 g	610193 620193
Tryptose Broth Broth for cultivating and maintaining of fastidious microorganisms, including <i>Brucella spp.</i>		500 g 100 g	610233 620233
Tryptose Phosphate Broth Broth for the cultivation of a wide variety of fastidious microorganisms.		500 g	610364
Urea Agar Base Basal medium for enterobacteria urease test. Urea 40% supplement	ISO 6785 ISO 6579	500 g 100 g 10x5 mL	610107 620107 80292
Urea Broth Broth for the detection of microorganisms which metabolize urea.		500 g	610311
Violet Red Bile Agar MUG Selective medium for isolation and counting of coliforms and <i>Escherichia coli</i> .		500 g 100 g	610187 620187
Violet Red Bile Glucose Agar Selective medium for enterobacteria isolation and enumeration.	ISO 21528-2 E.P.	500 g 100 g	610059 620059
Violet Red Bile Lactose Agar Selective medium for coliforms isolation and enumeration.	ISO 4832	500 g 100 g 5 kg	610058 620058 6100585
Vogel Johnson Agar Selective basal medium for <i>Staphylococcus aureus</i> isolation. Potassium Tellurite 1% supplement Potassium Tellurite 3.5% supplement	U.S.P.	500 g 100 g 5x10 mL 5x10 mL	610186 620186 80022 80291
V.R.B. Agar w Glucose / Lactose Medium recommended for the detection and enumeration of enterobacteria.	E.P.	500 g	610336
VRE Agar Base Medium for isolating Vancomycin-resistant enterococci. Meropenem Supplement		500 g 10 vials	610501 81083
Water Detect Coliforms Chromogenic medium for detection of <i>E. coli</i> and other coliforms in water		500 g 100 g	610382 620382
W.L. Nutrient Agar Medium for the enumeration and cultivation of yeasts and bacteria in the microbiological control carried out in the brewing and other fermentation industries.		500 g 100 g	610234 620234
Wort Agar Medium for the cultivation, isolation and enumeration or enrichment of fungi, especially of yeasts.		500 g	610323
Wort Broth w/o NaCl Basal Broth for isolation and enumeration of yeasts and moulds. Glycerol supplement (Glycerol)		500 g 4x50 mL	610080 80021
Wurtz Lactose Agar Medium for enterobacteria isolation.		500 g 100 g	611210 621210
XLD Agar Differential and selective medium for pathogenic enterobacteria isolation.	ISO 21567 ISO 6579	500 g 100 g 5 kg	610060 620060 6100605
XLT4 Agar Base Basal medium for selective isolation of <i>Salmonella non typhi</i> and <i>paratyphi</i> .		500 g 100 g	610092 620092
XLT (Tergitol 4)		2x50 mL	80010
Yeast Extract Agar Medium for microorganisms enumeration in water.	ISO 6222	500 g 100 g	611016 621016
Yeast Glucose Chloramphenicol Agar Selective medium for yeasts and moulds isolation and enumeration.	ISO 21527-1	500 g 100 g	610070 620070

Dehydrated Culture Media			
Description		Packaging	Ref.
Yersinia ITC Broth Base Enrichment broth for detecting <i>Yersinia enterocolitica</i> in food and water.	ISO10273	500 g	610369
Yersinia ITC Supplement	ISO10273	10 vials	81092
Yersinia PSB Broth Enrichment broth for detecting <i>Yersinia enterocolitica</i> in food and water.	ISO 10273	500 g	610370
Yersinia Selective Agar Base Basal medium for selective isolation of <i>Yersinia enterocolitica</i> .	ISO 10273	500 g 100 g	610111 620111
Yersinia supplement (Cefsulodin, Irgasan, Novobiocin)	ISO 10273	10 vials	81039

Ingredients for culture media		
Description	Packaging	Ref.
Acid Hydrolysate of Casein Product obtained by acid hydrolysis of casein.	500 g	610496
Agar Bacteriological Purified Agar for culture media preparation.	500 g 100 g 5 kg	611001 621001 6110015
Beef Extract Beef extract obtained by enzymatic hydrolysis of animal proteins.	500 g 100 g 5 kg	610497 620497 6104975
Bile Bacteriological Bacteriological bile obtained by bile purification.	500 g 100 g	611367 621367
Bile Salt No. 3 Purified extract obtained of bovine bile.	500 g 5 kg	611901 6119015
Brain Heart Infusion Infusion obtained by enzymatic hydrolysis of brain and heart.	500 g	610495
Gelatin Bacteriological Purified gelatin for bacteriological use.	500 g	611002
Gelatin Pancreatic Peptone Peptone obtained by pancreatic digest of gelatin.	500 g	611503
Glucose Glucose for bacteriological use.	500 g 100 g 5 kg	611601 621601 6116015
Hemoglobin Autoclavable hemoglobin from beef blood, for preparing microbiological culture media.	500 g	611706
Lactose Lactose for bacteriological use.	500 g	610498
Liver Extract Dehydrated extract of bovine liver.	500 g 5 kg	612001 6120015
Malt Extract Dehydrated extract obtained from sprouted grain.	500 g 100 g 5 kg	611006 621006 6110065
Maltose Maltose for bacteriological use.	500 g	611602
Peptone Bacteriological Peptone obtained by enzymatic hydrolysis of animal tissues.	500 g 100 g 5 kg	611701 621701 6117015
Peptone Mycological Mixture of peptones for yeast culture.	500 g	612101
Polypeptone Pancreatic digest of casein and peptic digest of animal tissue combined in equal parts.	500 g	610492
Proteose Peptone Special peptone obtained by enzymatic hydrolysis of animal tissues.	500 g	612201
Skim Milk Pasteurized product resulting from the removal of water from nonfat milk (APHA).	500 g	610491
Sodium Chloride Sodium Chloride for bacteriological use.	500 g	611000
Sodium Desoxycholate Sodium desoxycholate for bacteriological use.	500 g 5 kg	612401 6124015
Sodium Selenite Sodium selenite for bacteriological use.	500 g	611003
Soy Peptone Soy peptone obtained by enzymatic hydrolysis of soy bean meal.	500 g 5 kg	612501 6125015
Sucrose Sucrose for bacteriological use.	500 g	611801
Tryptone Product obtained by pancreatic hydrolysis of casein.	500 g 5 Kg	611004 6110045
Tryptose Product obtained by enzymatic hydrolysis of casein, yeasts and animal tissues.	500 g	611008
Yeast Extract Yeast extract obtained by autolysis of <i>Saccharomyces cerevisiae</i> .	500 g 100 g 5 Kg	611005 621005 6110055

Freeze-dried selective and growth supplements		
Description	Packaging	Ref.
		
Ampicillin supplement Selective supplement for <i>Aeromonas hydrophila</i> isolation.	10 vials	81001
Bacillus cereus supplement (Polymyxin B) Selective supplement for <i>Bacillus cereus</i> isolation (ISO 7932).	10 vials	81016
Bordetella supplement (Cephalexin) Selective supplement for Bordetella isolation.	10 vials	81013
Brucella supplement (Polymyxin B, Bacitracin, Cycloheximide, Nalidixic acid, Nystatin, Vancomycin) Selective supplement for <i>Brucella</i> spp. isolation.	10 vials	81003
Campylobacter growth supplement (Sodium pyruvate, Sodium metabisulphite, Ferrous sulphate) Growth supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81050
Campylobacter Blaser Wang supplement (Polymyxin B, Trimethoprim, Amphotericin B, Cephalexin, Vancomycin) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81051
Campylobacter Butzler supplement (Bacitracin, Cycloheximide, Colistin sulphate, Sodium cephalozin, Novobiocin) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81015
Campylobacter C.C.D.A. supplement (Cefoperazone, Amphotericin B) Selective supplement for <i>Campylobacter</i> spp. isolation (ISO 10272-1).	10 vials	81037
Campylobacter C.T.V.N. supplement (Cefoperazone, Trimethoprim, Vancomycin, Cycloheximide) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81038
Campylobacter Karmali supplement (Sodium pyruvate, Cefoperazone, Cycloheximide, Vancomycin) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81036
Campylobacter Preston supplement (polymyxin B, Rifampicin, Cycloheximide, Trimethoprim) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81004
Campylobacter Skirrow supplement (Polymyxin B, Vancomycin, Trimethoprim) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81055
C.F.C. (Pseudomonas) supplement (Cetramide, Fucidin, Cephalotin) Selective supplement for <i>Pseudomonas</i> isolation.	10 vials	81049
Cefoxitin MRSA Supplement Selective supplement for the detection of Methicillin resistant <i>Staphylococcus aureus</i> .	10 vials	81087
Chloramphenicol supplement Selective supplement for yeasts and moulds isolation.	10 vials	81017
Clostridium difficile supplement (D-Cycloserine, Cefoxitin) Selective supplement for <i>Clostridium difficile</i> isolation.	10 vials	81007
Clostridium perfringens (T.S.C.) supplement (D-Cycloserine) Selective supplement for <i>Clostridium perfringens</i> isolation (ISO 7937).	10 vials	81011
C.N. (Pseudomonas) supplement (Cetrimide, Nalidixic acid) Selective supplement for <i>Pseudomonas aeruginosa</i> isolation (ISO 12780).	10 vials	81006
C.N.A. (Staf / Strep) supplement (Colistin sulphate, Nalidixic acid) Selective supplement for staphylococci, streptococci, pneumococci isolation.	10 vials	81048
Dermatophyte supplement (Chlortetracycline) Selective supplement for dermatophytes fungi isolation.	10 vials	81025
Gardnerella vaginalis supplement (Gentamicin, Nalidixic acid, Amphotericin B) Selective supplement for <i>Gardnerella vaginalis</i> isolation.	10 vials	81040
Gentamicin supplement Selective supplement for yeasts and moulds isolation.	10 vials	81033
Haemophilus supplement (Bacitracin, Vancomycin, Clindamycin) Selective supplement for <i>Haemophilus</i> spp. isolation.	10 vials	81014
Kanamycin supplement Selective supplement for group D streptococci isolation (APHA).	10 vials	81030
Kanamycin / Polymyxin B supplement Selective supplement for <i>Clostridium perfringens</i> supplement.	10 vials	81031

Freeze-dried selective and growth supplements		
Description	Packaging	Ref.
L.C.A.T. supplement (Lincomycin, Colistin sulphate, Amphotericin B, Trimethoprim) Selective supplement for pathogenic neisseriae isolation.	10 vials	81012
Legionella AB supplement Selective supplement for enrichment of <i>Legionella</i> spp.	10 vials	81099
Legionella growth supplement (L-Cysteine, Ferric pyrophosphate) Growth supplement for <i>Legionella</i> spp. isolation.	10 vials	80056
Legionella (B.C.Y.E.) growth supplement (L-Cysteine, Ferric pyrophosphate, ACES, Potassium hydroxide, Bromothymol blue, α-ketoglutarate) Growth supplement for <i>Legionella</i> spp. isolation (ISO 11731).	10 vials	81056
Legionella (B.M.P.A.) supplement (Anysomycin, Cefamandole, Polymyxin B) Selective supplement for <i>Legionella</i> spp. isolation.	10 vials	81002
Legionella (G.V.P.C.) supplement (Cycloheximide, Glycine, Polymyxin B, Vancomycin) Selective supplement for <i>Legionella</i> spp. isolation (ISO 11731).	10 vials	81008
Legionella (M.W.Y.) supplement (Anysomycin, Glycine, Polymyxin B, Vancomycin, Bromothymol blue, Bromocresol purple) Selective supplement for <i>Legionella</i> spp. isolation.	10 vials	81019
Listeria Fraser supplement (Ferric ammonium citrate) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10 vials	81028
Listeria Oxford supplement (Cefotetan, Fosfomicin) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10 vials	81027
Listeria Oxford WITH CEF/FOS/COL/CYC/ACR supplement (Cefotetan, Fosfomicin, Colistin, Cycloheximide, Acriflavine) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10 vials	81629
Listeria Palcam supplement (Polymyxin B, Cefotaxime, Acriflavine HCl) Selective supplement for <i>Listeria monocytogenes</i> isolation (ISO 11290-1).	10 vials	81026
Listeria supplement for enrichment (FDA/IDF) Selective supplement for <i>Listeria</i> isolation.	10 vials	81045
m-CP supplement Selective supplement for <i>Clostridium perfringens</i> isolation (D.L.vo 2/2/2001 np.31).	10 vials	81034
MUG supplement (Methylumbelliferyl-β-D-glucosonide) Supplement for <i>Escherichia coli</i> identification.	10 vials	81020
Novobiocin supplement Selective supplement for <i>Salmonella</i> isolation (ISO 6579).	10 vials	81021
O.N.P.G. 1.5% supplement (α-nitrophenyl-β-D-galactopyronoside) Supplement for enzyme β -galactosidase determination.	10 vials	81032
Preston Supplement ISO 10272-1 Selective supplement for <i>Campylobacter</i> spp. isolation (ISO 10272-1).	10 vials	81100
Rosolic acid supplement Selective supplement for faecal coliforms isolation.	10 vials	81029
R.P.F. supplement (Rabbit plasm, Bovine fibrinogen, Tripsin inhibitor, Potassium tellurite) Selective supplement for coagulase positive staphylococci identification (ISO 6888-2).	8 vials	81057
Schaedler supplement (Kanamycin, Vancomycin) Selective supplement for anaerobes isolation.	10 vials	81054
V.C.A.T. supplement (Vancomycin, Colistin sulphate, Amphotericin B, Trimethoprim) Selective supplement for pathogenic neisseria isolation.	10 vials	81041
V.C.N. supplement (Vancomycin, Colistin sulphate, Nystatin) Selective supplement for pathogenic neisseria isolation.	10 vials	81022
V.C.N.T. supplement (Vancomycin, Colistin sulphate, Nystatin, Trimethoprim) Selective supplement for pathogenic neisseria isolation.	10 vials	81024
Vancomycin (5 mg) Supplement Selective supplement for the isolation of <i>Cronobacter</i> (formerly called <i>Enterobacter</i>) <i>sakazakii</i> .	10 vials	81064
VITALEX growth supplement (Glutamine, Adenine, NAD) Growth supplement for fastidious microorganisms isolation.	10 vials	81023
Yersinia supplement (Cefsulodin, Irgasan, Novobiocin) Selective supplement for <i>Yersinia</i> spp. isolation (ISO 10273).	10 vials	81039

Liquid Supplements		
Description	Packaging	Ref.
		
Egg Yolk emulsion Supplement for lecithinase producers bacteria (<i>Staphylococcus</i> , <i>Bacillus</i> and <i>Clostridium</i>) identification (ISO 7932).	4x50 mL 6x100 mL	80219 80124
Egg Yolk Tellurite emulsion Egg Yolk Tellurite emulsion 20% Supplement for coagulase positive staphylococci isolation (ISO 6888-1).	4x50 mL 4x50 mL	80122 80123
Eugon LT 100 Supplement Supplement for the enrichment and detection of aerobic bacteria in cosmetic products	4x10 mL	80048
Glycerol supplement Supplement for culture media preparation.	4x50 mL	80021
Iodine solution Supplement for culture media preparation.	5x10 mL	81009
Lecithin supplement Emulsion of lecithin.	4x50 mL	80007
Letheen supplement Sterile emulsion of lecithin for Letheen media.	4x50 mL	80008
Listeria Fraser supplement (Ferric ammonium citrate) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10x5 mL	80304
Middlebrook 7 H 10 supplement (Albumina Bovina, Glucosio) Supplement for Middlebrook 7 H 10 medium preparation.	4x50 mL	81035
Muller Kauffmann supplement (Iodine, Potassium iodide, Brilliant green) Selective supplement for Muller Kauffmann Broth.	2x50 mL Iodine 50 mL B.G. 0.1%	80047
O.A. Listeria supplement Supplement for selective enrichment of <i>Listeria monocytogenes</i> .	(4+4) x 500 mL	81074
Potassium tellurite 1% supplement Potassium tellurite 3.5% supplement Selective supplement for pathogenic staphylococci isolation.	5x10 mL 5x10 mL	80022 80291
Tergitol 0.2% supplement Tergitol 1.5% supplement Supplement for coliforms isolation (ISO 9308-1).	5x5 mL 4x40 mL	80035 80042
TTC 1% supplement Supplement for bacteria differentiation.	5x10 mL	80300
Tween 80 supplement Supplement for culture media preparation.	2x50 mL	80031
Urea 40% supplement Supplement for urease activity of bacteria (ISO 6785) (ISO 6579).	10x5 mL 6/pack x 100 mL	80292 80110
Vitamin K 1% supplement Growth supplement for anaerobes isolation.	5x5 mL	80053
XLT4 supplement (Tergitol 4) Selective supplement for X.L.T. 4 Agar Base.	2x50 mL	80010

Animal blood and derivates			# Product with minimum request
Description	Packaging	Ref.	
Bovine serum	50 mL	83830	
Defibrinated horse blood	50 mL	83395	
	100 mL	83396	
	500 mL	83394	
Lysed horse blood	100 mL	83397	
Horse serum	100 mL	83810	
Defibrinated sheep blood	25 mL	83298	
	50 mL	83296	
	100 mL	83297	
	500 mL	83294	
Sheep Red Cells 10%	25 mL	83701	

Saline solutions and sputum treatment		
Description	Packaging	Ref.
DECONTAM KIT Kit for fluidification and decontamination of sputum samples to submit to culture for mycobacteria detection.	3 bottles x 9	80060
RINGER'S Solution (1/4X) Freeze-dried product for preparation of Ringer (1/4X) isotonic solution for bacteriological use.	10 vials	81059

Incubators		
Description	Packaging	Ref.
Mini Incubator EM4 ATL	1 unit	96750
Digital Thermoblock Incubator for ampoules and tubes. Range of temperature: 5 to 150 °C.	1 unit	96752
		
Aluminum block FOR STERILTEST Aluminum block with 20 holes.	1 unit	96753
Aluminum block FOR STERIL CONTROL Aluminum block with 12 holes.	1 unit	96754
Aluminum block FOR STRIP CONTROL MiRA Test ,MeRA Test Aluminum block with 12 holes.	1 unit	96755

Air sampler		
Description	Packaging	Ref.
		
PORTABLE AIR SAMPLER "55" Instrument for microbiological monitoring of environment with 55 mm Contact Plates, complete of: carrying case, air sampler, spare aluminum sampling head, remote control, battery charger, pencil torch for head sterilization, calibration certificate.	1 unit	96760
PORTABLE AIR SAMPLER "90" Instrument for microbiological monitoring of environment with 90 mm Petri dishes, complete of: carrying case, air sampler, spare aluminum sampling head, remote control, battery charger, pencil torch for head sterilization, calibration certificate.	1 unit	96770
SAS Aluminum Head Settle Plate (90 mm Petri dish)	1 unit	96783
SAS Replacement Battery	1 unit	96784



Radial Immunodiffusion Plates

Description	Packaging	ref.
EASY RID h-IgG	1 x 12 Tests	93001
EASY RID h-IgA	1 x 12 Tests	93002
EASY RID h-IgM	1 x 12 Tests	93003
EASY RID h-C3c	1 x 12 Tests	93004
EASY RID h-C4	1 x 12 Tests	93005
EASY RID h- Transferrin	1 x 12 Tests	93006
EASY RID h- Albumin	1 x 12 Tests	93007
EASY RID h- Apolipoprotein A1	1 x 12 Tests	93008
EASY RID h- Apolipoprotein B	1 x 12 Tests	93009
EASY RID h- Alpha-1-acid Glycoprotein	1 x 12 Tests	93010
EASY RID h- Fibrinogen	1 x 12 Tests	93011
EASY RID h- Antithrombin III	1 x 12 Tests	93012
EASY RID h-IgG Kappa light chain	1 x 12 Tests	93013
EASY RID h-IgG Lambda light chain	1 x 12 Tests	93014
EASY RID h- Alfa 1 Antitrypsin	1 x 12 Tests	93015
EASY RID h- Ceruloplasmin	1 x 12 Tests	93016
EASY RID h- Haptoglobin	1 x 12 Tests	93018
Ruler for RID		E20316

Radial Immunodiffusion Multiplates

Description	Packaging	ref.
Multiplate h-IgG/IgA/IgM	3 sectors x 4 Test	93104
Multiplate h-C3c/C4	2 sectors x 6 Test	93106
Multiplate h-ApoA1/ApoB	2 sectors x 6 Test	93110
Multiplate h- Kappa chain / Lambda chain	2 sectors x 6 Test	93115

Double Immunodiffusion Plates

Description	Packaging	ref.
Bence Jones Test Urine	20 x 2 Test	93201

Protein Control Serum

Description	Packaging	ref.
RID CONTROL SERUM IgG + IgA + IgM + C3c + C4 + Transferrin + α1 Acid Glycoprotein + Ig k chain + Ig λ chain + Ceruloplasmine + α1 Antitrypsin + Albumin + Haptoglobin + Apo A + Apo B	4 x 0.5 mL	940010

Antisera for microbial ID		
Description	Packaging	ref.
Shigella Antiserum Lyophilized polyvalent serum for serological <i>Shigella</i> characterization.	50 tests	96148
Yersinia enterocolitica Antiserum Test serum for <i>Yersinia enterocolitica</i> identification.	50 tests	96147

Latex kit for antibody determination		
Description	Packaging	ref.
ASO KIT Qualitative and semi-quantitative determination of antibodies to O-Streptolysin.	5 mL suspension 0.5 mL positive control 0.5 mL negative control	96101
CRP KIT Qualitative and semi-quantitative determination of C reactive protein.	5 mL suspension 0.5 mL positive control 0.5 mL negative control	96105
RF KIT Qualitative and semi-quantitative determination of reumathoid factor.	5 mL suspension 0.5 mL positive control 0.5 mL negative control	96110
WAALER ROSE KIT Qualitative and semi-quantitative determination of reumathoid factor.	5 mL suspension 0.5 mL positive control 0.5 mL negative control	96111
RPR KIT Qualitative and semi-quantitative determination of reagin joinable to Lue.	4 mL suspension 0.5 mL positive control 0.5 mL negative control	96115
TPHA KIT Qualitative and semi-quantitative determination of Treponema pallidum antibodies.	100 tests	96117
MONONUCLEOSI KIT Qualitative and semi-quantitative determination of eterophil antibodies joined to infectious mononucleosis.	2.5 mL red cells test 0.5 mL positive control 0.5 mL negative control	96120
ANTI n-DNA KIT Determination of antibodies to n-DNA joined to Systemic Lupus Eythematosus.	2.5 mL suspension 0.5 mL positive control 0.5 mL negative control	96125
CRP (C Reactive Protein) Latex	5 mL latex vial, liquid	96181
RF (Rheumatoid Factor) Latex	5 mL latex vial, liquid	96182

Pregnancy		
Description	Packaging	ref.
HCG URINE CARD 20 mUI/mL Rapid test for qualitative determination of human chorionic gonadotropin in urine.	50 Cards	96460

Rapid tests on card		
Description	Packaging	ref.
FECAL H.pylori ANTIGEN Card Rapid test for qualitative detection of <i>Helicobacter pylori</i> antigen in human fecal specimens.	20 Cards	97803
FECAL OCCULT BLOOD Card Research occult blood in feces (without diet).	20 Cards	96415/20
CHLAMYDIA Card	20 Cards	96450
FLAVISCREEN-HCV Card (CE mark)	25 Cards	96435
ONE STEP ROTA-ADENOVIRUS COMBO	20 cards	97802
ONE STEP ROTAVIRUS TEST	20 cards	97800
CEA Card	20 Cards	96485
PSA Card	40 Cards	96470
Ig E Tot. Card	20 Cards	96480
AFP – ALFA FETOCARD	20 Cards	96465
STREPTO A Card	30 Cards	96418
STREPTO B Card	25 Cards	97809
TREPONEMA PALLIDUM (Syphilis)	20 Cards	96495
HIV 1 & 2 Card (CE mark)	25 Cards	96440
Microalbumin Card (Urine)	20 Cards	96462
Tuberculosis Card	20 Cards	96468
TSH Card	20 Cards	96486
Myoglobin Card	20 Cards	96487
TROPONIN Card	20 Cards	96488
FERRITIN Card	20 Cards	96490
IM Card (mononucleosis infection)	20 Cards	96498
ONE STEP ADENOVIRUS TEST	20 Cards	97807
Clostridium difficile GDH	10 Cards	96316
Clostridium difficile Toxine A+B	10 Cards	96317
Giardia	10 Cards	96318
E. coli O157	10 Cards	96321
Gonorrhea	10 Cards	96441
Gardnerella vaginalis	20 Cards	96442
Trichomonas vaginalis	10 Cards	96443
HCV	50 Cards	96435/50
HIV 1 & 2 Card (CE mark)	50 Cards	96440/50
BJ Free k/λ dipstick urine	25 Cards	96444

Vidal-Wright kit

Description	Packaging	ref.
FEBRILE MULTITEST KIT Kit for <i>Salmonella</i> and <i>Brucella</i> disease determination.	8 x 5 mL	96015

MACRO bacterial suspensions for tube test

Description	Packaging	ref.
Salmonella typhi H Macro	20 mL	96001
Salmonella typhi O Macro	20 mL	96002
Salmonella typhi Total Macro	20 mL	96009
Salmonella paratyphi AH Macro	20 mL	96003
Salmonella paratyphi AO Macro	20 mL	96004
Salmonella paratyphi A Total Macro	20 mL	96010
Salmonella paratyphi BH Macro	20 mL	96005
Salmonella paratyphi BO Macro	20 mL	96006
Salmonella paratyphi B Total Macro	20 mL	96018
Salmonella paratyphi CH Macro	20 mL	96019
Salmonella paratyphi CO Macro	20 mL	96020
Salmonella paratyphi C Total Macro	20 mL	96021
Brucella Total Macro (only for human diagnostics)	20 mL	96007
Brucella abortus Macro (only for human diagnostics)	20 mL	96008
Brucella melitensis Macro (only for human diagnostics)	20 mL	96022
Brucella suis Macro (only for human diagnostics)	20 mL	96023
Proteus OX2 Macro	20 mL	96011
Proteus OX19 Macro	20 mL	96013
Proteus OXK Macro	20 mL	96012

SLIDE bacterial suspensions for slide test

Description	Packaging	ref.
Salmonella typhi H Slide	5 mL	96031
Salmonella typhi O Slide	5 mL	96032
Salmonella typhi Total Slide	5 mL	96033
Salmonella paratyphi AH Slide	5 mL	96034
Salmonella paratyphi AO Slide	5 mL	96035
Salmonella paratyphi A Total Slide	5 mL	96036
Salmonella paratyphi BH Slide	5 mL	96037
Salmonella paratyphi BO Slide	5 mL	96038
Salmonella paratyphi B Total Slide	5 mL	96039
Salmonella paratyphi CH Slide	5 mL	96040
Salmonella paratyphi CO Slide	5 mL	96041
Salmonella paratyphi C Total Slide	5 mL	96042
Brucella Total Slide (only for human diagnostics)	5 mL	96043
Brucella abortus Slide (only for human diagnostics)	5 mL	96044
Brucella melitensis Slide (only for human diagnostics)	5 mL	96045
Brucella Bengal Rose Slide (only for human diagnostics)	5 mL	96046
Proteus OX2 Slide	5 mL	96047
Proteus OX19 Slide	5 mL	96048
Proteus OXK Slide	5 mL	96049

Controls for Bacterial Suspensions

Description	Packaging	ref.
Positive control for Salmonella	0.5 mL	96096
Positive control for Proteus	0.5 mL	96097
Positive control for Brucella	0.5 mL	96098
Negative control	0.5 mL	96093

Dextrose Drink (not available in the European Union)



Description	Packaging	ref.
Dextrose Drink 50 orange 50 g glucose, orange flavor	12 x 200 mL	412512
Dextrose Drink 50 lemon 50 g glucose, lemon flavor	12 x 200 mL	412517
Dextrose Drink 75 orange 75 g glucose, orange flavor	12 x 200 mL	412511
Dextrose Drink 75 lemon 75 g glucose, lemon flavor	12 x 200 mL	412516
Dextrose Drink 100 orange 100 g glucose, orange flavor	12 x 200 mL	412510
Dextrose Drink 100 lemon 100 g glucose, lemon flavor	12 x 200 mL	412515
Dextrose Drink 50% 50% glucose solution	6 X 150 mL	402510

Swabs		
Description	Packaging	Ref.
STERILE DRYSWAB, FINE TIP, NARROW PLASTIC SHAFT Shaft: Flexible plastic for narrow orifices (urethral, pediatric). Rayon fiber: non-inhibitory to fastidious microorganisms.	100 units	85283
STERILE DRYSWAB, STRAIGHT ALUMINUM WIRE URETHRAL ENT Shaft: Controlled sampling from narrow orifices, ideal for male urethral sampling (urethral, ear). Rayon fiber: non-inhibitory to fastidious microorganisms.	100 units	85281
STERILE DRYSWAB, STANDARD PLASTIC SHAFT Shaft: standard plastic, suitable for most specimens (wound, nasal, throat, vaginal, skin). Alginate fiber: bud can be dissolved in weak organic acids to give total release of collected organisms, ideal for recovering organisms present in low numbers.	100 units	85282
STERILE DRYSWAB, ULTRA FINE FLEXIBLE TWIST-WIRE NASOPHARYNGEAL PEDIATRIC Shaft: ultra fine, twisted wire, for sampling nasopharyngeal region at back of nose (especially for whopping cough). Rayon fiber: non-inhibitory to fastidious microorganisms.	100 units	85280

Miscellaneous sterile disposable labware

Description	Packaging	ref.
Petri dishes 90 mm	600 units	96706
Petri dishes 60 mm	500 units	96705
Contact petri dishes 55 mm	500 units	96707
Two Sectors Petri dishes 90 mm	360 units	96709
Three Sectors Petri dishes 90 mm	500 units	96716
Disposable Plastic Loops 1 µL	100 units	96702
Disposable Plastic Loops 10 µL	100 units	96701
Sterile container with sampler for stool collection	120 units	96708
Tablets for CO ₂ generation	50 units	96740
Disposable Plastic Pasteur Pipettes	100 units	96703
Glass Test Tubes 16 mm x 160 mm	299 units	96717
Aluminum cups for test tubes	300 units	96718
L-Shape Spreaders	200 units	96763

© Liofilchem® Product Catalogue 2023

May 22nd, 2023
ref. 6553866

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.



Printed in Italy

Liofilchem® s.r.l.

Via Scozia, Zona Industriale
64026, Roseto degli Abruzzi (Te)
Italy
tel +39 0858930745
fax +39 0858930330
www.liofilchem.com